Prolactin

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INTENDED USE

The Prolactin kit is intended for the quantitative determination of prolactin in human serum.

Prolactin is a single-chain protein hormone with a molecular weight of about 23 000 Da. It is synthesized in the frontal lobe of pituitary gland. The determination of prolactin is a valuable tool in diagnostics of testicular and ovarian dysfunctions. Hyperprolactinemia in women may lead to galactorrhea, amenorrhea and other distortions in menstrual cycle. In men it may result in impotence or loss of libido. Measurement of circulating prolactin is used as a primary test for barrenness.

Pathological hyperprolactinemia takes place in the case of hypothyroidism, renal insufficiency and in patients with a pituitary gland tumor - prolactinoma. Physiological increase of prolactin level takes place during the gestation, lactation, in sleep and after physical and emotional stress.

PRINCIPLE OF TEST

The Prolactin kit is a “sandwich” type of solid-phase enzyme immunoassay, based on two monoclonal antibodies that are specific for different epitopes of prolactin molecule. One of these antibodies is conjugated with horseradish peroxidase; the other is immobilized on the inner surface of microwells. Prolactin molecules from the serum sample bind to both immobilized antibody and anti-prolactin-peroxidase conjugate. Then the wells are washed with wash buffer to remove any material not bound on the inner surface of the wells. Quantity of the bound conjugate is directly proportional to the prolactin level in sample.

During the incubation with substrate colour development occurs in microwells. Optical density of the solution in the wells is directly proportional to the quantity of the bound conjugate. The standard curve is plotted by using the prolactin concentrations in the calibrators (x-axis) and their corresponding OD values (y-

axis). The prolactin concentration of the specimen is directly read off from the standard curve.

PATIENTS' SERUM SAMPLES

Specimen collection and storage:

Blood is taken aseptically by venipuncture. After clotting, the serum is separated by centrifugation. Do not use plasma, hemolysed or lipemic serum and serum with sodium azide added as preservative.

Store specimens at +2...8°C for no more than 2 days; for longer storage it is recommended to aliquote and freeze them at -20°C or below. Avoid repeated freezing.

Preparation before use:

Prior to assay, allow the samples to reach room temperature (+18...25°C). Take care to agitate specimens gently in order to ensure homogeneity.

MATERIAL PROVIDED

MP	Microtiter plate: 12 breakable strips of 8 wells each (total 96 wells), coated with anti-prolactin monoclonal antibodies	1 plastic bag
CONJ	Conjugate. Contains anti-prolactin monoclonal antibodies conjugated with HRP	14 ml, ready to use
0-5 CAL	Prolactin calibrators (protein-based buffer containing known prolactin concentrations). Calibrators were standardized against 3 rd WHO International Standard 84/500. For exact prolactin concentrations, see vial labels.	6 vials , 0,5 ml each. Ready to use
CONTROL	Prolactin control (protein-based buffer containing known prolactin concentration). For exact range of prolactin concentration, see vial label.	0,5 ml, ready to use
SUB	TMB Substrate: 3,3',5,5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide	14 ml, ready to use
WASHB	Concentrated wash buffer, enough for 280 ml of solution	14 ml, concentrated
STOP	Stop reagent: 1N HCl solution	14 ml, ready to use

MATERIALS REQUIRED BUT NOT SUPPLIED

- a set of digital variable pipettes that cover volume range from 5 μ l to 5 ml, with appropriate disposable tips;
- 8-channel digital variable pipette that covers volume range up to 0.3 ml, with appropriate disposable tips;
- microplate shaker-incubator, which is able to maintain temperature +37°C and room temperature (+18...25°C) and shaking speed 500 to 800 rpm;
- microplate washer;
- automatic microplate reader with optical filter for 450 nm;
- volumetric cylinder, 200 ml;
- volumetric beaker, 300 ml;
- distilled or deionized water;
- gloves (rubber or plastic);
- trays for conjugate, substrate and stop-reagent.

APPLICATION AND STORAGE

Prolactin kit is designed for 96 determinations. This is sufficient for 40 unknowns, 6 calibrators, 1 control and 1 substrate (control of OD of TMB solution) in duplicates, provided that all the strips are used simultaneously.

Please take into consideration that calibrators should be measured in each separate run. It is also recommended to measure each time prolactin concentration in the control.

NOTE: If used partially, kit should be used within a month after opening.

The expiry date of the kit (12 months) is reported on the box label, expiry date for each component is indicated on the respective label.

Upon receipt, Prolactin kit should be stored at +2...8°C, preferably in the original kit box. Kit storage is permissible at the most +25°C for no more than 5 days.

If used for separate experiments, kit contents should be stored as follows:

- unused strips: in the original plastic bag with ziplock closure and close firmly. Store at +2...8°C until expiry date;
- opened vials with conjugate: at +2...8°C for no more than 1 month;
- calibrators and control: at +2...8°C for no more than 1 month after opening;
- opened vials with substrate: at +2...8°C for no more than 1 month;
- concentrated wash buffer: at +2...8°C until expiry date;
- wash buffer prepared for use: at room temperature (+18...25°C) for no more than 5 days, in a firmly closed container;
- stop-reagent: at +2...8°C until expiry date.

PREPARATION OF REAGENTS

Before the assay, allow all the kit components to reach room temperature and stir thoroughly.

MP The **microtiter plate** consists of a holder and breakable strips, packed in the plastic bag. Before opening, keep the bag at room temperature for at least 30 minutes. Open the bag and place required number of strips on holder. Put remaining strips back in the original plastic bag with ziplock closure and close firmly. Keep at +2...8°C until expiry date stated on the label.

WASHB Prepare the necessary volume of **wash buffer** by dilution the 20-fold concentrate with distilled or deionized water.

For example:

5 ml of concentrate + 95 ml of distilled water.

Mix thoroughly, avoiding foaming. Store at room temperature for no more than 5 days. Keep the prepared wash buffer firmly closed.

The rest of the concentrated wash buffer should be stored firmly closed at +2...8°C until expiry date.

CONJ Transfer conjugate into the tray at 1 ml per strip.

SUB Transfer substrate into the tray at 1 ml per strip directly before dispensing.

Substrate is sensitive to light.

STOP Transfer stop-reagent into the tray at 1 ml per strip before dispensing.

ASSAY PROCEDURE

All samples should be tested in duplicates.

1. Bring all the reagents to room temperature (+18...25°C) before use. Stir gently, avoid foaming. Assay scheme is on page 27.

2. Dispense **100 µl** of conjugate **CONJ** into each well, **except two wells «Blank»**.

3. Dispense

20 µl of prolactin **calibrators CAL** (0 to 5);

20 µl of prolactin **control CONTROL**;

20 µl of patients' samples into the respective wells.

Note: Total time of dispensing should not exceed 15 minutes, otherwise the test result may be unreliable, because the time of incubation will substantially vary for different samples.

4. Incubate strips for **1 hour while shaking (500–800 rpm)** at **+37°C**.

5. Decant, then wash each well 5 times with **300 µl** of **wash buffer (prepared from WASHB)**. Shake strips during 5–10 seconds each time before decanting. Make sure that after the last washing cycle the residual buffer is thoroughly aspirated from the wells.

6. Immediately add **100 µl** of **substrate SUB** into each well.

7. Incubate strips at **room temperature (+18...25°C)** in the **dark for 15–30 minutes**, depending on the color intensity.

8. Add **100 µl** of **stop-reagent STOP** to all the wells and shake well for **1–2 minutes** at a room temperature.


9. Read the optical density at a wavelength of **450 nm**. Colouring is stable at a room temperature in the course of 20 minutes.

DATA PROCESSING

Data processing is done by a computer-assisted analysis by plotting the mean OD values of the calibrators at 450 nm versus their respective prolactin concentrations using 4PL or 5PL fit (see typical standard curve). Arithmetical mean OD of the two «Blank» wells is used for zero setting of device.

Any extrapolation of the standard curve to prolactin concentration above the nominal value of calibrator 5 is not permitted.

TYPICAL EXAMPLE OF STANDARD CURVE *Do not use for evaluation of real assay data!*

EXPECTED VALUES

Serum samples of 120 apparently healthy people were assayed with Prolactin kit. The results are shown in the table. These limits should be considered as guidelines only.

Group	n	Prolactin concentration range (mIU/L)
Female	120	67-726
Male		105-540

It is highly recommended for each laboratory to determine its own reference range of prolactin concentrations.

PERFORMANCE CHARACTERISTICS OF THE ASSAY

Calibration:

Prolactin test kit was calibrated against the 3rd WHO International Standard 84/500.

Specificity:

No cross-interaction of both anti-prolactin monoclonal antibodies with LH, FSH, TSH, placental lactogen, and hGH was detected

Prolactin kit shows no **high dose hook effect** for prolactin concentrations up to 100 000 mIU/L. **High dose hook effect** was determined by spiking Calibrator 0 matrix with antigen.

Analytical sensitivity:

Analytical sensitivity of Prolactin assay, or the lowest detectable concentration that can be distinguished from zero calibrator, is 50 mIU/L. It is defined as mean OD of 10 replicates of Calibrator 0 plus two standard deviations.

Intra- and inter-assay variation:

To determine **intra-assay CV** 7 serum samples were assayed in 9 replicates each. The results are shown below.

Sample	Mean prolactin concentration, mIU/L	Intra-assay CV	
		SD	CV
HS 1	155	21,7	14,0%
HS 2	432	26,8	6,2%
HS 3	513	29,2	5,7%
HS 4	660	33,7	5,1%
HS 5	1234	55,5	4,5%
HS 6	1546	94,3	6,1%
HS 7	2567	179,7	7,0%

To determine **inter-assay CV** 8 serum samples were assayed 3 times by different operators with 1-week interval. Each specimen was assayed in 9 replicates. The results are shown below.

Sample	Mean prolactin concentration, mIU/L			Inter-assay precision	
	1 assay	2 assay	3 assay	SD	CV
HS 1	154	139	144	7,6	5,4%
HS 2	218	202	195	11,8	5,9%
HS 3	302	278	262	20,1	7,5%
HS 4	324	301	281	21,5	7,4%
HS 5	348	311	308	22,3	7,2%
HS 6	470	432	431	22,2	5,2%
HS 7	958	840	851	65,2	7,7%
HS 8	2410	2409	2049	208,1	9,3%

LIMITATION OF THE METHOD

Any clinical diagnosis should not be based on the results of in vitro diagnostic methods alone. To state a diagnosis, the physician is supposed to consider all the available clinical and laboratory findings.

SAFETY PRECAUTIONS

- **This kit is for in vitro diagnostic use only.** To obtain the reliable data an operator should thoroughly follow the instruction for use. This instruction for use is valid only for the present kit with the listed composition. Any exchange of kit components is not allowed by CE regulations.

- Do not use kits or components after expiry date stated on the label. Take into consideration stability period for reconstituted reagents.

- Do not mix or use together reagents from different lots, except substrate, stop-reagent and wash buffer.

- Stop-reagent is 1N HCl solution. Avoid contacts with skin and mucosa. In case of contact rinse affected region thoroughly with plenty of water and seek medical advice.

- Source materials of human origin that were used  in preparation of this kit were tested and found negative  for HBsAg, anti-HIV and anti-HCV antibodies. However, no known laboratory test guarantees the absence of these viral agents. Therefore, the kit components and patient's samples should be handled as potentially hazardous.

- As the kit contains potentially hazardous material,  the following precautions should be observed:

- wear disposable gloves when handling reagents and specimens;

- do not pipette by mouth;

- in the case of spilling, wipe up the spills promptly and wash the affected area thoroughly with decontaminant

- do not smoke, eat or drink while performing the assay.

- GLP with all general and individual regulations should be applied to the use of the kit.

ASSAY SCHEME

Wells	«Blank»	<div style="border: 1px solid black; display: inline-block; padding: 2px;">CAL</div> <div style="border: 1px solid black; display: inline-block; padding: 2px;">CONTROL</div>	Samples
Reagents			
CONJ	–	100 µl	100 µl
CAL CONTROL	–	20 µl	–
Samples	–	–	20 µl
Incubation №1	60 min, +37°C, 500–800 rpm		
WASHB (diluted)	5 x 300 µl		
SUB	100 µl	100 µl	100 µl
Incubation №2	15–30 min, +18...25°C, in the dark		
STOP	100 µl	100 µl	100 µl
Stirring	1–2 min, +18...25°C		
OD measuring	450 nm		
Calculations	Corresponding software		

December, 15, 2010

Astra Biotech GmbH
Im Biotechnologiepark, TGZ I
14943 Luckenwalde, Germany
Telefon: +49 (0) 3371 / 681-450
E-Mail: info@astrabiotech.de
www.astrabiotech.de