# Thrombosis kit

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### INTENDED USE

The **Thrombosis** diagnostic kit is intended for detection of polymorphisms in *F5* gene encoding coagulation factor V (+1691G>A polymorphism – Leiden mutation), in *F2* gene encoding prothrombin or coagulation factor II (+20210G>A polymorphism), and in the *MTHFR* gene encoding methylenetetrahydrofolate reductase enzyme (+677C>T polymorphism).

The most general factor of hereditary thrombophilia is +1691G/A mutation in *F5* gene (Leiden mutation) encoding coagulation factor V which is the cofactor at formation thrombin from prothrombin.

Polymorphism +1691G>A is located in exon 10 of the *F5* gene.

Prothrombin (coagulation factor II) is one of the main components of blood coagulation system. Mutation +20210G/A of the prothrombin gene increases the expression of this gene that leads to increase of the prothrombin concentration in blood in 1.5-2 times.

Polymorphism +20210G>A is located in 3'-untranslated region of the *F2* gene.

Hyperhomocysteinemia i.e. increase of the homocysteine concentration in blood is one of the risk factors of venous thrombosis. The most widespread polymorphism influencing the homocysteine level in blood is the +677C>T polymorphism in the *MTHFR* gene encoding methylenetetrahydrofolate reductase enzyme.

Polymorphism +677C>T is located in exon 4 of the *MTHFR* gene.

### PRINCIPLE OF THE TEST

The **Thrombosis** kit contains a set of primers for amplification of the specific sequence within the human genes *F5*, *F2* and *MTHFR*. Amplicons of the varying genotypes are determined by a subsequent specific restriction enzyme digestion.

In the case of +1691G>A (Leiden mutation): The major variant (G, protective) possesses a restriction site for the endonuclease whereas rare variant (A, pathogen) does not.

In the case of +20210G>A polymorphism: The major variant (G, protective) possesses a restriction site for the endonuclease whereas rare variant (A, pathogen) does not.

In the case of +677C>T polymorphism: The rare variant (T, pathogen) possesses a restriction site for the endonuclease whereas the major variant (C, protective) does not.

The products of enzymatic digestion are further detected by a gel electrophoresis.

### PATIENT'S SAMPLES

It is recommended to use patient's total DNA samples with a concentration range 10÷500 ng/μl. Genomic DNA can be extracted from blood, tissue samples or epithelium.

### Material Provided

<b>M1</b>	<b>PCR Mix M1 (+1691G&gt;A <i>F5</i>)</b>	Micro test-tube with green cover – 1.5ml	585 μl, ready to use
<b>M2</b>	<b>PCR Mix M2 (+20210G&gt;A <i>F2</i> and +677C&gt;T <i>MTHFR</i>)</b>	Micro test-tube with orange cover – 1.5ml	585 μl, ready to use
<b>P</b>	<b>DNA Polymerase enzyme (purple)</b>	Micro test-tube with yellow cover – 0.5ml	45 μl, ready to use
<b>C</b>	<b>Positive control for M1 and M2</b>	Micro test-tube with brown cover – 0.5ml	19 μl, ready to use
<b>E</b>	<b>Restriction enzyme</b>	Micro test-tube with purple cover – 0.5ml	39 μl, ready to use
<b>B</b>	<b>Buffer for restriction enzyme</b>	Micro test-tube with purple cover – 0.5ml	210 μl, ready to use
	Holder for test-tubes	1	

## MATERIALS REQUIRED BUT NOT SUPPLIED WITH THE KIT:

### Reagents:

- mineral oil (optional, for thermal cycler without heated lid);
- ethidium bromide (10 mg/ml);
- reagents for electrophoresis in agarose or acrylamide gels.

### Instruments:

- PCR-box;
- thermal cycler;
- centrifuge;
- vortex;
- pipettes (0.5 - 1000 µl) and sterile pipette tips;
- sterile micro tubes suitable for the thermal cycler in use;
- thermoblock;
- instruments for gel electrophoresis.

## SIZE AND STORAGE

The kit is designed for 112 test reactions that will be sufficient for testing of 40 clinical samples. 112 test reactions include: 80 mixes for patient samples (40 for M1 and 40 for M2), 16 for positive and 16 for negative controls.

The kit should be stored at -18...-22°C. The reagents are stable in the unopened micro tubes until the expiration date indicated on the package.

## ASSAY PROCEDURE

### Manufacturer's recommendations:

We recommend setting up blank protocols of PCR and enzyme digestion stages before the analysis.

We recommend using the diagnostic kit for 7 reactions in each PCR mix: Mix1 (analysis of +1691G/A polymorphism in the *F5* gene) and Mix2 (analysis of +20210G/A polymorphism in *F2* and +677C/T polymorphism in *MTHFR*).

Composition of single recommended series of experiments: 5 test samples + 1 positive control of PCR and enzyme digestion + 1 negative control. Each series is to be duplicated: 2x5 test samples + 2x1 reactions with positive controls + 2x1 reactions with negative controls according to two PCR mixes. In total: 14 reactions in one series of experiments.

The complete procedure consists of three steps:

1. Amplification;
2. Digestion of the PCR products with restriction enzyme;
3. Detection of the digested DNA by gel electrophoresis.

### 1. Amplification

#### 1.1. PCR carrying

- Defrost the M1, M2, C micro tubes, while keeping the rest of the kit in freezer. Stir and centrifuge for 10 seconds all the components before use.
- Take micro tubes twice as much as a number of samples to be analysed,
  - + 2 for control DNA samples (C1 and C2),
  - + 2 for negative controls of PCR (deionized water (O1 and O2); is not provided).
 and place them in a holder.
- Take P micro tube and place it in a holder. Put the kit in refrigerator.
- M1, M2, C and P are ready to use.
- Prepare Master Mixes according to the table below (where N is a number of reactions in each mix):

<i>PCR reagents</i>	<i>Master Mix (M1) volume</i>	<i>Master Mix (M2) volume</i>
PCR Mix (M1 or M2)	8.7 µl x (N+1)	8.7 µl x (N+1)
Polymerase P	0.3 µl x (N+1)	0.3 µl x (N+1)
<b>Total volume:</b>	<b>9.0 µl x (N+1)</b>	<b>9.0 µl x (N+1)</b>

- Stir micro tubes with Master Mixes and centrifuge for 10 seconds.
- Aliquot Master Mixes (**9 µl** for M1 and M2) into sterile micro tubes suitable for the thermal cycler, add **1 µl** of the **extracted DNA** (or control samples: DNA from C micro tubes into C1 and C2 micro tubes, and deionized water into O1 and O2 micro tubes) and centrifuge.
- If using thermal cycler without a heated lid overlay the mixes with 30 µl of mineral oil.
- Put the micro tubes into the thermal cycler.
- Perform the following amplification protocol:

PCR		
Stage	Temperature	Time
Initial Hold:	96 <sup>0</sup> C	2 min
30 cycles:	96 <sup>0</sup> C	10 s
	62 <sup>0</sup> C	15 s
	72 <sup>0</sup> C	20 s
Final Hold:	72 <sup>0</sup> C	10 s

## 2. Digestion of the Amplified DNA

**2.1.** Take B micro tube and place it in a holder. Keep the component at room temperature until final defrosting. Mix the component and centrifuge for 10 seconds before use.

**2.2.** Take E micro tube and place it on ice.

**2.3.** B and E are ready to use.

**2.4.** Prepare Master Mix according to the table below (where N is a number of reactions in each mix):

Digestion reagents	Master Mix volume
Buffer (B)	1.7µl x (2N+2)
Enzyme (E)	0.3 µl x (2N+2)
<b>Total volume:</b>	<b>2 µl x (2N+2)</b>

**2.5.** Stir Master Mix and centrifuge for 10 seconds.

**2.6.** Add **2 µl** of Master Mix into all PCR samples.

**2.7.** Stir Master Mixes and centrifuge for 10 seconds.

**2.8.** Put the micro tubes into a thermoblock and incubate for 2 hours (or night) at 65<sup>0</sup>C.

## 3. Detection of the Digested DNA and Interpretation of the Results

### 3.1. Gel electrophoresis

Electrophoresis can be run either in acrylamide or agarose gel:


- In case of using acrylamide, 8% gel is recommended. Load **5 µl** of each sample. Run the gel at voltage 15 V/cm. The dye band should pass 7 cm from the start. Stain the gel in the ethidium bromide solution (1 µg/ml).

- In case of using agarose, 2-3% gel is recommended. Load **5 µl** of each sample. The gel and the running buffer should contain ethidium bromide (0,1 µg/ml). Run the gel at voltage 15 V/cm. The dye band should pass 3-4 cm from the start.

**Note.** It's not necessary to add loading dye to the samples before electrophoresis.

### 3.2. Registration of the Digestion results


There are appropriate bands in the gel lanes with positive digestion control (C) (Pic.1 and 2).

**Picture 1.** Example of the digestion results for M1 Mix (polymorphism +1691G/A in F5) (All genotype variants are present on the picture).

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**Picture 2.** Example of the digestion results for M2 Mix (polymorphism +20210G/A in *F2* and +677C/T in *MTHFR*) (All genotype variants are present on the picture).

Explanation of the Picture 2:

Lane	Genotype: <i>F2</i> +20210G/A	Genotype: <i>MTHFR</i> +677C/T	Fragment length (bp)
1	pathogen:A / pathogen:A	protective:C / protective:C	311 / 125
2	pathogen:A / pathogen:A	protective:C / pathogen:T	311 / 257 / 125 / 54
3	pathogen:A / pathogen:A	pathogen:T / pathogen:T	/ 257 / 125 / 54
4	protective:G / pathogen:A	protective:C / protective:C	311 / 125 / 99
5	protective:G / pathogen:A	protective:C / pathogen:T	311 / 257 / 125 / 99 / 54
6	protective:G / pathogen:A	pathogen:T / pathogen:T	/ 257 / 125 / 99 / 54
7	protective:G / protective:G	protective:C / protective:C	311 / 99
8	protective:G / protective:G	protective:C / pathogen:T	311 / 257 / 99 / 54
9(C)	protective:G / protective:G	pathogen:T / pathogen:T	/ 257 / 99 / 54
10(O)	-	-	-

## PERFORMANCE CHARACTERISTICS OF THE ASSAY

### Specificity:

- Specificity of the kit is achieved by the selection of specific primers flanking the analyzed region of DNA molecule. Positive control is designed to validate quality and functionality of all the PCR and restriction components. The respective digested PCR products of positive control must appear on the gel as bands of the expected size (Pic.1 and 2).
- Negative control is supposed to represent absence of a DNA contamination. The respective gel lanes must not contain any PCR products (Pic.1 and 2).

### Sensitivity:

Sensitivity of the Thrombosis kit, i.e. lowest detectable concentration of total genomic DNA in a sample that can be analyzed, is 10 ng/μl.

### **LIMITATION OF THE METHOD**

Any clinical diagnosis should not be based on the results of in vitro diagnostic methods alone. To state a diagnosis, the physician is supposed to consider all the available clinical and laboratory findings.

All the resulting digestion products must correspond to the indicated sizes and match the positive controls. If this is not the case, the sample must be tested one more time or the complete analysis must be repeated with freshly isolated DNA. If the positive control digestion products are not present, the analysis is to be considered incorrect.

### SAFETY PRECAUTIONS

- **This kit is for in vitro diagnostic use only.** The operator should thoroughly follow the manual to obtain the reliable data. This instruction manual is valid only for the present kit with the listed contents. Any exchange of the kit components is not allowed by CE regulations.
- Source materials of human origin that were used in preparation of this kit were tested and found negative for HBsAg, anti-HIV and anti-HCV antibodies. However, no known laboratory test guarantees the absence of these viral agents. Therefore, the kit components and patient's samples should be handled as potentially hazardous. 
- Test should only be performed by skilled persons considering GLP (Good Laboratory Practice) guidelines.
- Don't use the kit after its expiration date.
- PCR technology is extremely sensitive. The amplification of a single DNA molecule generates millions of identical copies. Therefore, set up three separate working areas for a) sample preparation, b) PCR reagent preparation and c) restriction and DNA detection. For each working area a different set of pipettes should be reserved.
- Wear separate coats and gloves in each working area.
- Use sterile filter tips for pipetting and use special PCR pipettes for aerosol free pipetting.
- Routinely decontaminate your pipettes and the laboratory benches.
- Avoid aerosols.

April, 13, 2011

Astra Biotech GmbH  
Im Biotechnologiepark, TGZ I  
14943 Luckenwalde, Germany  
Telefon: +49 (0) 3371 / 681-450  
E-Mail: [info@astrabiotech.de](mailto:info@astrabiotech.de)  
[www.astrabiotech.de](http://www.astrabiotech.de)