17-OH-Progesterone kit

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2. INTENDED USE

The **17-OH-Progesterone** kit is provided for the **quantitative** determination of **17-OH-progesterone (17-OHP)** in human serum and plasma.

This test has 2 complete sets:

REF 20-03 for manual use;

REF 20-03 A must be used with ELISA automatic instrument “Alisei” manufactured by NEXT Level S.r.l. hereinafter referred to as analyser “Alisei”.

Instructions for use are described in:

section 13 for manual test,

section 14 for automated kit for “Alisei”.

Note 1: *Take into account that calibrators' nominals can be different for manual and automatic test kits.*

Note 2: *We guarantee applications of the test*

REF 20-03 A only on analyser “Alisei”,

REF 20-03 only for manual use.

While using a non-predefined method of use, it is under end user responsibility to make sure that it was appropriately tested for ELISA kits.

17-OHP is a steroid hormone with a molecular mass of 330.5 Da. It is secreted in the suprarenal gland cortex and sex glands.

A level of 17-OHP demonstrates circadian oscillations in parallel with cortisol secretion. Its maximum is detected in samples taken in the morning.

For non-pregnant females 17-OHP level in blood circulation depends on a menstrual period phase.

17-OHP is produced by mature ovarian follicle and by *corpus luteus*, as progesterone is. 17-OHP level increases essentially during pregnancy due to its buildup by adrenal cortex of foetus and by placenta.

Quantitative measurement of serum 17-OHP is a valuable tool for monitoring the activity of 21-hydroxylase of the suprarenal gland cortex. Usually lack of 21-hydroxylase is a result of congenital hyperplasia of suprarenal gland cortex and leads to over secretion of 17-OHP, the level of which rises in the peripheral blood. Here the deficiency of 11-hydroxylase leads only to moderate increase in concentration of 17-OHP. Thus, detection of this steroid hormone plays a very important role in differential diagnostics of congenital hyperplasia of suprarenal gland cortex.

3. PRINCIPLE OF THE TEST

The **17-OH-Progesterone** kit is a competitive solid phase enzyme immunoassay. During the incubation period 17-OHP of the tested samples and horseradish peroxidase (HRP) labeled 17-OHP binds to the antibodies coated onto the inner surface of the microplate wells until balance between them occurs. Separation of free and bound to antibodies 17-OHP and conjugate 17-OHP - peroxidase occurs while extracting the contents of the wells. The amount of bound conjugate is inversely proportional to the quantity of 17-OHP in the sample (Fig. 1).

During the incubation with TMB solution the colour is developing. The intensity of the colour is inversely proportional to the 17-OHP concentration in specimens.


Fig. 1 Assay scheme

The 17-OHP concentration in the patient sample is read from a standard curve that is processed in each assay.

4. STORAGE AND STABILITY OF THE KIT

The expiration date of the kit is printed on the box label; the expiration date for each component is printed on the respective label.

17-OH-Progesterone kit should be stored at +2...+8 °C upon receipt, preferably in the original kit box, until the expiration date. Storage at +25 °C is allowed but for no more than 5 days.

Shelf life of the kit is 18 months.

After initial opening the kit is stable for 12 months if stored at +2...+8 °C.

If used in several separate experiments, after initial opening kit contents should be stored as follows but never used longer than the expiration date:

- unused strips: in a firmly closed resealable zipper bag at +2...+8 °C until the expiration date;
- vials with conjugate, calibrators and control (ready-to-use): at +2...+8 °C for 12 months; vials with calibrators and control (reconstituted): at +2...+8 °C for no more than 1 month after opening;
- vial with substrate: at +2...+8 °C until the expiration date, protected from light;
- vials with concentrated Trial, concentrated wash solution and stop solution: at +2...+8 °C until the expiration date;
- wash solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days or at +2...+8°C for no more than 4 weeks, in a firmly closed bottle;
- Trial solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days, in a firmly closed bottle.

Damaged Test Kits

In case of any severe damage of the test kit or components, Astra Biotech GmbH has to be informed in writing, during one week after receiving the kit. Usage of severely damaged single components for a test run is not recommended.

5. SAMPLE COLLECTION AND STORAGE

Collect blood by venipuncture in a tube without anticoagulants (for serum) or evacuated tube with EDTA or

heparin (for plasma). Possibility to use other anticoagulants for the **17-OH-Progesterone kit** is not established.

Allow blood to clot for serum samples. Centrifuge the specimens to separate serum or plasma from the blood corpuscles.

Do not use haemolyzed (bright red) or lipaemic (milky) blood samples as well as samples containing sodium azide as preservative.

Store the serum and plasma samples at +2...+8 °C for no more than 5 days. Aliquot and freeze the samples for longer storage (≤ -20 °C). Avoid repeated freezing.

6. EXPECTED VALUES

Serum samples collected between 9 and 11 a.m. from 120 apparently healthy people, both males and females, between the ages of 21–45, were assayed with **17-OH-Progesterone kit**. The results are listed below. These limits should be considered as guidelines only.

Group	No	17-OHP concentration range (nmol/L)
<i>Female</i>	120	
Follicular phase		<0.3-2.06
Luteal phase		1.42-6.91
<i>Male</i>		0.4-8.3

It is highly recommended each laboratory to determine its own reference range of 17-OHP concentrations.

7. QUALITY CONTROL

It is recommended to use control samples according to the state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

8. REAGENT PREPARATION

Allow all the reagents to reach room temperature (+18...+25 °C), and then thoroughly stir.

MP Microplate

Keep the **microplate** at room temperature for at least 30 minutes before opening the bag. Place required number of strips onto strip holder. Place unused strips into the resealable zipper bag and reseal duly.

CAL **CONTROL** Calibrators and Control

Liquid **calibrators** and **control** are ready to use.

Prepare lyophilized calibrators and control as follows. Gently tap on the vial caps to knock off all the dry matter. Open the vials and carefully place the caps upside down on the clean dry surface. Add 0.5 mL of distilled or deionized water to each vial with lyophilized calibrators and control, close vials with the corresponding caps and leave for 10 min at room temperature without stirring. Then stir gently avoiding foaming, until the dry matter is completely dissolved. Leave for another 10 minutes at room temperature stirring gently periodically. Make sure that no dry matter is left on the caps and walls of the vials.

WASH P Wash solution

Prepare required volume of **wash solution** by dilution of the concentrate 20-fold with distilled or deionized water. For example:

5 mL of **WASH P 20X** + 95 mL of water

Mix thoroughly, avoid foaming.

SUB Substrate

Protect **substrate** from direct light.

9. SAMPLE PREPARATION

Allow samples to reach room temperature (+18...+25 °C). Stir samples gently in order to ensure homogeneity.

10. PERFORMANCE CHARACTERISTICS OF THE ASSAY**10.1. Calibration-Traceability:**

The **17-OH-Progesterone kit** was calibrated against the Working Standard, which had been manufactured by gravimetric method based on weighing purified synthetic 17-OHP into analyte-free matrix.

10.2. Specificity:

Cross-reaction of anti-17-OHP polyclonal antibodies with different steroids is shown below:

Steroid	Cross-reaction, %
C21-steroids	
17-OH-Progesterone	100
Pregnenolone	0.008
17-OH-Pregnenolone	0.626

17-OH-Pregnenolone sulfate	0.586
Progesterone	0.340
Deoxycorticosterone	0.018
11-Deoxycortisol	0.586
21-Deoxycortisol	1.690
Corticosterone	0.003
Aldosterone	< 0.001
Cortisol	0.003
C19-steroids	
Androstendione	0.002
Testosterone	0.004
5- α -Dihydrotestosterone	< 0.001
DHEA	< 0.001
DHEAS	< 0.001
C18-steroids	
17- β -Estradiol	< 0.001
17- α -Estradiol	< 0.001
Estrone	< 0.001
Estriol	< 0.001

10.3. Analytical Sensitivity:

Analytical sensitivity of **17-OH-Progesterone kit**, i.e. concentration, that can be distinguished from zero calibrator is 0.3 nmol/L (for manual and Alisei kit). It was defined as mean OD of 10 replicates of calibrator 0 minus 2 SD.

10.4. Measurement Range:

17-OH-Progesterone kit (manual and Alisei kit) was validated for measurement of 17-OHP concentration within the concentration diapason of 0.3 - 60 nmol/L.

10.5. Measurement Units:

In **17-OH-Progesterone kit** the concentrations of calibrators are specified in nmol/L. To convert into ng/mL, multiply the concentration in nmol/L by 0.33.

10.6. Intra- and Inter-Assay Variation (Precision):

For **Intra-Assay CV** determination, 6 serum samples were run, each in 2 replicates. The results are shown below.

Manual kit

Sample	Mean 17-OHP concentration, nmol/L	Intra-Assay CV	
		SD	CV, %
1	2.8	0.1	2.6
2	8.2	0.2	2.7
3	12.0	0.3	2.5
4	15.6	0.3	1.9
5	23.2	0.9	3.6
6	34	2.0	5.9

Alisei kit

Sample	Mean 17-OH-progesterone concentration, nmol/L	Intra-Assay CV	
		SD	CV, %
1	2.7	0.06	2.4
2	8.4	0.07	0,8
3	13.7	0.28	2.1
4	18.1	0.35	1.9
5	25.7	0.21	0.8
6	43.2	0.42	1.0

For **Inter-Assay CV** determination, 6 serum samples were assayed 3 times by different operators with 1-week interval. Each specimen was run in 2 replicates. The results are shown below.

Manual kit

Sample	Mean 17-OHP concentration, nmol/L			Inter-Assay CV	
	Assay 1	Assay 2	Assay 3	SD	CV, %
1	1.3	1.4	1.4	0.03	2.2
2	2.8	2.9	2.7	0.1	1.6
3	8.3	8.1	8.2	0.1	0.8
4	15.2	15.5	15.1	0.2	1.3
5	24.3	24.8	24.6	0.3	1.0
6	34.2	34.6	36.2	1.1	3.0

Alisei kit

Sample	Mean 17-OHP concentration, nmol/L			Inter-Assay CV	
	Assay 1	Assay 2	Assay 3	SD	CV, %
1	1.5	1.5	1.4	0.06	3.9
2	3.0	3.2	3.1	0.10	3.2
3	7.7	7.9	8.0	0.15	1.9
4	17.6	18.2	16.9	0.65	3.7
5	22.3	22.0	22.5	0.25	1.1
6	38.5	39.0	39.2	0.36	0.9

11. LIMITATION OF THE METHOD

Any clinical diagnosis should not be based on the results of in vitro diagnostic methods alone. For diagnosis establishment, a physician is supposed to consider all available clinical and laboratory findings.

12. SAFETY PRECAUTIONS

- **This kit is for in vitro diagnostic use only.** Operator should follow the manual closely in order to ensure reliable data. The manual is valid for the present kit only, within the

listed composition. Any substitution of kit components is not allowed by CE regulations.

- Do not use kit or its components after the expiration date indicated on the label. Take into account stability period for reconstituted reagents.
- Do not mix or use together reagents from different lots of the kit except substrate, stop solution and wash solution.
- Do not use substrate, stop solution and wash solution supplied by other vendors.
- Use only “P”-labeled Astra Biotech wash solution.
- Note that stop solution is 1 N HCl solution. Avoid contacts with skin and mucosa. In case of contact rinse affected area thoroughly with plenty of water and seek medical advice.
- Take into account the following common procedural notes:
 - always pipette reagents into wells immediately after washing procedure;
 - avoid contamination of the solutions;
 - in case of partial use of the kit, dispense only required volume of the reagent into the tray;
 - do not pour unused reagents back into the original vials;
 - avoid exposure to direct sunlight during incubations;
 - always pipette reagents in the same order to minimize reaction time differences between wells; the total dispensing time for the calibrators, control and samples must not exceed 15 min;
 - the incubation temperature for all the immunological reactions must be kept at +37 °C;
 - do not touch the bottom of the wells;

- calibrators should be measured in each separate assay. It is also recommended to measure each time 17-OHP concentration in the control;
- remove the adhesive foil carefully to avoid contamination and don't use the adhesive foil repeatedly.
- TMB solution should be colourless. Light colouring of solution is admissible. Avoid direct exposure of substrate to sunlight.
-  After usage strips, calibrators, control, specimens and all consumables which contacted with specimens during handling, storage or assay (tubes, vials, gloves, pipette tips etc.) should be collected separately and sterilized by autoclaving. Instead of autoclaving pipette tips may be sterilized by disinfectant treatment. After sterilization all components and expendable materials may be utilized as non-dangerous garbage. Other components of the kit should be discarded into conventional garbage.
-  During manual washing procedure do not discard the contents of the wells directly to drainage. Use a container with disinfectant solution.
-  As the kit contains potentially hazardous material, the following precautions should be taken:
 - do not smoke, eat or drink while performing the assay;
 - always use protective gloves;
 - never pipette material by mouth;
 - in case of spilling, wipe up the spills promptly and wash affected area thoroughly using decontaminant.
- GLP including all general and individual regulations should be applied for the kit usage.

As the kit contains irritant (**CONJ**, **CAL**, **CONTROL**), the following precautions should be observed:

- P261 - Avoid breathing spray;
- P272 - Contaminated work clothing should not be allowed out of the workplace;
- P280 - Wear protective gloves/protective clothing/eye protection;
- P302+P352 - IF ON SKIN: Wash with plenty of soap and water;
- P333+P313 - If skin irritation or rash occurs: Get medical advice/attention;
- P363 - Wash contaminated clothing before reuse;
- P501 - Dispose of contents/container in accordance with national regulation.

Precautionary statements according to Regulation EC No 1272/2008.

13. MANUAL TEST (REF 20-03)

13.1. Material Provided

MP	Microplate: 12 breakable 8-well strips (total 96 wells) coated with anti-17-OHP polyclonal antibodies	1 pcs
CONJ	Conjugate: solution containing 17-OHP conjugated with HRP	9 mL, ready to use
0-5 CAL	17-OHP calibrators: protein-based solutions or lyophilized preparations containing known 17-OHP concentrations – 0; 1; 2; 6; 10; 60 nmol/L (approximate values - do not use for evaluation of real assay data). The concentrations of calibrators may be different for Assay Protocol 1 and 2. For lot specific 17-OHP concentrations see values for respective Assay Protocol provided in the Quality Control Sheet.	6 vials, 0.5 mL each; ready to use or lyophilized
CONTROL	17-OHP control: protein-based solution or lyophilized preparation containing known 17-OHP concentration. The range of 17-OHP concentration may be different for Assay Protocol 1 and 2. For lot specific ranges of 17-OHP concentrations see values for respective Assay Protocol provided in the Quality Control Sheet.	0.5 mL, ready to use or lyophilized
SUB	Substrate (TMB solution): 3,3',5,5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide	14 mL, ready to use
WASH P	Wash solution P, 20X concentrated: surfactant in buffered saline, sufficient for	2x14 mL, concentrated

20X	preparation of 560 mL of solution	
STOP	Stop solution: 1 N HCl solution	14 mL, ready to use
	Adhesive foil	2x1 foil (optional)

13.2. Equipment and Materials Required but not provided

- 1-channel calibrated variable precision pipettes, with disposable tips;
- 8-channel calibrated variable precision pipette, with disposable tips;
- microplate incubator (+37 °C) or microplate incubator/shaker (+37 °C, shaking speed 650–800 rpm);
- manual or automatic equipment for rinsing wells;
- calibrated microplate reader (450 nm);
- vortex tube mixer;
- deionized or distilled water;
- graduated beaker and cylinder of appropriate volume;
- latex or plastic gloves;
- trays for pipetting reagents with 8-channel pipette;
- disinfectant;
- absorbent material (for manual wash).

13.3. Test Procedure

The **17-OH-Progesterone kit** is designed for 96 tests. This is sufficient for the **quantitative assay** of 40 unknowns, 6 calibrators, 1 control and 1 blank (OD of TMB solution) in duplicates, provided that all the strips are used simultaneously.

13.3.1. Assay Procedure

Note: *Please check in the Quality Control Sheet supplied with the kit for which Assay Protocol your kit has been validated. If your kit has been validated for Assay Protocol 1, please follow the Assay Protocol 1 described below. If your kit has been additionally validated for Assay Protocol 2, please consider the Supplement 1 to the Instructions for use provided with the kit.*

Assay Protocol 1

(see also Assay scheme to Assay Protocol 1, section 13.5)

Note: *Consider concentrations of calibrators and control range for Assay Protocol 1 provided in the Quality Control Sheet.*

- A.** Pipette **50 µL calibrators** **CAL** (0-5), **control** **CONTROL** and **patient's samples in duplicates** into the respective wells; **Leave wells A1-A2 empty for blank!**
- B.** Pipette **50 µL conjugate** **CONJ** into each well, **except wells A1-A2.**

Note: *total time of dispensing must not exceed 15 minutes, otherwise the test result may be unreliable, because the time of incubation will substantially vary for different samples.*

Note: *If an adhesive foil is used, put it now onto the plate to seal the cavities.*

C. Incubate strips for **30 minutes while shaking (650–800 rpm) at +37 °C.**

Note: *If an adhesive foil is used, remove it now from the plate.*

D. Wash 4 times, as described in section 13.3.2.

E. Pipette **100 µL substrate** **SUB** into each well (including blank); incubate **at room temperature in the dark for 15–30 minutes**, depending on the colour intensity, or **10 minutes while shaking (650–800 rpm) at +37 °C.**

F. Pipette **100 µL stop solution** **STOP** into each well (including blank) in the same sequence and at the same speed as used for dispensing TMB substrate. Shake for 1–2 min at room temperature.

G. Read **OD at 450 nm within 20 min.**

13.3.2. Wash Procedure

It is advisable to use an automatic microplate washer set at 4 wash cycles and a volume of 300 µL of wash solution per well per cycle.

If an automatic washer is not available, the wash procedure can be carried out manually as follows:

- remove the contents of the wells into a container with disinfectant;
- dispense 300 µL of wash solution (section 8) into each well, shake the plate carefully for 5–10 sec and remove the contents of the wells; repeat 4 times;
- strike the wells sharply on absorbent material to remove any liquid residue.

13.4. Data Processing

If the reader cannot be adjusted to zero using the substrate blank in wells A1-A2, subtract mean OD value of wells A1-A2 from all OD values before further calculations.

Example:

OD (Cal 0) measured = 2.28 and OD (blank) = 0.06;

OD (Cal 0) calculated = $2.28 - 0.06 = \underline{2.22}$

13.4.1. Data Reliability (for OD measured at 450 nm)

The data should meet the following criteria:

- average blank OD (in wells A1-A2) ≤ 0.100 ;
- average OD of Cal 0 ≥ 1.500 (after blank subtraction);
- control's concentration must fall within the acceptability range that is shown on the Quality Control Sheet.

If the data obtained do not meet the criteria, the results are considered unreliable and the test should be repeated.

13.4.2. Quantitative Determination

Specialized software for quantitative determination is recommended. Mean OD of calibrators are plotted versus their respective 17-OHP concentrations using **4PL fit** (see typical standard curve in Fig. 2). Calculate concentration of 17-OHP in samples using standard curve.

Any extrapolation of the standard curve to 17-OHP concentration above the nominal value of the calibrator 5 (approximately 60 nmol/L) is forbidden.

47/52


Fig. 2 Example of typical standard curve
Do not use for evaluation of real assay data!

13.5. Assay scheme to Assay Protocol 1

Reagents	Wells	«Blank»	CAL CONTROL	Samples
CAL CONTROL		–	50 µL	–
Samples		–	–	50 µL
CONJ		–	50 µL	50 µL
Incubation No.1	30 min, +37 °C, 650–800 rpm			
WASH P (diluted)	4 x 300 µL			
SUB		100 µL	100 µL	100 µL
Incubation No.2	15-30 min, +18...+25 °C, in the dark			
	10 min, +37 °C, 650–800 rpm			
STOP		100 µL	100 µL	100 µL
Stirring	1–2 min, +18...+25 °C			
OD measuring	450 nm			
Calculations	Corresponding software			

14. AUTOMATIC TEST (REF 20-03A)

14.1. Material Provided

MP	Microplate: 12 breakable 8-well strips (total 96 wells) coated with anti-17-OHP polyclonal antibodies	1 pcs
CONJ	Conjugate: solution containing 17-OHP conjugated with HRP	9 mL, ready to use
0-5 CAL	17-OHP calibrators: protein-based solution or lyophilized preparations containing known 17-OHP concentrations - 0; 1; 2; 6; 10; 60 nmol/L (approximate values - do not use for evaluation of real assay data). For lot specific 17-OHP concentrations see values provided in the Quality Control Sheet	8 vials, 0.5 mL each: CAL 1, 3 – 2 × 0.5 mL; CAL 0, 2, 4, 5 – 1 × 0.5 mL; ready to use or lyophilized
CONTROL	17-OHP control: protein-based solution or lyophilized preparation containing known 17-OHP concentration. For lot specific ranges of 17-OHP concentration see Quality Control Sheet	2 × 0.5 mL; ready to use or lyophilized
WASH P 20X	Wash solution P, 20X concentrated: surfactant in buffered saline, sufficient for preparation of 1000 mL of solution	50 mL, concentrated
SUB	Substrate (TMB solution): 3,3',5,5'-tetramethyl benzidine solution in citrate buffer containing hydrogen peroxide	14 mL, ready to use
STOP	Stop solution: 1 N HCl solution	19 mL, ready to use
TRIAL 5000X	Trial solution, 5000X concentrated: solution of detergent	It is delivered by separate order.

Note: *Extra vials of **CAL** 1, 3 and **CONTROL** are provided for recalibration on reference calibration curve. For detailed information see instructions for use for analyser “Alisei”.*

14.2. Equipment and Materials required but not provided

- Automatic analyser for ELISA kits on microplate (“Alisei”);
- multipurpose polypropylene tubes 12x75 volume 5.5 mL;
- **TRIAL** **5000X** Trial solution, 5000X concentrated²;
- deionized or distilled water;
- latex or plastic gloves.

14.3. Test Procedure

The **17-OH-Progesterone kit** is designed for 96 tests. This is sufficient for 88 unknowns, 6 calibrators, 1 control and 1 blank (OD of TMB solution) in monoplicates, provided that all the strips are used simultaneously.

14.3.1. Reagent Preparation

Allow all the reagents to reach room temperature, and then thoroughly stir.

Note: *Preparation of other reagents see in section 8.*

TRIAL Working solution of Trial for washing internal hydraulic lines and needles of automatic analyser

Prepare required volume of Trial solution by dilution of the concentrate 5000-fold with distilled or deionized water before analysis:

2 mL **TRIAL** **5000X** + 9998 mL water

Mix thoroughly, avoid foaming.

² Reagent is not included in the kit, it is delivered by separate order.

14.3.2. Assay Procedure for Automatic Test

While using for the procedure automatic instrument for microplates “Alisei”, refer to its relative manual.

The analysis on EIA analyser “Alisei” is completely automatic: pipetting of reagents, washing, incubation, OD measurement, analysis of results. The program of calculation 17-OHP concentration is hardwired in analyser memory.

14.4. Data Processing

14.4.1. Data Reliability (OD measurement at 450 nm)

See criteria in section 13.4.1.

14.5. Safety precautions

- If kit for automatic analyser “Alisei” is used in several separate experiments it is necessary to take reagents from automatic analyser “Alisei” immediately after pipetting them in the wells of all plates because liquid evaporates from vials. Put the reagents in refrigerator;
- Recalibration using calibration curve, obtained with a kit of any other lot, is not permitted.

August, 22, 2019

Astra Biotech GmbH
Rudower Chaussee 29,
12489 Berlin, Germany
Telefon: +49 (0)30 74696509
E-Mail: info@astrabiotech.de
www.astrabiotech.de