AFP kit

ELISA Testkit zur quantitativen Bestimmung
 von Alpha-Fetoprotein in humanem Blutserum

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2. INTENDED USE

The **AFP kit** is provided for the **quantitative** determination of **alpha-fetoprotein (AFP) in human serum**.

This test has 2 complete sets:

REF 22-02 for manual use,
REF 22-02 A must be used with ELISA automatic instrument "Alisei" manufactured by NEXT Level S.r.l. hereinafter referred to as analyser "Alisei"

Instructions for use:

Section 13 for manual test,
Section 14 for *analyser*.

Note 1: *Take into account that calibrators' nominals can be different for manual and automatic test kits.*

Note 2: *We guarantee applications of test:*

REF 22-02 A only on analyser "Alisei",
REF 22-02 only for manual use.

While using non predefined methods of use, it is under end user responsibility, to make sure that it was appropriately tested for ELISA kits.

AFP is a glycoprotein with a molecular mass of about 70 000 Da. During the ontogenesis AFP is produced mostly in yolk-sac and embryonal liver and, in less extent, in gastrointestinal tract. AFP is the main serum protein of early fetus. AFP concentration in fetal serum reaches its climax of up to 240 000 IU/mL between 12 and 15 weeks of gestation.

Then its concentration decreases and at the end of second year reaches the normal level (0-16 IU/mL). Normally, AFP concentration maintains at this level for life.

AFP is excreted through the fetus' kidneys into the amniotic fluid. From there it penetrates into maternal blood through placenta and umbilical cord. AFP level in serum of pregnant women continuously increases and reaches the maximum (up to 400 IU/mL) at the middle of the third trimester. Serum AFP measurement in pregnant women is an important method of early diagnostics of some inborn diseases. Besides, this method is widely used in obstetrics for diagnostics of multiple pregnancy, prenatal death and threatening abortion.

AFP measurement can be also used for the diagnostics and monitoring of different forms of cancer. For instance, high and maintaining AFP level (800...80 000 IU/mL) is often connected with primary hepatomas, testis teratomas and ovary tumors.

3. PRINCIPLE OF TEST

The **AFP kit** is a “sandwich” type of solid-phase enzyme immunoassay, based on two monoclonal antibodies that are specific for different epitopes of AFP molecule. One of these antibodies is conjugated with horseradish peroxidase (HRP); the other is coated onto the inner surface of microwells. The AFP molecules from the serum sample bound to both immobilized antibody and anti-AFP-peroxidase conjugate (Fig. 1). Then the wells are washed with wash solution to remove any material not bound to the inner surface of the wells. Quantity of the bound conjugate is directly proportional to AFP level in the sample.

During incubation with substrate the colour is developing. The intensity of the colour correlates to the concentration of AFP in specimens or calibrators. AFP concentration in the patient sample is read from a standard curve that is processed in each assay.


Fig. 1 Assay scheme

4. STORAGE AND STABILITY OF THE KIT

The expiration date of the kit is printed on the box label; expiration date for each component is printed on the respective label.

The **AFP kit** should be stored at +2...+8 °C upon receipt, preferably in the original kit box, until the expiration date. Storage at +25 °C is allowed but for no more than 5 days.

Shelf life of the kit is 18 months.

After initial opening the kit is stable for 12 months if stored at +2...+8 °C. If used in several separate experiments, after initial opening kit contents should be stored as follows but never used longer than the expiration date:

- unused strips: in a firmly closed resealable zipper bag at +2...+8 °C until the expiration date;
- vials with conjugate, calibrators and control (ready-to-use): at +2...+8 °C for 12 months; vials with calibrators and

control (reconstituted): at +2...+8 °C for no more than 1 month after opening;

- vial with substrate: at +2...+8 °C until the kit expiration date, protected from light;
- vials with concentrated Trial, sample diluent, concentrated wash solution and stop solution: at +2...+8 °C until the expiration date;
- wash solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days or at +2...+8°C for no more than 4 weeks, in a firmly closed bottle;
- Trial solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days, in a firmly closed bottle.

Damaged Test Kits

In case of any severe damage of the test kit or components, Astra Biotech GmbH has to be informed in writing, during one week after receiving the kit. Usage of severely damaged single components for a test run is not recommended.

5. SAMPLE COLLECTION AND STORAGE

Collect blood by venipuncture. After clotting, the serum is separated by centrifugation.

Do not use plasma, haemolyzed (bright red) or lipaemic (milky) serum samples as well as samples containing sodium azide as preservative. Store serum and plasma samples at +2...+8 °C for no more than 2 days. Aliquot and freeze the samples for a longer storage (≤ -20 °C). Avoid repeated freezing.

6. EXPECTED VALUES

Serum samples of 200 apparently healthy people (both males and females) were assayed with **AFP kit**, the AFP range was 0–14.4 IU/mL. These limits should be considered as guidelines only.

22 058 blood samples of pregnant women have been measured with the Astra Biotech AFP kit (2009-2010).

Pregnancy term on FDLMP, quantity of completed weeks	AFP			
	Median, ng/mL	Tolerant limits, ng/mL	Median, IU/mL	Tolerant limits, IU/mL
14	35	17-59	28	13-47
15	39	22-87	31	17-70
16	44	24-77	35	19-62
17	48	28-86	38	22-69
18	56	31-76	45	25-60
19	60	35-115	48	28-92
20	68	37-130	54	29-104

FDLMP – the first day of the last menstrual period

Tolerant limit – a limit in which with 95 % probability all physiological values of a marker lie.

It is highly recommended for each laboratory to determine its own reference ranges of AFP concentrations.

7. QUALITY CONTROL

It is recommended to use control samples according to the state and federal regulations. The use of control samples is advised to assure the day to day validity of results

8. REAGENT PREPARATION

Allow all the reagents to reach room temperature, and then thoroughly stir.

MP Microplate

Keep the **microplate** at room temperature (+18...+25 °C) for at least 30 minutes before opening the bag. Place required number of strips onto strip holder. Place unused strips into the resealable zipper bag and reseal duly.

CAL CONTROL Calibrators and Control

Liquid calibrators and control are ready to use.

Prepare lyophilized calibrators and control as follows. Gently tap on the vial caps to knock off all the dry matter. Open the vials and carefully place the caps upside down on the clean dry surface. Add 0.5 mL of distilled or deionized water to each vial with lyophilized calibrators and control, close vials with the corresponding caps and leave for 10 min at room temperature without stirring. Then stir gently avoiding foaming, until the dry matter is completely dissolved. Leave for another 10 minutes at room temperature stirring gently periodically. Make sure than no dry matter is left on the caps and walls of the vials.

WASH P Wash solution

Prepare required volume of **wash solution** by dilution of the concentrate 20-fold with distilled or deionized water. For example:

5 mL **WASH P 20X** + 95 mL water

Mix thoroughly, avoid foaming.

SUB Substrate

Protect **substrate** from direct light.

9. SAMPLE PREPARATION

Allow samples to reach room temperature. Stir samples gently in order to ensure homogeneity.

If expected AFP concentration in the serum sample is more than in calibrator 5, the sample should be additionally diluted **20-fold** (Sample 1) or **400-fold** (Sample 2) **DIL** in concordance with instructions for use for analyser “Alisei” or manually before analysis.

The example of manual sample dilution as follows:

Sample 1 (**20-fold dilution**):

380 μ L sample diluent **DIL** + 20 μ L serum sample

Sample 2 (**400-fold dilution**):

380 μ L sample diluent **DIL** + 20 μ L Sample 1

Vortex or mix thoroughly.

10. PERFORMANCE CHARACTERISTICS OF THE ASSAY

10.1. Calibration-Traceability:

The **AFP kit** was calibrated against the WHO 1st International Standard preparation 72/225.

10.2. Specificity:

No cross-reaction was detected between anti-AFP monoclonal antibodies used in the assay and serum albumin, chorionic gonadotropin, human placental lactogen and human immunoglobulin.

10.3. Analytical Sensitivity:

Analytical sensitivity of **AFP kit** i.e. concentration, that can be distinguished from zero calibrator, is 0.9 IU/mL. It was defined as mean OD of 10 replicates of calibrator 0 plus two standard deviations.

10.4. Measurement Units:

In **AFP kit** concentrations of calibrators are specified in IU/mL. To convert into ng/mL, multiply the concentration in IU/mL by 1.21.

10.5. Measurement Range:

The AFP kits (manual and Alisei kit) were validated for measurement of AFP within the concentration diapason (without dilution) of 0.9-300 IU/mL.

10.6. Hook Effect:

For AFP assay no high dose hook effect was detected for AFP concentration up to 12500 IU/mL.

10.7 Intra- and Inter-Assay Variation:

To determine **intra-assay CV** 8 serum samples were assayed in 9 replicates each. The results are shown below.

Manual kit

Sample	Mean AFP concentration, IU/mL	Intra-assay CV	
		SD	CV, %
1	4.8	0.47	9.8
2	24.3	1.85	7.6
3	46.8	2.46	5.2
4	49.2	2.19	4.4
5	61.3	4.12	6.7
6	72.5	3.14	4.3
7	86.3	6.75	7.8
8	167.6	4.99	2.9

Alisei kit

Sample	Mean AFP concentration, IU/mL	Intra-assay CV	
		SD	CV, %
1	5.1	0.21	4.1
2	10.7	0.83	7.8
3	24.5	1.67	6.8
4	57.1	3.26	5.7
5	72.2	2.59	3.6
6	93.4	3.21	3.4
7	125.0	6.40	5.1
8	189.9	9.94	5.2

To determine **inter-assay CV** 8 serum samples were assayed 3 times by different operators with 1-week interval. Each specimen was assayed in 9 replicates. The results are shown below.

Manual kit

Sample	Mean AFP concentration, IU/mL			Inter-assay CV	
	1 assay	2 assay	3 assay	SD	CV, %
1	5.0	4.3	5.2	0.46	9.6
2	26.2	22.0	24.4	2.09	8.6
3	48.0	43.8	48.6	2.61	5.6
4	49.3	46.9	51.3	2.22	4.5
5	63.5	56.5	63.7	4.1	6.7
6	72.7	69.3	75.0	2.89	4.0
7	89.4	79.3	90.4	6.14	7.1
8	166.2	167.0	172.2	3.28	1.9

Alisei kit

Sample	Mean AFP concentration, IU/mL			Inter-assay precision	
	1 assay	2 assay	3 assay	SD	CV, %
1	5.2	4.7	5.4	0.36	7.1
2	11.5	12.0	10.3	0.87	7.7
3	25.7	23.8	26.9	1.56	6.1
4	61.1	54.0	55.7	3.70	6.5
5	78.0	69.9	72.1	4.19	5.7
6	100.0	93.1	89.9	5.16	5.5
7	119.9	125.9	130.3	5.22	4.2
8	201.6	194.3	188.5	6.56	3.4

10.8. Dilution Parallelism of Serum Samples:

Serial dilutions of three human serum samples with predetermined AFP concentration were assayed with **AFP kit** with the following result:

Manual kit

Sample	Dilution with sample diluent	Measured concentration, IU/mL	Expected concentration, IU/mL	Measured/expected concentration ratio, %
1	Undiluted	273.2		
	1:2	133.4	136.6	98
	1:4	70.2	68.3	103
	1:8	36.7	34.15	107
	1:16	17.0	17.075	100
	1:32	8.90	8.54	104
2	Undiluted	179.4		
	1:2	83.9	89.7	93.5
	1:4	44.5	44.9	99.2
	1:8	21.7	22.4	96.7
	1:16	9.80	11.2	87.3
	1:32	5.41	5.6	96.7
3	Undiluted	288.0		
	1:2	135.1	144.0	93.8
	1:4	70.7	72.0	98.2
	1:8	38.6	36.0	107
	1:16	20.2	18.0	112

Alisei kit

Sample	Dilution	Measured concentration, IU/mL	Expected concentration, IU/mL	Observed/expected concentration ratio, %
1	Undiluted	270.5		
	1:2	140.9	136.6	103
	1:4	72.7	68.3	106
	1:8	32.1	34.15	94
	1:16	18.6	17.075	109
	1:32	7.8	8.54	91
2	Undiluted	177.1		
	1:2	81.6	89.7	91

Sample	Dilution	Measured concentration, IU/mL	Expected concentration, IU/mL	Observed/expected concentration ratio, %
	1:4	44.9	44.9	100
	1:8	22.9	22.4	102
	1:16	11.8	11.2	105
	1:32	6.1	5.6	109
3	Undiluted	284.6		
	1:2	148.2	144	103
	1:4	76.9	72	107
	1:8	35.4	36	98
	1:16	19.6	18	109

11. LIMITATION OF THE METHOD

Any clinical diagnosis should not be based on the results of in vitro diagnostic methods alone. For diagnosis establishment, a physician is supposed to consider all available clinical and laboratory findings.

12. SAFETY PRECAUTIONS

- **This kit is for in vitro diagnostic use only.** Operator should follow the manual closely in order to ensure reliable data. The manual is valid for the present kit only, within the listed composition. Any substitution of kit components is not allowed by CE regulations.

- Do not use kit or its components after the expiration date indicated on the label. Take into account stability period for reconstituted reagents.
- Do not mix or use together reagents from different lots of the kit except substrate, stop solution and wash solution.
- Do not use substrate, stop solution and wash solution supplied by other vendors.
- Use only “P”-labeled Astra Biotech wash solution.
- Note that stop solution is 1 N HCl solution. Avoid contacts with skin and mucosa. In case of contact rinse affected area thoroughly with plenty of water and seek medical advice.
- Take into account the following common procedural notes:
 - always pipette reagents into wells immediately after washing procedure;
 - avoid contamination of the solutions;
 - in case of partial use of the kit, dispense only required volume of the reagent into the tray;
 - do not pour unused reagents back into the original vials;
 - avoid exposure to direct sunlight during incubations;
 - always pipette reagents in the same order to minimize reaction time differences between wells; the total dispensing time for the calibrators, control and samples must not exceed 15 min;
 - the incubation temperature for all the immunological reactions must be kept at +37 °C;
 - do not touch the bottom of the wells;
 - calibrators should be measured in each separate assay. It is also recommended to measure each time AFP concentration in the control.

- remove the adhesive foil carefully to avoid contamination and don't use the adhesive foil repeatedly.
- TMB solution should be colourless. Light colouring of solution is admissible. Avoid direct exposure of substrate to sunlight.
-  Source materials of human origin used for kit components preparation were tested and found negative for HBsAg, anti-HIV and anti-HCV antibodies. However, none of known laboratory test guarantees absence of these viral agents. Therefore, all kit components and patient's samples should be handled as potentially hazardous.
-  After usage strips, calibrators, control, specimens and all consumables which contacted with specimens during handling, storage or assay (tubes, vials, gloves, pipette tips etc.) should be collected separately and sterilized by autoclaving. Instead of autoclaving pipette tips may be sterilized by disinfectant treatment. After sterilization all components and expendable materials may be utilized as non-dangerous garbage. Other components of the kit should be discarded into conventional garbage.
-  During manual washing procedure do not discard the contents of the wells directly to drainage. Use a container with disinfectant solution.
-  As the kit contains potentially hazardous material, the following precautions should be taken:
 - do not smoke, eat or drink while performing the assay;
 - always use protective gloves;

- never pipette material by mouth;
- in case of spilling, wipe up the spills promptly and wash affected area thoroughly using decontaminant.
- GLP including all general and individual regulations should be applied for the kit usage.
-  As the kit contains irritant (**CONJ**, **DIL**, **CAL**, **CONTROL**), the following precautions should be observed:
 - P261 - Avoid breathing spray;
 - P272 - Contaminated work clothing should not be allowed out of the workplace;
 - P280 - Wear protective gloves/protective clothing/eye protection;
 - P302+P352 - IF ON SKIN: Wash with plenty of soap and water;
 - P333+P313 - If skin irritation or rash occurs: Get medical advice/attention;
 - P363 - Wash contaminated clothing before reuse;
 - P501 - Dispose of contents/container in accordance with national regulation.

Precautionary statements according to Regulation EC No 1272/2008.

13. MANUAL TEST (REF 22-02)

13.1. Material Provided

MP	Microplate: 12 breakable 8-well strips (total 96 wells) coated with anti-AFP monoclonal antibodies	1 pcs
CONJ	Conjugate: solution containing anti-AFP monoclonal antibodies conjugated with HRP	14 mL, ready to use
0-5 CAL	AFP calibrators: protein-based solution or lyophilized preparations containing known AFP concentrations – 0; 10; 50; 100; 200; 300 IU/mL (approximate values - do not use for evaluation of real assay data). The concentrations of calibrators may be different for Assay Protocol 1 and 2. For lot specific AFP concentrations see values for respective Assay Protocol provided in the Quality Control Sheet.	6 vials, 0.5 mL each; ready to use or lyophilized
CONTROL	AFP control: protein-based solution or lyophilized preparation containing known AFP concentration. The range of AFP concentration may be different for Assay Protocol 1 and 2. For lot specific ranges of AFP concentrations see values for respective Assay Protocol provided in the Quality Control Sheet.	0.5 mL, ready to use or lyophilized
SUB	Substrate (TMB solution): 3,3',5,5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide	14 mL, ready to use
WASH P 20X	Wash solution P, 20X concentrated: surfactant in buffered saline, sufficient for preparation of 560 mL of solution	2x14 mL, concentrated

STOP	Stop solution: 1 N HCl solution	14 mL, ready to use
DIL	Sample diluent	10 mL, ready to use
	Adhesive foil	2x1 foil (optional)

13.2. Equipment and Materials Required but not provided

- 1-channel calibrated variable precision pipettes, with disposable tips;
- 8-channel calibrated variable precision pipette, with disposable tips;
- microplate incubator (+37 °C) or microplate incubator-shaker (37°C, shaking speed 500–800 rpm);
- manual or automatic equipment for rinsing wells;
- calibrated microplate reader (450 nm);
- vortex tube mixer;
- deionized or distilled water;
- graduated beaker and cylinder of appropriate volume;
- latex or plastic gloves;
- trays for pipetting reagents with 8-channel pipette;
- disinfectant;
- absorbent material (for manual wash).

13.3. Test Procedure

The **AFP kit** is designed for 96 tests. This is sufficient for the **quantitative assay** of 40 unknowns, 6 calibrators, 1 control and 1 blank (OD of TMB solution) in duplicates, provided that all the strips are used simultaneously.

13.3.1. Assay Procedure

Note: Please check in the *Quality Control Sheet* supplied with the kit for which Assay Protocol your kit has been validated. If your kit has been validated for Assay Protocol 1, please follow the Assay Protocol 1 described below. If your kit has been additionally validated for Assay Protocol 2, please consider the Supplement 1 to the Instructions for use provided with the kit.

Assay Protocol 1

(see also Assay scheme to Assay Protocol 1, section 13.5)

Note: Consider concentrations of calibrators and control range for Assay Protocol 1 provided in the *Quality Control Sheet*.

- A. Pipette **100 µL conjugate** **CONJ** into each well, **except wells A1-A2 (blank)**.
- B. Pipette **20 µL** of **calibrators** **CAL**, **control** **CONTROL** and prepared **patient's samples** in **duplicates** into the respective wells.

Leave wells A1-A2 empty for blank!

Note: *total time of dispensing must not exceed 15 minutes, otherwise the test result may be unreliable, because the time of incubation will substantially vary for different samples.*

Note: *If an adhesive foil is used, put it now onto the plate to seal the cavities.*

- C. Incubate for **60 minutes at +37 °C** while shaking (**500–800 rpm**).

Note: *If an adhesive foil is used, remove it now from the plate.*

- D. Wash 5 times, as described below.

- E.** Pipette **100 µL substrate** **SUB** into each well (including blank).
- F.** Incubate strips **at room temperature** (+18...+25 °C) **in the dark for 15-30 minutes**, depending on the colour intensity, or **10 minutes while shaking (500-800 rpm) at +37 °C**.
- G.** Pipette **100 µL stop solution** **STOP** to all the wells (including blank) in the same sequence and at the same speed as used for dispensing TMB substrate. **Shake for 1–2 min at room temperature**.
- H.** Read **OD at 450 nm within 20 min**.

13.3.2 Wash Procedure

It is advisable to use an automatic microplate washer set at 5 wash cycles and a volume of 300 µL of wash solution per well per cycle.

If an automatic washer is not available, the wash procedure can be carried out manually as follows:

- remove the contents of the wells into a container with disinfectant;
- dispense 300 µL of wash solution (section 8) into each well, shake the plate carefully for 5–10 sec and remove the contents of the wells; repeat 5 times;
- strike the wells sharply on absorbent material to remove any liquid residue.

13.4. Data Processing

If the reader cannot be adjusted to zero using the substrate blank in wells A1-A2, subtract mean OD value of wells A1-A2 from all OD values before further calculations.

Example:

OD (Cal 5) measured = 2.28 and OD (blank) = 0.06;

OD (Cal 5) calculated = $2.28 - 0.06 = \underline{2.22}$

13.4.1. Data Reliability (for OD measured at 450 nm)

The data should meet the following criteria:

- average blank OD (in wells A1-A2) ≤ 0.100 ;
- average OD of Cal 5 ≥ 1.0 (after blank subtraction).
- control's concentration must fall within the acceptability range that is shown on the Quality Control Sheet.

If the data obtained do not meet the criteria, the results are considered unreliable and the test should be repeated.

13.4.2 Quantitative Determination

Specialized software for quantitative determination is recommended. Mean OD of calibrators are plotted versus their respective AFP concentrations using **4PL or 5PL fit** (Fig. 2 typical standard curve).

Calculate concentration of AFP in samples using standard curve.

Any extrapolation of the standard curve to AFP above the nominal value of the calibrator 5 (approximately 300 IU/mL) is forbidden. In this case the sample should be diluted 50- or 2000-fold with sample diluent and re-tested. Multiply the measured concentration of pre-diluted samples by corresponding dilution factor.


Fig. 2 Example of typical standard curve
Do not use for evaluation of real assay data!

13.5 Assay scheme to Assay Protocol 1

Reagents	Wells	«Blank»	CAL CONTROL	Samples
CONJ	–	–	100 µL	100 µL
CAL CONTROL	–	–	20 µL	–
Samples	–	–	–	20 µL
Incubation No. 1	60 min, +37 °C, 500–800 rpm			
WASH P (diluted)	5 x 300 µL			
SUB	100 µL	100 µL	100 µL	100 µL
Incubation No. 2	15–30 min, +18...+25 °C, in the dark			
	10 min, +37 °C, 500–800 rpm			
STOP	100 µL	100 µL	100 µL	100 µL
Stirring	1–2 min, +18...+25 °C			
OD measuring	450 nm			
Calculations	Corresponding software			

14. AUTOMATIC TEST (REF 22-02 A)

14.1. Material Provided

MP	Microplate: 12 breakable 8-well strips (total 96 wells) coated with anti-AFP monoclonal antibodies	1 pcs
CONJ	Conjugate: solution containing anti-AFP monoclonal antibodies conjugated with HRP	14 mL, ready to use
0-5 CAL	AFP calibrators: protein-based solution or lyophilized preparations containing known AFP concentrations- 0; 10; 50; 100; 200; 300 IU/mL (approximate values - do not use for evaluation of real assay data). For lot specific AFP concentrations see values provided in the Quality Control Sheet. 	8 vials, 0.5 mL each: CAL 1, 4 – 2 × 0.5 mL; CAL 0, 2, 3, 5 – 1 × 0.5 mL; ready to use or lyophilized
CONTR OL	AFP control: protein-based solution or lyophilized preparation containing known AFP concentration. For exact range of AFP concentration see Quality Control Sheet	2 × 0.5 mL; ready to use or lyophilized
WASH P 20X	Wash solution P, 20X concentrated: surfactant in buffered saline, sufficient for preparation of 1000 mL of solution	50 mL, concentrated
SUB	Substrate (TMB solution): 3,3',5,5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide	14 mL, ready to use
STOP	Stop solution: 1 N HCl solution	19 mL, ready to use
TRIAL 5000X	Trial solution, 5000X concentrated: solution of detergent	It is delivered by separate order.
DIL	Sample diluent: protein-based solution	2 × 10 mL, ready to use

Note: Extra vials of **CAL** 1, 4 and **CONTROL** are provided for recalibration on reference calibration curve. For detailed information see instructions for use for analyser “Alisei”.

14.2. Equipment and Materials Required but not provided

- Automatic analyser for ELISA kits on microplate (analyser “Alisei”);
- multipurpose polypropylene tubes 12x75 volume 5.5 mL;
- **TRIAL 5000X** Trial solution, 5000X concentrated;
- deionized or distilled water;
- latex or plastic gloves.

14.3. Test Procedure

The **AFP kit** is designed for 96 tests. This is sufficient for 88 unknowns, 6 calibrators, 1 control and 1 blank (OD of TMB solution) in monoplicates, provided that all the strips are used simultaneously.

14.3.1. Reagent Preparation

Allow all the reagents to reach room temperature (+18...+25 °C), and then thoroughly stir.

Note: Preparation of other reagents see in section 8.

TRIAL Working solution of Trial for washing internal hydraulic lines and needles of automatic analyser.

Prepare required volume of **Trial solution** by dilution of the concentrate 5000-fold with distilled or deionized water before analysis:

2 mL **TRIAL 5000X** + 9998 mL water

Mix thoroughly, avoid foaming.

14.3.2. Assay Procedure for Automatic Test

While using for the procedure analyser “Alisei”, refer to its relative manual.

The analysis on analyser “Alisei” is completely automatic: pipetting of reagents, washing, incubation, OD measurement, analysis of results. The program of calculation AFP concentration is hardwired in analyser memory.

14.4. Data Processing

14.4.1. Data Reliability (OD measured at 450 nm)

See criteria in section 13.4.1.

14.5. Safety precautions

- If kit for analyser “Alisei” is used in several separate experiments it is necessary to take reagents from analyser “Alisei” immediately after pipetting them in the wells of all plates because liquid evaporates from vials. Put the reagents in refrigerator;
- Recalibration using calibration curve, obtained with a kit of any other lot, is not permitted.

August, 21, 2019

Astra Biotech GmbH
Rudower Chaussee 29,
12489 Berlin, Germany
Telefon: +49 (0)30 74696509
E-Mail: info@astrabiotech.de
www.astrabiotech.de