Anti-TPO kit

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2. INTENDED USE

The **Anti-TPO kit** is provided for the **quantitative** determination of **autoantibodies to thyroid peroxidase (Anti-TPO)** in human serum.

This test has 2 complete sets:

REF 24-06 for manual use;

REF 24-06 A must be used with ELISA automatic instrument “Alisei” manufactured by NEXT Level S.r.l. hereinafter referred to as analyser “Alisei”.

Instructions for use are described in:

section 13 for manual test,

section 14 for analyser

Note 1: *Take into account that calibrators' nominals can be different for manual and automatic test kits.*

Note 2: *We guarantee applications of the test*

REF 24-06 A *only on analyser “Alisei”,*

REF 24-06 *only for manual use.*

While using a non-predefined method of use it is under end user responsibility to make sure that it was appropriately tested for ELISA kits.

TPO is a key enzyme in the biosynthesis of thyroid hormones. Anti-TPO are mainly IgG and, less frequently, IgM. In healthy individuals, concentrations of Anti-TPO in serum can reach 30 U/mL. A level exceeding this limit is the sign of an autoimmune process that usually leads to thyroid dysfunction, mainly due to reduced thyroxin secretion.

Quantitative determination of Anti-TPO in serum is used for diagnostics of autoimmune thyroid diseases, such as idiopathic mixoedema, Hashimoto thyroiditis and Graves' disease.

3. PRINCIPLE OF TEST

The **Anti-TPO kit** detects Anti-TPO of IgG class. During the first incubation, Anti-TPO from test samples and calibrators binds to TPO coated onto the inner surface of the microplate wells. During the second incubation immobilized Anti-TPO reacts with conjugate (HRP-labeled anti-human IgG). Quantity of the bound conjugate is directly proportional to Anti-TPO concentration in tested sample (Fig. 1).

During the incubation with TMB solution the colour is developing. The intensity of the colour is directly proportional to Anti-TPO concentration in specimens. Anti-TPO concentration in the patient sample is read from a standard curve that is processed in each assay.


Fig. 1 Assay scheme

4. STORAGE AND STABILITY OF THE KIT

The expiration date of the kit is printed on the box label; the expiration date for each component is printed on the respective label. The **Anti-TPO kit** should be stored at +2...+8 °C upon receipt, preferably in the original kit box, until the expiration date. Storage at +25 °C is allowed, but for no more than 5 days.

Shelf life of the kit is 18 months.

After initial opening the kit is stable for 12 months, if stored at +2...+8 °C.

If used in several separate experiments, after initial opening kit contents should be stored as follows but never used longer than the expiration date:

- unused strips: in a firmly closed resealable zipper bag at +2...+8 °C until the expiration date;
- vials with conjugate, calibrators and control (ready-to-use): at +2...+8 °C for 12 months; vials with calibrators and control (reconstituted): at +2...+8 °C for no more than 1 month after opening;
- vial with substrate: at +2...+8 °C until expiration date, protected from light;
- vials with concentrated Trial, sample diluent, concentrated wash solution, stop solution and assay buffer: at +2...+8 °C until the expiration date;
- wash solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days or at +2...+8°C for maximal 4 weeks in a firmly closed bottle;
- Trial solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days in a firmly closed bottle.

Damaged Test Kits

In case of any severe damage of the test kit or components, Astra Biotech GmbH has to be informed in writing, during one week after receiving the kit. Usage of severely damaged single components for a test run is not recommended.

5. SAMPLE COLLECTION AND STORAGE

Collect blood by venipuncture. After clotting, the serum is separated by centrifugation.

Do not use plasma, haemolyzed (bright red) or lipaemic (milky) serum samples as well as samples containing sodium azide as preservative.

Store serum samples at +2...+8 °C for no more than 5 days. Aliquot and freeze the samples for a longer storage (≤ -20 °C). Avoid repeated freezing.

6. EXPECTED VALUES

324 serum samples collected between 9–11 a.m. from healthy individuals at the age 21-45 were assayed with **Anti-TPO kit**. In 99.4% of patients Anti-TPO concentration did not exceed 30 U/mL. This limit should be considered as guidelines only.

It is highly recommended for each laboratory to determine its own reference range of Anti-TPO concentrations.

7. QUALITY CONTROL

It is recommended to use control samples according to the state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

8. REAGENT PREPARATION

Allow all the reagents to reach room temperature (+18...+25 °C), and then thoroughly stir.

MP Keep the **microplate** at room temperature for at least 30 minutes before opening the bag. Place required number of strips onto strip holder. Place unused strips into the resealable zipper bag and reseal duly.

CAL CONTROL Calibrators and Control

Liquid calibrators and control are ready to use.

Prepare lyophilized calibrators and control as follows. Gently tap on the vial caps to knock off all the dry matter. Open the vials and carefully place the caps upside down on the clean dry surface. Add 0.5 mL of distilled or deionized water to each vial with lyophilized calibrators and control, close vials with the corresponding caps and leave for 10 min at room temperature without stirring. Then stir gently avoiding foaming, until the dry matter is completely dissolved. Leave for another 10 minutes at room temperature stirring gently periodically. Make sure that no dry matter is left on the caps and walls of the vials.

WASH P Prepare required volume of **wash solution** by dilution of the concentrate 20-fold with distilled or deionized water. For example:

5 mL **WASH P 20X** + 95 mL water

Mix thoroughly, avoid foaming.

SUB Protect **substrate** from direct light.

9. SAMPLE PREPARATION

Allow samples to reach room temperature (+18...+25 °C). Stir samples gently in order to ensure homogeneity.

Before the analysis all the serum samples (except calibrators and control) should be diluted 100-fold with sample diluent **DIL** in concordance with instructions for use for manual test or automatic test. Example of manual sample dilution as follows (Sample 1):

Sample 1 (100-fold dilution):

990 µL sample diluent **DIL** + 10 µL serum sample

If expected Anti-TPO concentration in the sample 1 is higher than in calibrator 5, it should be diluted 20-fold with sample diluent **DIL** (Sample 2).

Sample 2 (additional 20-fold dilution):

190 µL sample diluent **DIL** + 10 µL sample 1

10. PERFORMANCE CHARACTERISTICS OF THE ASSAY

10.1. Calibration-Traceability:

Anti-TPO test kit was calibrated against internationally approved NIBSC Research Standard 66/387.

10.2. Specificity:

High specificity of the assay is ensured by high purity TPO that was used for manufacture of the kit.

10.3. Analytical Sensitivity:

Analytical sensitivity of **Anti-TPO assay**, i.e. concentration that can be distinguished from calibrator 0, is 4 U/mL. It is defined as mean OD of 10 replicates of calibrator 0 plus two SD standard deviations.

10.4. Measurement Range:

The **Anti-TPO kit** was validated for measurement of Anti-TPO concentration within the concentration diapason (without dilution) of 4-500 U/mL (for manual test) and 4-1000 U/mL (for analyser “Alisei”).

10.5. Intra- and Inter-Assay Variation:

To determine **Intra-Assay CV**, 8 serum samples were assayed in 9 replicates each. The results are shown below.

Manual kit

Sample	Mean concentration of Anti-TPO, U/mL	Intra-Assay CV	
		SD	CV, %
1	22.945	0.4	1.6
2	51.285	3.6	7.0
3	100.91	0.8	0.9
4	234.27	12.7	5.4
5	383.27	9.1	2.4
6	78.033	2.6	3.3
7	34.335	2.1	6.1
8	444.85	28.2	6.3

Alisei kit

Sample	Mean concentration of Anti-TPO, U/mL	Intra-Assay CV	
		SD	CV, %
1	67.345	2.5	6.4
2	44.794	1.5	7.6
3	13.108	13.1	7.9
4	278.902	15.2	5.9
5	95.711	7.3	1.7
6	383.27	11.3	3.9
7	843.81	27.2	6.9
8	354.007	16.8	5.8

To determine **Inter-Assay CV** 8 serum samples were assayed 3 times by different operators with 1-week interval. Each specimen was assayed in 9 replicates. The results are shown below.

Manual kit

Sample	Mean concentration of Anti-TPO, U/mL			Inter-Assay CV	
	1 assay	2 assay	3 assay	SD	CV, %
1	34.3	33.8	34.3	0.3	0.9
2	23.1	23.5	26.0	1.6	6.5
3	96.0	93.4	103.3	5.1	5.2
4	134.3	123.3	133.3	6.0	4.6
5	155.5	146.8	156.1	5.2	3.4
6	157.8	146.8	147.1	6.3	4.2
7	101.2	92.0	102.1	5.6	5.7
8	318.0	315.7	307.0	5.8	1.8

Alisei kit

Sample	Mean concentration of Anti-TPO, U/mL			Inter-Assay precision	
	1 assay	2 assay	3 assay	SD	CV, %
1	14.9	15.2	14.8	1.9	5.3
2	34.6	34.5	37.8	0.6	1.1

3	55.2	56.1	55.0	2.6	3.5
4	77.8	75.3	72.6	1.1	1.2
5	90.1	89.3	91.4	4.5	3.4
6	137.7	128.9	131.7	22.3	6.6
7	344.6	312.4	355.1	31.5	4.2
8	722.1	754.7	785.0	0.2	1.4

11. LIMITATION OF THE METHOD

Any clinical diagnosis should not be based on the results of in vitro diagnostic methods alone. For diagnosis establishment, a physician is supposed to consider all available clinical and laboratory findings.

12. SAFETY PRECAUTIONS

- This kit is for in vitro diagnostic use only. Operator should follow the manual closely in order to ensure reliable data. The manual is valid for the present kit only, within the listed composition. Any substitution of kit components is not allowed by CE regulations.
- Do not use kit or its components after the expiration date indicated on the label. Take into account stability period for reconstituted reagents.
- Do not mix or use together reagents from different lots of the kit except substrate, stop solution and wash solution.
- Do not use substrate, stop solution and wash solution supplied by other vendors.
- Use only “P”-labeled Astra Biotech wash solution.

- Note that stop solution is 1 N HCl solution. Avoid contacts with skin and mucosa. In case of contact rinse affected area thoroughly with plenty of water and seek medical advice.
- Take into account the following common procedural notes:
 - always pipette reagents into wells immediately after washing procedure;
 - avoid contamination of the solutions;
 - in case of partial use of the kit, dispense only required volume of the reagent into the tray;
 - do not pour unused reagents back into the original vials;
 - avoid exposure to direct sunlight during incubations;
 - always pipette reagents in the same order to minimize reaction time differences between wells; the total dispensing time for the calibrators, control and samples must not exceed 15 min;
 - the incubation temperature for all the immunological reactions must be kept at +37 °C;
 - do not touch the bottom of the wells;
 - calibrators should be measured in each separate assay. It is also recommended to measure each time anti-TPO concentration in the control.
 - Carefully remove the protective adhesive film to avoid contamination and do not use the same protective adhesive film again.
- TMB solution should be colourless. Light colouring of solution is admissible. Avoid direct exposure of substrate to sunlight.
-  Source materials of human origin used for kit components preparation were tested and found negative for HBsAg, anti-HIV and anti-HCV antibodies. However, none of known laboratory test guarantees absence of these viral agents. Therefore, all kit components and patient's samples should be handled as potentially hazardous.

-  After usage strips, calibrators, control, specimens and all consumables which contacted with specimens during handling, storage or assay (tubes, vials, gloves, pipette tips etc.) should be collected separately and sterilized by autoclaving. Instead of autoclaving pipette tips may be sterilized by disinfectant treatment.

After sterilization all components and expendable materials may be utilized as non-dangerous garbage. Other components of the kit should be discarded into conventional garbage.

-  During manual washing procedure do not discard the contents of the wells directly to drainage. Use a container with disinfectant solution.
-  As the kit contains potentially hazardous material, the following precautions should be taken:
 - do not smoke, eat or drink while performing the assay;
 - always use protective gloves;
 - never pipette material by mouth;
 - in case of spilling, wipe up the spills promptly and wash affected area thoroughly using decontaminant.
- GLP including all general and individual regulations should be applied for the kit usage.

As the kit contains irritant (**CONJ**, **CAL**, **DIL**, **CONTROL**), the following precautions should be observed:

- P261 - Avoid breathing spray;
- P272 - Contaminated work clothing should not be allowed out of the workplace;
- P280 - Wear protective gloves/protective clothing/eye protection;

- P302+P352 - IF ON SKIN: Wash with plenty of soap and water;
- P333+P313 - If skin irritation or rash occurs: Get medical advice/attention;
- P363 - Wash contaminated clothing before reuse;
- P501 - Dispose of contents/container in accordance with national regulation.

Precautionary statements according to Regulation EC № 1272/2008.

13. MAUAL TEST (REF 24-06)

13.1. Material Provided

MP	Microplate: 12 breakable 8-well strips (total 96 wells) coated with human TPO	1 pcs
CONJ	Conjugate: solution containing monoclonal antibodies against human IgG conjugated with HRP	16 mL, ready to use
0-5 CAL	Anti-TPO calibrators: protein-based solution containing known Anti-TPO concentrations – 0; 25; 50; 100; 250; 500 U/mL (approximate values - do not use for evaluation of real assay data). The concentrations of calibrators may be different for Assay Protocol 1 and 2. For lot specific Anti-TPO-concentrations see values for respective Assay Protocol provided in the Quality Control Sheet.	6 vials, 0.5 mL each; ready to use or lyophilized
CONTROL	Anti-TPO control: protein-based solution or lyophilized preparation containing known Anti-TPO concentration. The range of Anti-TPO-concentration may be different for Assay Protocol 1 and 2. For lot specific ranges of Anti-TPO-concentration see values for respective Assay Protocol provided in the Quality Control Sheet.	0.5 mL, ready to use or lyophilized
SUB	Substrate (TMB solution): 3,3',5,5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide	14 mL, ready to use
WASH P 20X	Wash solution P, 20X concentrated: surfactant in buffered saline, sufficient for preparation of 560 mL of solution	2 x 14 mL, concentrated

STOP	Stop reagent: 1 N HCl solution	14 mL, ready to use
ASSAYB	Assay buffer	14 mL, ready to use
DIL	Sample diluent	50 mL, ready to use
	Adhesive foil	2x1 foil (optional)

13.2. Equipment and Materials Required but not provided

- 1-channel calibrated variable precision pipettes, with disposable tips;
- 8-channel calibrated variable precision pipette, with disposable tips;
- microplate incubator (+37 °C) or microplate incubator-shaker (+37 °C, 500–800 rpm);
- manual or automatic equipment for rinsing wells;
- calibrated microplate reader (450 nm);
- vortex tube mixer;
- deionized or distilled water;
- graduated beaker and cylinder of appropriate volume;
- latex or plastic gloves;
- trays for pipetting reagents with 8-channel pipette;
- disinfectant;
- absorbent material (for manual wash).

13.3. Test Procedure

The **Anti-TPO kit** is designed for 96 tests. This is sufficient for the **quantitative assay** of 40 unknowns, 6 calibrators, 1 control and 1 blank (OD of TMB solution) in duplicates, provided that all the strips are used simultaneously.

13.3.1. Assay Procedure

Note: *Please check in the Quality Control Sheet supplied with the kit for which Assay Protocol your kit has been validated. If your kit has been validated for Assay Protocol 1, please follow the Assay Protocol 1 described below. If your kit has been additionally validated for Assay Protocol 2, please consider the Supplement 1 to the Instructions for use provided with the kit.*

Assay Protocol 1

(see also Assay scheme to Assay Protocol 1, section 13.5)

Note: *Consider concentrations of calibrators and control range for Assay Protocol 1 provided in the Quality Control Sheet.*

All samples should be tested in duplicates.

A. Pipette **100 μ L assay buffer** **ASSAYB** into each well, **Leave wells A1-A2 empty for blank!**

B. Pipette **50 μ L of calibrators** **CAL** (0-5), **control** **CONTROL** and diluted **patient's samples in duplicates** into the respective wells.

Note: *total time of dispensing must not exceed 15 minutes, otherwise the test result may be unreliable, because the time of incubation will substantially vary for different samples.*

Note: *If an adhesive foil is used, put it now onto the plate to seal the cavities.*

C. Incubate strips for **30 minutes while shaking (500–800 rpm) at +37 °C.**

Note: *If an adhesive foil is used, remove it now from the plate.*

D. Wash 5 times, as described below.

E. Pipette **120 µL conjugate** **CONJ** into each well, **except wells A1-A2.**

F. Incubate strips for **30 minutes while shaking (500–800 rpm) at +37 °C.**

G. Wash 5 times, as described below.

H. Pipette **100 µL substrate** **SUB** into each well (including blank); incubate **at room temperature (+18...+25 °C) in the dark for 15-30 minutes** (depending on the colour intensity) or **10 minutes while shaking (500-800 rpm) at +37 °C.**

I. Pipette **100 µL stop solution** **STOP** into each well (including blank) in the same sequence and at the same speed as used for dispensing TMB substrate. **Shake for 1–2 min at room temperature.**

J. Read OD at **450 nm within 20 min.**

13.3.2. Wash Procedure

It is advisable to use an automatic microplate washer set at 5 wash cycles and a volume of 300 µL of wash solution per well per cycle. If an automatic washer is not available, the wash procedure can be carried out manually as follows:

- remove the contents of the wells into a container with disinfectant;
- dispense 300 µL of wash solution (section 8) into each well, shake the plate carefully for 5–10 sec and remove the contents of the wells; repeat 5 times;

- strike the wells sharply on absorbent material to remove any liquid residue.

13.4. Data Processing

If the reader cannot be adjusted to zero using the substrate blank in wells A1-A2, subtract mean OD value of wells A1-A2 from all OD values before further calculations.

Example:

OD (Cal 5) measured = 2.28 and OD (blank) = 0.06;

OD (Cal 5) calculated = $2.28 - 0.06 = \underline{2.22}$

13.4.1. Data Reliability (OD 450 nm)

The data should meet the following criteria:

- average blank OD (in wells A1-A2) ≤ 0.100 ;
- average OD of Cal 5 ≥ 1.0 (after blank subtraction);
- control's concentration must fall within the acceptability range that is shown in the Quality Control Sheet.

If the data obtained do not meet the criteria, the results are considered unreliable and the test should be repeated.

13.4.2. Quantitative Determination

Specialized software for quantitative determination is recommended. Mean OD of calibrators are plotted versus their respective Anti-TPO concentrations using **4PL** or **5PL fit** (Fig. 2 Typical standard curve). Calculate concentration of Anti-TPO in samples using standard curve.


Fig. 2 Example of typical standard curve.

Do not use for evaluation of real assay data!

Any extrapolation of the standard curve to Anti-TPO concentration above the nominal value of **calibrator 5** is forbidden. In this case, the sample should be additionally diluted 20-fold with sample diluent and re-tested. Multiply the measured concentration of additionally pre-diluted samples by additional dilution factor (20-fold).

13.5. Assay scheme to Assay Protocol 1

Wells	«Blank»	CAL CONTROL	Samples
Reagents			
ASSAYB	–	100 µL	100 µL
CAL CONTROL	–	50 µL	–
Samples	–	–	50 µL
Incubation No. 1	30 min, +37 °C, 500–800 rpm		
WASH P (diluted)	5 x 300 µL		
CONJ	–	120 µL	120 µL
Incubation No. 2	30 min, +37 °C, 500–800 rpm		
WASH P (diluted)	5 x 300 µL		
SUB	100 µL	100 µL	100 µL
Incubation No. 3	15–30 min, +18...+25 °C, in the dark		
	10 min, +37 °C, 500–800 rpm		
STOP	100 µL	100 µL	100 µL
Stirring	1–2 min, +18...+25 °C		
OD measuring	450 nm		
Calculations	Corresponding software		

14. AUTOMATIC TEST (REF 24-06 A)

14.1. Material Provided

MP	Microplate: 12 breakable 8-well strips (total 96 wells) coated with human TPO	1 pcs
CONJ	Conjugate: solution containing monoclonal antibodies against human IgG conjugated with HRP	17 mL, ready to use
0-5 CAL	Anti-TPO calibrators: protein-based solution containing known Anti-TPO concentrations – 0; 25; 100; 250; 500; 1000 U/mL approximate values - do not use for evaluation of real assay data). For lot specific Anti-TPO concentrations see values provided in the Quality Control Sheet	8 vials, 0.5 mL each: CAL 1,4- 2×0.5 mL; CAL 0, 2, 3, 5 – 1 × 0.5 mL; ready to use or lyophilized
CONTRO L	Anti-TPO control: protein-based solution or lyophilized preparation containing known Anti-TPO concentration. For lot specific ranges of Anti-TPO- concentration see Quality Control Sheet	2 × 0.5 mL; ready to use or lyophilized
WASH P 20X	Wash solution P, 20X concentrated: surfactant in buffered saline, sufficient for preparation of 1000 mL of solution	50 mL, concentrated
SUB	Substrate (TMB solution): 3,3',5,5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide	14 mL, ready to use
STOP	Stop reagent: 1 N HCl solution	19 mL, ready to use
TRIAL 5000X	Trial solution, 5000X concentrated: solution of detergent	It is delivered by separate order.
ASSAYB	Assay buffer	14 mL, ready to use
DIL	Sample diluent	2x 50 mL, ready to use

Note: Extra vials of **CAL** 1, 4 and **CONTROL** are provided for recalibration on reference calibration curve. For detailed information see instructions for use for analyser “Alisei”.

14.2. Equipment and Materials Required but not provided

- Automatic analyser for ELISA kits on microplate (analyser “Alisei”);
- multipurpose polypropylene tubes 12x75 volume 5.5 mL;
- **TRIAL** **5000X** Trial solution, 5000X concentrated;
- deionized or distilled water;
- latex or plastic gloves.

14.3. Test Procedure

The **Anti-TPO kit** is designed for 96 tests. This is sufficient for 88 unknowns, 6 calibrators, 1 control and 1 blank (OD of TMB solution) in monoplicates, provided that all the strips are used simultaneously.

14.3.1. Reagent Preparation

Allow all the reagents to reach room temperature (+18...+25 °C), and then thoroughly stir.

Note: Preparation of other reagents see in section 8.

TRIAL Working solution of Trial for washing internal hydraulic lines and needles of automatic analyser “Alisei”.

Prepare required volume of **Trial solution** by dilution of the concentrate 5000-fold with distilled or deionized water before analysis:

2 mL **TRIAL** **5000X** + 9998 mL water

Mix thoroughly, avoid foaming.

14.3.2. Assay Procedure for Automatic Test

While using for the analyser “Alisei”, refer to its relative manual.

The analysis on analyser “Alisei” is completely automatic: pipetting of reagents, washing, incubation, OD measurement, analysis of results. The program of calculation anti-TPO concentration is hardwired in analyser memory.

14.4. Data Processing

14.4.1. Data Reliability (OD 450 nm)

See criteria in section **13.4.1.**

14.5. Safety precautions

- If kit for analyser “Alisei” is used in several separate experiments it is necessary to take reagents from analyser “Alisei” immediately after pipetting them in the wells of all plates because liquid evaporates from vials. Put the reagents in refrigerator;
- Recalibration using calibration curve, obtained with kit of any other lot, is not permitted.

August 26, 2019

Astra Biotech GmbH
Rudower Chaussee 29,
12489 Berlin, Germany
Telefon: +49 (0) 30 74696509
E-Mail: info@astrabiotech.de
www.astrabiotech.de