# Ferritin kit

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## 2. INTENDED USE

The **Ferritin kit** is provided for the **quantitative** determination of **ferritin** in human serum.

This test has 2 complete sets:

**REF 40-01** for manual use;

**REF 40-01 A** must be used with ELISA automatic instrument “Alisei” manufactured by NEXT Level S.r.l. hereinafter referred to as analyser “Alisei”

Instructions for use are described in:

section 13 for manual test,

section 14 for automated kit for “Alisei”.

**Note 1:** *Take into account that calibrators’ nominals can be different for manual and automatic test kits.*

**Note 2:** *We guarantee applications of the test*

**REF 40-01 A** only on analyser “Alisei”,

**REF 20-01** only for manual use.

*While using a non-predefined method of use it is under end user responsibility to make sure that it was appropriately tested for ELISA kits.*

Ferritin is an intracellular protein. It makes soluble iron pool that is necessary for successful erythropoiesis. At the same time ferritin protects the organism from the toxic effect of iron (connected with catalysis of free radicals production).

Ferritin molecule consists of iron-containing core and protein shell (apoferritin) with a molecular mass of about 450 kDa. Apoferritin consists of 24 subunits of “light” (L) and “heavy” (H) types. Each ferritin molecule can bind up to 4500 iron atoms in hydroxide and phosphate complexes.

Ferritin is located mainly in spleen, liver and red bone marrow blood corpuscles; in small quantities it is also present in plasma. Ferritin concentration in plasma usually permits to estimate total iron reserve in the organism adequately.

Ferritin concentration at the moment of birth is high (up to 600 ng/mL). In the first months of life it falls and remains at the level of about 30 ng/mL till pubescence. Then ferritin concentration slowly begins to grow, the level reached at the age of 24–25 remains for life. Normal concentration of ferritin in serum is 50–250 ng/mL for men and 20–150 ng/mL for women.

Ferritin concentration below 10 ng/mL is a sign of asiderotic anemia. This parameter permits successfully differentiate asiderotic anemia from another types of anemia. Measurement of serum ferritin is also useful for monitoring of iron reserve in pregnant women, donors and patients subjected to regular haemodialysis. In the case of iron overloading ferritin concentration exceeds 400–500 ng/mL. It may reach several thousands ng/mL in the case of distinct hemochromatosis.


Concentration of serum ferritin increases in the case of infection, inflammatory processes (such as osteomyelitis and rheumatoid arthritis), acute and chronic liver disorders, leukemia, Hodgkin disease, breast cancer and some other oncopathologies.

### 3. PRINCIPLE OF THE TEST

The **Ferritin kit** is a “sandwich” type of solid-phase enzyme immunoassay, based on two monoclonal antibodies that are specific for different epitopes of ferritin molecule. One of these antibodies is conjugated with horseradish peroxidase (HRP); the other coated onto on the inner surface of microwells.

Ferritin molecules from the serum sample are bound to both immobilized antibody and anti-ferritin-peroxidase conjugate. Then the wells are washed with wash solution to remove any material not bound to the inner surface of the wells. Quantity of the bound conjugate is directly proportional to ferritin level in the sample (Fig.1).

During the incubation with TMB solution the colour is developing. The intensity of the colour is directly proportional to ferritin concentration in specimens. Ferritin concentration in the patient sample is read from a standard curve that is processed in each assay.

#### 4. STORAGE AND STABILITY OF THE KIT

The expiration date of the kit is printed on the box label; the expiration date for each component is printed on the respective label.

The **Ferritin kit** should be stored at +2...+8 °C upon receipt, preferably in the original kit box, until the expiration date. Storage at +25 °C is allowed but for no more than 5 days.

Shelf life of the kit is 18 months.

After initial opening, the kit is stable for 12 months if stored at +2...+8 °C.

If used in several separate experiments, after initial opening kit contents should be stored as follows but never used longer than the expiration date:

- unused strips: in a firmly closed resealable zipper bag at +2...+8 °C until the expiration date;
- vials with conjugate, calibrators and control (ready-to-use): at +2...+8 °C for 12 months; vials with calibrators and control (reconstituted): at +2...+8 °C for no more than 1 month after opening;
- vial with substrate: at +2...+8 °C until the expiration date, protected from light;
- vials with concentrated Trial, concentrated wash solution and stop solution: at +2...+8 °C until the expiration date;
- wash solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days or at +2...+8°C for no more than 4 weeks in a firmly closed bottle;
- Trial solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days in a firmly closed bottle.

### **Damaged Test Kits**

In case of any severe damage of the test kit or components, Astra Biotech GmbH has to be informed in writing, during one week after receiving the kit. Usage of severely damaged single components for a test run is not recommended.

## 5. SAMPLE COLLECTION AND STORAGE

Collect blood by venipuncture. After clotting, the serum is separated by centrifugation.

Do not use plasma, haemolyzed (bright red) or lipaemic (milky) serum samples as well as samples containing sodium azide as preservative.

Store serum samples at +2...+8 °C for no more than 2 days. Aliquot and freeze the samples for longer storage ( $\leq -20$  °C). Avoid repeated freezing.

## 6. EXPECTED VALUES

Serum samples taken from 184 healthy men (age 21–45); 140 healthy women (age 19–43) and 52 pregnant women (age 20–34) were assayed with **ferritin assay**. The results are shown in the table below. These limits should be considered as guidelines only.

<b>Study Group</b>	<b>Ferritin concentration range (ng/mL)</b>	<b>Mean ferritin concentration (ng/mL)</b>
Healthy men	22-346	158
Healthy women	10-147	73
Pregnant women (I .trimester)	55-90	61
Pregnant women (II. trimester)	25-74	32
Pregnant women (III. trimester)	10-16	12

It is highly recommended for each laboratory to determine its own reference range of ferritin concentrations.

## 7. QUALITY CONTROL

It is recommended to use control samples according to the state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

## 8. REAGENT PREPARATION

Allow all the reagents to reach room temperature, and then thoroughly stir.

### **MP** Microplate

Keep the **microplate** at room temperature (+18...+25 °C) for at least 30 minutes before opening the bag. Place required number of strips onto strip holder. Place unused strips into the resealable zipper bag and reseal duly.

### **CAL** **CONTROL** Calibrators and Control

Liquid calibrators and control are ready to use.

Prepare lyophilized calibrators and control as follows. Gently tap on the vial caps to knock off all the dry matter. Open the vials and carefully place the caps upside down on the clean dry surface. Add 0.5 mL of distilled or deionized water to each vial with lyophilized calibrators and control, close vials with the corresponding caps and leave for 10 min at room temperature (+18...+25 °C) without stirring. Then stir gently avoiding foaming, until the dry matter is completely dissolved. Leave for another 10 minutes at room temperature stirring gently periodically. Make sure that no dry matter is left on the caps and walls of the vials.

**WASH P Wash solution**

Prepare required volume of **wash solution** by dilution of the concentrate 20-fold with distilled or deionized water. For example:

5 mL **WASH P 20X** + 95 mL water

Mix thoroughly, avoid foaming.

**SUB SUB Substrate**

Protect **substrate** from direct light.

## 9. SAMPLE PREPARATION

Allow samples to reach room temperature. Stir samples gently in order to ensure homogeneity.

## 10. PERFORMANCE CHARACTERISTICS OF THE ASSAY

### 10.1 Calibration – Traceability:

The **Ferritin kit** was calibrated against the WHO 2<sup>nd</sup> International Standard 80/578.

### 10.2 Specificity:

According to Product specification provided by the supplier, no cross-interaction of both anti-ferritin monoclonal antibodies with serum albumin, haemoglobin, transferrin, alpha-fetoprotein and iron chloride was detected.

### 10.3 Analytical sensitivity:

Analytical sensitivity of **Ferritin kit** i.e. concentration that can be distinguished from zero calibrator, is 5 ng/mL. It is defined as mean OD of 10 replicates of Calibrator 0 plus two standard deviations.

### 10.4 Measurement Range:

The **Ferritin kits** (manual and Alisei kits) were validated for measurement of ferritin concentration within the concentration diapason of 5-1000 ng/mL.

### 10.5 Hook Effect:

For **Ferritin kit** high dose hook effect was not detected for concentrations up to 10 000 ng/mL. High dose hook effect was determined by spiking calibrator 0 matrix with antigen.

### 10.6 Intra- and Inter-Assay Variation:

To determine **Intra-Assay CV** 8 serum samples were assayed in 9 replicates each. The results are shown below.

#### Manual kit

Sample	Mean ferritin concentration, ng/mL	Intra-Assay CV	
		SD	CV, %
1	16.4	0.62	3.8
2	18.1	0.64	3.5
3	59.9	1.85	3.1
4	113	2.4	2.1
5	151	3.7	2.4
6	343	21.3	6.2
7	450	19.8	4.4
8	627	12.5	2.0

#### Alisei kit

Sample	Mean Ferritin concentration, ng/mL	Intra-Assay CV	
		SD	CV, %
1	11.1	0.78	7.0
2	21.4	1.13	5.3
3	75.7	2.45	3.2
4	137	4.9	3.6
5	225	5.3	2.4
6	450	9.8	2.2
7	664	26.2	3.9
8	897	38.9	4.3

To determine Inter-Assay CV 8 serum samples were assayed 3 times by different operators with 1-week interval. Each specimen was assayed in 9 replicates. The results are shown below.

#### Manual kit

Sample	Mean ferritin concentration, ng/mL			Inter-Assay CV	
	1 assay	2 assay	3 assay	SD	CV, %
1	15.2	13.7	14.1	0.79	5.5
2	44.0	44.8	43.2	0.90	2.0
3	51.3	49.1	54.2	2.52	5.0
4	81.0	80.5	80.1	0.50	0.6
5	120	124	120	2.1	1.8
6	277	275	283	4.2	1.5
7	366	375	382	8.0	2.2
8	583	588	609	13.8	2.4

#### Alisei kit

Sample	Mean ferritin concentration, ng/mL			Inter-Assay precision	
	1 assay	2 assay	3 assay	SD	CV,%
1	11.1	12.3	10.6	0.87	7.7
2	21.4	20.1	22.7	1.30	6.1
3	75.7	88.2	85.3	6.54	7.9

4	137	145	149	6.11	4.2
5	225	258	244	16.6	6.8
6	450	498	462	25.0	5.3
7	664	623	651	20.9	3.2
8	897	923	915	13.3	1.5

## 11. LIMITATION OF THE METHOD

Any clinical diagnosis should not be based on the results of in vitro diagnostic methods alone. For diagnosis establishment, a physician is supposed to consider all available clinical and laboratory findings.

## 12. SAFETY PRECAUTIONS

- **This kit is for in vitro diagnostic use only.** Operator should follow the manual closely in order to ensure reliable data. The manual is valid for the present kit only, within the listed composition. Any substitution of kit components is not allowed by CE regulations.
- Do not use kit or its components after the expiration date indicated on the label. Take into account stability period for reconstituted reagents.
- Do not mix or use together reagents from different lots of the kit except substrate, stop solution and wash solution.
- Do not use substrate, stop solution and wash solution supplied by other vendors.
- Use only “P”-labeled Astra Biotech wash solution.

- Note that stop solution is 1 N HCl solution. Avoid contacts with skin and mucosa. In case of contact rinse affected area thoroughly with plenty of water and seek medical advice.
- Take into account the following common procedural notes:
  - always pipette reagents into wells immediately after washing procedure;
  - avoid contamination of the solutions;
  - in case of partial use of the kit, dispense only required volume of the reagent into the tray; do not pour unused reagents back into the original vials;
  - avoid exposure to direct sunlight during incubations;
  - always pipette reagents in the same order to minimize reaction time differences between wells; the total dispensing time for the calibrators, control and samples must not exceed 15 min;
  - the incubation temperature for all the immunological reactions must be kept at +37 °C;
  - do not touch the bottom of the wells;
  - calibrators should be measured in each separate assay. It is also recommended to measure each time ferritin concentration in the control;
  - remove the adhesive foil carefully to avoid contamination and don't use the adhesive foil repeatedly.
- TMB solution should be colourless. Light colouring of solution is admissible. Avoid direct exposure of substrate to sunlight.
-  Source materials of human origin used for kit components preparation were tested and found negative for HBsAg, anti-HIV and anti-HCV antibodies. However, none of known laboratory test guarantees absence of these viral

agents. Therefore, all kit components and patient's samples should be handled as potentially hazardous.

-  After usage strips, calibrators, control, specimens and all consumables which contacted with specimens during handling, storage or assay (tubes, vials, gloves, pipette tips etc.) should be collected separately and sterilized by autoclaving. Instead of autoclaving pipette tips may be sterilized by disinfectant treatment.

After sterilization all components and expendable materials may be utilized as non-dangerous garbage. Other components of the kit should be discarded into conventional garbage.

-  During manual washing procedure do not discard the contents of the wells directly to drainage. Use a container with disinfectant solution.
-  As the kit contains potentially hazardous material, the following precautions should be taken:
  - do not smoke, eat or drink while performing the assay;
  - always use protective gloves;
  - never pipette material by mouth;
  - in case of spilling, wipe up the spills promptly and wash affected area thoroughly using decontaminant.
- GLP including all general and individual regulations should be applied for the kit usage.

As the kit contains irritant (**CONJ**, **CAL**, **CONTROL**), the following precautions should be observed:

- P261 - Avoid breathing spray;

- P272 - Contaminated work clothing should not be allowed out of the workplace;
- P280 - Wear protective gloves/protective clothing/eye protection;
- P302+P352 - IF ON SKIN: Wash with plenty of soap and water;
- P333+P313 - If skin irritation or rash occurs: Get medical advice/attention;
- P363 - Wash contaminated clothing before reuse;
- P501 - Dispose of contents/container in accordance with national regulation.

Precautionary statements according to Regulation EC No 1272/2008.

## 13. MANUAL TEST (REF 40-01)

### 13.1 Material Provided

<b>MP</b>	<b>Microplate:</b> 12 breakable 8-well strips (total 96 wells) coated with anti-ferritin monoclonal antibodies.	1 pcs.
<b>CONJ</b>	<b>Conjugate:</b> solution contains anti-ferritin monoclonal antibodies conjugated with HRP.	14 mL, ready to use
<b>0-5 CAL</b>	<b>Ferritin calibrators:</b> protein-based solution or lyophilized preparations containing known ferritin concentrations: 0;10;30;100;300;1000 ng/mL (approximate values - do not use for evaluation of real assay data). The concentrations of calibrators may be different for Assay Protocol 1 and 2. For lot specific ferritin concentrations see values for respective Assay Protocol provided in the Quality Control Sheet.	6 vials, 0.5 mL each ready to use or lyophilized
<b>CONTROL</b>	<b>Ferritin control:</b> protein-based solution or lyophilized preparation containing known ferritin concentration. The range of ferritin concentration may be different for Assay Protocol 1 and 2. For lot specific ranges of ferritin concentration, see values for respective Assay Protocol provided in the Quality Control Sheet.	0.5 mL, ready to use or lyophilized
<b>SUB</b>	<b>Substrate (TMB solution):</b> 3,3',5,5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide.	14 mL, ready to use
<b>WASH P 20X</b>	<b>Wash solution P, 20X concentrated:</b> surfactant in buffered saline, sufficient for	2x14 mL, concentrated

	preparation of 560 mL of solution.	
<b>STOP</b>	<b>Stop solution:</b> 1 N HCl solution.	14 mL, ready to use
	<b>Adhesive foil</b>	2x1 foil (optional)

### 13.2 Equipment and Materials Required but not provided

- 1-channel calibrated variable precision pipettes, with disposable tips;
- 8-channel calibrated variable precision pipette, with disposable tips;
- microplate incubator (+37 °C) or microplate incubator\_shaker (+37 °C, 500–800 rpm);
- manual or automatic equipment for rinsing wells;
- calibrated microplate reader (450 nm);
- vortex tube mixer;
- deionized or distilled water;
- graduated beaker and cylinder of appropriate volume;
- latex or plastic gloves;
- trays for pipetting reagents with 8-channel pipette;
- disinfectant;
- absorbent material (for manual wash).

### 13.3 Test Procedure

The **Ferritin kit** is designed for 96 tests. This is sufficient for 40 unknowns, 6 calibrators, 1 control and 1 blank (OD of TMB solution) in duplicates, provided that all the strips are used simultaneously.

#### 13.3.1 Assay Procedure

**Note:** *Please check in the Quality Control Sheet supplied with the kit for which Assay Protocol your kit has been validated. If your kit has been validated for Assay Protocol 1, please follow the Assay Protocol 1 described below. If your kit has been additionally validated for Assay Protocol 2, please consider the Supplement 1 to the Instructions for use provided with the kit.*

#### Assay Protocol 1

(see also Assay scheme to Assay Protocol 1, section 13.5)

**Note:** *Consider concentrations of calibrators and control range for Assay Protocol 1 provided in the Quality Control Sheet.*

**All samples should be tested in duplicates.**

- A.** Pipette **20  $\mu$ L** of calibrators **CAL**, control **CONTROL** and patient's samples in duplicates into the respective wells; **Leave wells A1-A2 empty for blank!**
- B.** Pipette **100  $\mu$ L conjugate CONJ** into each well, **except wells A1-A2.**

**Note:** *total time of dispensing must not exceed 15 minutes, otherwise the test result may be unreliable, because the time of incubation will substantially vary for different samples.*

**Note:** *If an adhesive foil is used, put it now onto the plate to seal the cavities.*

- C.** Incubate for **30 minutes at +37 °C while shaking (500–800 rpm)**.
- D.** Wash 5 times, as described below.
- E.** Pipette **100 µL substrate** **SUB** into each well (including blank); incubate strips at room temperature (+18...+25 °C) **in the dark for 15-30 minutes**, depending on the colour intensity, or **10 minutes while shaking (500-800 rpm) at +37 °C**.
- F.** Pipette **100 µL stop solution** **STOP** into each well (including blank) in the same sequence and at the same speed as used for dispensing TMB substrate. **Shake for 1–2 min at room temperature**.
- G.** Read OD at **450 nm within 20 min**.

### 13.3.2 Wash Procedure

It is advisable to use an automatic microplate washer set at 5 wash cycles and a volume of 300 µL of wash solution per well per cycle.

If an automatic washer is not available, the wash procedure can be carried out manually as follows:

- remove the contents of the wells into a container with disinfectant;
- dispense 300 µL of wash solution (section 8) into each well, shake the plate carefully for 5–10 sec and remove the contents of the wells; repeat 5 times;
- strike the wells sharply on absorbent material to remove any liquid residue.

## 13.4 . Data Processing

If the reader cannot be adjusted to zero using the substrate blank in wells A1-A2, subtract mean OD value of wells A1-A2 from all OD values before further calculations.

*Example:*

OD (Cal 5) measured = 2.28 and OD (blank) = 0.06;

OD (Cal 5) calculated =  $2.28 - 0.06 = \underline{2.22}$

### 13.4.1. Data Reliability (for OD measured at 450 nm)

The data should meet the following criteria:


- average blank OD (in wells A1-A2)  $\leq 0.100$ ;
- average OD of Cal 5  $\geq 1.5$  (after blank subtraction).
  - control's concentration must fall within the acceptability range that is shown in the Quality Control Sheet.

If the data obtained do not meet the criteria, the results are considered unreliable and the test should be repeated.

### 13.4.2. Quantitative Determination

Specialized software for quantitative determination is recommended. Mean OD of calibrators are plotted versus their respective ferritin concentrations using **4PL fit** (Fig. 2 Typical standard curve). Calculate concentration of ferritin in samples using standard curve.

Any extrapolation of the standard curve to ferritin above the nominal value of the calibrator 5 (approximately 1000 ng/mL) is forbidden.

*Fig. 2 Example of typical standard curve*  
**Do not use for evaluation of real assay data!**

### 13.5 Assay scheme to Assay Protocol 1

	Wells	«Blank»	<b>CAL</b> <b>CONTROL</b>	Samples
Reagents				
<b>CAL</b> <b>CONTROL</b>		–	20 µL	–
Samples		–	–	20 µL
<b>CONJ</b>		–	100 µL	100 µL
Incubation No.1	30 min, +37 °C, 500–800 rpm			
<b>WASH P</b> (diluted)	5 x 300 µL			
<b>SUB</b>	100 µL	100 µL	100 µL	100 µL
Incubation No.2	15–30 min, +18...+25 °C, in the dark			
	10 min, +37 °C, 500–800 rpm			
<b>STOP</b>	100 µL	100 µL	100 µL	100 µL
Stirring	1–2 min, +18...+25 °C			
OD measuring	450 nm			
Calculations	Corresponding software			

## 14. AUTOMATIC TEST (REF 40-01 A)

### 14.1. Material Provided

<b>MP</b>	<b>Microplate:</b> 12 breakable 8-well strips (total 96 wells) coated with anti-ferritin monoclonal antibodies.	1 pcs
<b>CONJ</b>	<b>Conjugate:</b> solution contains anti-ferritin monoclonal antibodies conjugated with HRP.	14 mL, ready to use
<b>0-5 CAL</b>	<b>Ferritin calibrators:</b> protein-based solution or lyophilized preparations containing known ferritin concentrations – 0;10;30;100;300;1000 ng/mL (approximate values - do not use for evaluation of real assay data). For lot specific ferritin concentrations see values provided in the Quality Control Sheet	<b>8 vials,</b> 0.5 mL each: <b>CAL</b> 1, 4 – 2 × 0.5 mL; <b>CAL</b> 0, 2, 3, 5 – 1 × 0.5 mL; ready to use or lyophilized
<b>CONTROL</b>	<b>Ferritin control:</b> protein-based solution or lyophilized preparation containing known ferritin concentration. . For lot specific ranges of ferritin concentration see Quality Control Sheet	2 × 0.5 mL; ready to use or lyophilized
<b>WASH P 20X</b>	<b>Wash solution P, 20X concentrated:</b> surfactant in buffered saline, sufficient for preparation of 1000 mL of solution	50 mL, concentrated
<b>SUB</b>	<b>Substrate (TMB solution):</b> 3,3',5,5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide	14 mL, ready to use
<b>STOP</b>	<b>Stop solution:</b> 1 N HCl solution	19 mL, ready to use
<b>TRIAL 5000X</b>	<b>Trial solution, 5000X concentrated:</b> solution of detergent	It is delivered by separate order.

**Note:** Extra vials of **CAL** 1, 4 and **CONTROL** are provided for recalibration on reference calibration curve. For detailed information see instructions for use for analyser "Alisei".

## 14.2. Equipment and Materials Required but not provided

- Automatic analyser for ELISA kits on microplate (analyser “Alisei”);
- multipurpose polypropylene tubes 12x75 volume 5.5 mL;
- **TRIAL 5000X** Trial solution, 5000X concentrated<sup>\*</sup>;
- deionized or distilled water;
- latex or plastic gloves.

## 14.3. Test Procedure

The **Ferritin kit** is designed for 96 tests. This is sufficient for 88 unknowns, 6 calibrators, 1 control and 1 blank (OD of TMB solution) in monoplicates, provided that all the strips are used simultaneously.

### 14.3.1. Reagent Preparation

Allow all the reagents to reach room temperature (+18...+25 °C), and then thoroughly stir.

**Note:** *Preparation of other reagents see in section 8.*

**TRIAL** Working solution of Trial for washing internal hydraulic lines and needles of automatic analyzer.

Prepare required volume of **Trial solution** by dilution of the concentrate 5000-fold with distilled or deionized water before analysis:

2 mL **TRIAL 5000X** + 9998 mL water

Mix thoroughly, avoid foaming.

---

<sup>\*</sup> reagent is not included in the kit, it is delivered by separate order

### **14.3.2. Assay Procedure for Automatic Test**

While using for the procedure analyser “Alisei”, refer to its relative manual.

The analysis on analyser “Alisei” is completely automatic: pipetting of reagents, washing, incubation, OD measurement, analysis of results. The program of calculation Ferritin concentration is hardwired in analyser memory.

## **14.4. Data Processing**

### **14.4.1. Data Reliability (OD measured at 450 nm)**

See criteria in section **13.4.1.**

## **14.5. Safety precautions**

- If kit for analyser “Alisei” is used in several separate experiments it is necessary to take reagents from analyser “Alisei” immediately after pipetting them in the wells of all plates because liquid evaporates from vials. Put the reagents in refrigerator;
- Recalibration using calibration curve, obtained with a kit of any other lot, is not permitted.

August 27, 2019

Astra Biotech GmbH  
Rudower Chaussee 29,  
12489 Berlin, Germany  
Telefon: +49 (0)30 74696509  
E-Mail: [info@astrabiotech.de](mailto:info@astrabiotech.de)  
[www.astrabiotech.de](http://www.astrabiotech.de)