Free T4 kit

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2. INTENDED USE

The **Free T₄** kit is provided for the **quantitative** determination of **free thyroxin (free T₄)** in human serum.

This test has 2 complete sets:

REF 24-03 for manual use,
REF 24-03 A must be used with ELISA automatic instrument “Alisei” manufactured by NEXT Level S.r.l. hereinafter referred to as analyser “Alisei”.

Instructions for use are described in:

Section 13 for manual test,
Section 14 for automated kit for “Alisei”.

Note 1: *Take into account that calibrators’ nominals can be different for manual and automatic test kits.*

Note 2: *We guarantee applications of test*

REF 24-03 A only on analyzer “Alisei”,
REF 24-03 only for manual use.

While using non predefined methods of use, it is under end user responsibility, to make sure that it was appropriately tested for ELISA kits.

Thyroxin (T₄) is one of the thyroid gland hormones with a molecular mass of 777 Da. Quantitative determination of serum free thyroxin is significant for evaluation of thyroid gland function.

The bulk of serum T₄ persists in protein-bound form. The main T₄-binding proteins are thyroxin-binding globulin (TBG), thyroxin-binding prealbumin and albumin.

Concentration of non-bound form of T_4 (free T_4) is very low, approximately 0.03% of the total concentration of circulating hormone.

Namely free T_4 is considered to be responsible for biological activity of the hormone.

Free T_4 concentration does not depend on the concentration of binding proteins and remains normal when TBG levels are either increased (e.g. due to congenital pathology, pregnancy, intake of estrogens and oral contraceptives) or decreased (e.g. due to congenital pathology, intake of androgens and salicylates, renal insufficiency) and when the binding capacity of the proteins is diminished due to salicylates treatment. So, free T_4 level adequately reflects the actual thyroid status. Concentration of free T_4 is above normal in patients with thyroid hyperfunction or after intake of thyroxin preparations; it is below normal in the case of thyroid hypofunction.

3. PRINCIPLE OF THE TEST

Free T_4 kit is a competitive solid phase enzyme immunoassay. During the incubation free T_4 of tested samples and horseradish peroxidase (HRP) labeled T_4 bind to the antibodies coated onto the inner surface of the microplate wells until balance between them occurs. Separation of free and bound to antibodies T_4 and conjugate T_4 -peroxidase occurs while extracting the contents of the wells. The amount of bound conjugate is inversely proportional to the quantity of free T_4 in the sample (Fig. 1).

During the incubation with TMB solution the colour is developing. The intensity of the colour is inversely proportional to the free T_4 concentration in specimens.

The free T_4 concentration in the patient sample is read from a standard curve that is processed in each assay.


Fig. 1 Assay scheme

4. STORAGE AND STABILITY OF THE KIT

The expiration date of the kit is printed on the box label; the expiration date for each component is printed on the respective label.

The Free T_4 kit should be stored at $+2...+8$ °C upon receipt, preferably in the original kit box, until the expiration date. Storage at $+25$ °C is allowed but for no more than 5 days.

Shelf life of the kit is 12 months.

After initial opening, the kit is stable until the expiration date if stored at $+2...+8$ °C.

If used in several separate experiments, after initial opening kit contents should be stored as follows:

- unused strips: in a firmly closed resealable zipper bag at +2...+8 °C until the expiration date;
- vials with conjugate, calibrators and control (ready-to-use): at +2...+8 °C until the expiration date; vials with calibrators and control (reconstituted): at +2...+8 °C for no more than 1 month after opening;
- vial with substrate: at +2...+8 °C until the expiration date, protected from light;
- vials with concentrated Trial, concentrated wash solution and stop solution: at +2...+8 °C until the expiration date;
- wash solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days or at +2...+8°C for no more than 4 weeks in a firmly closed bottle;
- Trial solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days in a firmly closed bottle

Damaged Test Kits

In case of any severe damage of the test kit or components, Astra Biotech GmbH has to be informed in writing, during one week after receiving the kit. Usage of severely damaged single components for a test run is not recommended.

5. SAMPLE COLLECTION AND STORAGE

Collect blood by venipuncture in a tube without anticoagulants. Allow blood to clot. Centrifuge the specimens to separate serum from the blood corpuscles.

Do not use plasma, haemolyzed (bright red) or lipaemic (milky) serum samples as well as samples containing sodium azide as preservative.

Store serum samples at +2...+8 °C for no more than 2 days. Aliquot and freeze the samples for longer storage (≤ -20 °C). Avoid repeated freezing.

6. EXPECTED VALUES

Serum samples collected between 9 and 11 a.m. from 200 apparently healthy people, both males and females, between the ages of 18–45, were assayed with **Free T4 kit**. Free T₄ concentration range in the samples was 10.0–23.2 pmol/L. These limits should be considered as guidelines only.

It is highly recommended each laboratory to determine its own reference range of free T₄ concentrations.

7. QUALITY CONTROL

It is recommended to use control samples according to the state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

8. REAGENT PREPARATION

Allow all the reagents to reach room temperature, and then thoroughly stir.

MP Keep the **microplate** at room temperature (+18...+25 °C) for at least 30 minutes before opening the bag. Place required number of strips onto strip holder. Place unused strips into the resealable zipper bag and reseal duly.

CAL CONTROL Calibrators and Control

Liquid calibrators and control are ready to use.

Prepare lyophilized calibrators and control as follows. Gently tap on the vial caps to knock off all the dry matter. Open the vials and carefully place the caps upside down on the clean dry surface. Add 0.5 mL of distilled or deionized water to each vial with lyophilized calibrators and control, close vials with the corresponding caps and leave for 10 min at room temperature without stirring. Then stir gently avoiding foaming, until the dry matter is completely dissolved. Leave for another 10 minutes at room temperature stirring gently periodically. Make sure that no dry matter is left on the caps and walls of the vials.

WASH P Prepare required volume of **wash solution** by dilution of the concentrate 20-fold with distilled or deionized water. For example:

5 mL **WASH P** **20X** + 95 mL water.

Mix thoroughly, avoid foaming.

SUB Protect **substrate** from direct light.

9. SAMPLE PREPARATION

Allow samples to reach room temperature. Stir samples gently in order to ensure homogeneity.

10. PERFORMANCE CHARACTERISTICS OF THE ASSAY

10.1. Calibration-Traceability

Free T4 kit was calibrated via commercially available diagnostic kit, which had been standardized using equilibrium dialysis.

10.2. Specificity

Cross-reaction of anti-T₄ monoclonal antibodies with different thyroids is shown below.

Thyroid	Cross-reaction, %
L-Thyroxin	100.0
L-Triiodothyronine	8.0
L-Diiodotyrosine	0.001

10.3. Analytical Sensitivity

Analytical sensitivity of **Free T4 kit**, i.e. concentration, which can be distinguished from zero calibrator is 1.0 pmol/L. It was defined as mean OD of 10 replicates of calibrator 0 minus 2 SD.

10.4. Measurement Range

The **Free T4 kits** (manual and Alisei kit) were validated for measurement of free T₄ concentration within the concentration diapason of 1.0 – 100 pmol/L.

Measurement Units

In **Free T4 kit** the concentrations of calibrators are specified in pmol/L. To convert into ng/L, multiply the concentration in pmol/L by 0.777.

10.5. Intra- and Inter-Assay Variation (Precision)

For **intra-assay CV** determination, 8 serum samples were run, each in 9 replicates. The results are shown below.

Manual kit

Sample	Mean free T ₄ concentration, pmol/L	Intra-assay CV	
		SD	CV, %
1	11.1	0.22	2.0
2	12.8	0.51	4.0
3	12	0.46	3.8
4	18.1	0.39	2.1
5	18	0.19	1.0
6	26.2	1.15	4.4
7	50.2	3.71	7.4
8	95.0	7.70	7.3

Alisei kit

Sample	Mean concentration of free T ₄ , pmol/L	Intra-assay CV	
		SD	CV, %
1	10.2	0.65	6.4
2	15.4	0.54	3.8
3	23.5	1.08	4.6
4	33.4	1.46	4.4
5	43.2	1.86	4.3
6	53.8	3.53	6.6
7	75.5	2.05	2.7
8	91.4	3.68	4.0

For **inter-assay CV** determination, 8 serum samples were assayed 3 times by different operators with 1-week interval. Each specimen was run in 9 replicates. The results are shown below.

Manual kit

Sample	Mean concentration of free T ₄ , pmol/L			Inter-assay precision	
	1 assay	2 assay	3 assay	SD	CV, %
1	11.1	11.1	10.2	0.52	4.8
2	12.8	12.4	11.5	0.66	5.4
3	12	12.5	11.2	0.66	5.5
4	18.1	18.3	16.9	0.76	4.3
5	18	17.8	17.1	0.48	2.7
6	26.2	25.1	24.2	1.01	4
7	50.2	50.2	44.8	3.10	6.4
8	95	96.9	88.1	6.22	6.7

Alisei kit

Sample	Mean concentration of free T ₄ , pmol/L			Inter-assay precision	
	1 assay	2 assay	3 assay	SD	CV, %
1	10.2	11.6	10.1	0.84	7.9
2	15.4	16.7	15.9	0.66	4.1
3	23.5	23.0	23.9	0.45	1.9
4	33.4	30.9	32.8	1.30	4.0
5	43.2	44.5	49.8	3.50	7.6
6	53.8	53.8	55.1	0.75	1.4
7	75.5	80.4	78.3	2.46	3.1
8	91.4	80.5	83.6	5.62	6.6

11. LIMITATION OF THE METHOD

Any clinical diagnosis should not be based on the results of in vitro diagnostic methods alone. For diagnosis establishment, a physician is supposed to consider all available clinical and laboratory findings.

12. SAFETY PRECAUTIONS

- **This kit is for in vitro diagnostic use only.** Operator should follow the manual closely in order to ensure reliable data. The manual is valid for the present kit only, within the listed composition. Any substitution of kit components is not allowed by CE regulations.
- Do not use kit or its components after the expiration date indicated on the label. Take into account stability period for reconstituted reagents.
- Do not mix or use together reagents from different lots of the kit except substrate, stop solution and wash solution.
- Do not use substrate, stop solution and wash solution supplied by other vendors.
- Use only “P”-labeled Astra Biotech wash solution.
- Note that stop solution is 1 M HCl solution. Avoid contacts with skin and mucosa. In case of contact rinse affected area thoroughly with plenty of water and seek medical advice.
- Take into account the following common procedural notes:
 - always pipette reagents into wells immediately after washing procedure;
 - avoid contamination of the solutions;

- in case of partial use of the kit, dispense only required volume of the reagent into the tray;
- do not pour unused reagents back into the original vials;
- avoid exposure to direct sunlight during incubations;
- always pipette reagents in the same order to minimize reaction time differences between wells; the total dispensing time for the calibrators, control and samples must not exceed 15 min;
- the incubation temperature for all the immunological reactions must be kept at +37 °C;
- do not touch the bottom of the wells;
- calibrators should be measured in each separate assay; It is also recommended to measure each time free T₄ concentration in the control.
- TMB solution should be colourless. Light colouring of solution is admissible. Avoid direct exposure of substrate to sunlight.
-  Source materials of human origin used for kit components preparation were tested and found negative for HBsAg, anti-HIV and anti-HCV antibodies. However, none of known laboratory test guarantees absence of these viral agents. Therefore, all kit components and patient's samples should be handled as potentially hazardous.
-  After usage strips, calibrators, control, specimens and all consumables which contacted with specimens during handling, storage or assay (tubes, vials, gloves, pipette tips etc.) should be collected separately and sterilized by autoclaving. Instead of autoclaving pipette tips may be sterilized by disinfectant treatment. After sterilization all

components and expendable materials may be utilized as non-dangerous garbage. Other components of the kit should be discarded into conventional garbage.

-  During manual washing procedure do not discard the contents of the wells directly to drainage. Use a container with disinfectant solution.
-  As the kit contains potentially hazardous material, the following precautions should be taken:
 - do not smoke, eat or drink while performing the assay;
 - always use protective gloves;
 - never pipette material by mouth;
 - in case of spilling, wipe up the spills promptly and wash affected area thoroughly using decontaminant.
- GLP including all general and individual regulations should be applied for the kit usage.

As the kit contains irritant (**CONJ**, **CAL**, **CONTROL**), the following precautions should be observed:

- P261 - Avoid breathing spray;
- P272 - Contaminated work clothing should not be allowed out of the workplace;
- P280 - Wear protective gloves/protective clothing/eye protection;
- P302+P352 - IF ON SKIN: Wash with plenty of soap and water;
- P333+P313 - If skin irritation or rash occurs: Get medical advice/attention;
- P363 - Wash contaminated clothing before reuse;

- P501 - Dispose of contents/container in accordance with national regulation.

Precautionary statements according to Regulation EC № 1272/2008.

MANUAL TEST (REF 24-03)**13.1. Material Provided**

MP	Microplate: 12 breakable 8-well strips (total 96 wells) coated with anti-T ₄ monoclonal antibodies	1 pcs
CONJ	Conjugate: solution containing T ₄ conjugated with HRP	18 mL, ready to use
0-5 CAL	Free T₄ calibrators: protein-based solutions or lyophilized preparations containing known free T ₄ concentrations – 0; 3; 6; 12; 30; 100 pmol/L. For exact free T ₄ concentrations see vial labels	6 vials, 0.5 mL each; ready to use or lyophilized
CONTROL	Free T₄ control: protein-based solution or lyophilized preparation containing known free T ₄ concentration. For exact range of free T ₄ concentration see vial label.	0.5 mL, ready to use or lyophilized
SUB	Substrate (TMB solution): 3,3',5,5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide	14 mL, ready to use
WASH P 20X	Wash solution P, 20X concentrated: surfactant in buffered saline, sufficient for preparation of 280 mL of solution	14 mL, concentrated
STOP	Stop solution: 1 M HCl solution	14 mL, ready to use

13.2. Equipment and Materials Required but not provided

- 1-channel calibrated variable precision pipettes, with disposable tips;
- 8-channel calibrated variable precision pipette, with disposable tips;
- microplate incubator/shaker (+37 °C, 500–800 rpm);
- manual or automatic equipment for rinsing wells;
- calibrated microplate reader (450 nm);
- vortex tube mixer;
- deionized or distilled water;
- graduated beaker and cylinder of appropriate volume;
- latex or plastic gloves;
- trays for pipetting reagents with 8-channel pipette;
- disinfectant;
- absorbent material (for manual wash).

13.3. Test Procedure

The **Free T4 kit** is designed for 96 tests. This is sufficient for the **quantitative assay** of 40 unknowns, 6 calibrators, 1 control and 1 blank (OD of TMB solution) in duplicates, provided that all the strips are used simultaneously.

13.3.1. Assay procedure

(Assay scheme, section 13.5)

- A.** Pipette **20 µL** calibrators **CAL** (0-5), control **CONTROL** and **patient's samples in duplicates** into the respective wells; **Leave wells A1-A2 empty for blank!**
- B.** Pipette **150 µL** of conjugate **CONJ** into each well, **except wells A1-A2.**
- C.** Incubate strips for **60 minutes while shaking (500–800 rpm) at +37 °C.**
- D.** Wash 4 times, as described below.
- E.** Pipette **100 µL** substrate **SUB** into each well (including blank); incubate **at room temperature in the dark for 15-30 minutes** (depending on the colour intensity) or **10 minutes while shaking (500-800 rpm) at +37 °C.**
- F.** Pipette **100 µL** stop solution **STOP** into each well (including blank) in the same sequence and at the same speed as used for dispensing TMB substrate. **Shake for 1–2 min at room temperature.**
- G.** Read **OD at 450 nm within 20 min.**

13.3.2. Wash Procedure

It is advisable to use an automatic microplate washer set at 4 wash cycles and a volume of 300 µL of wash solution per well per cycle.

If an automatic washer is not available, the wash procedure can be carried out manually as follows:

- remove the contents of the wells into a container with disinfectant;
- dispense 300 μ L of wash solution (section 8) into each well, shake the plate carefully for 5–10 sec and remove the contents of the wells; repeat 4 times;
- strike the wells sharply on absorbent material to remove any liquid residue.

13.4. Data Processing

If the reader cannot be adjusted to zero using the substrate blank in wells A1-A2, subtract mean OD value of wells A1-A2 from all OD values before further calculations.

Example:

OD (Cal 0) measured = 2.28 and OD (blank) = 0.06;

OD (Cal 0) calculated = $2.28 - 0.06 = \underline{2.22}$

13.4.1. Data Reliability (OD measured at 450 nm)

The data should meet the following criteria:

- average blank OD (in wells A1-A2) ≤ 0.100 ;
- average OD of Cal 0 ≥ 1.0 (after blank subtraction);
- control's concentration must fall within the acceptability range that is shown on the vial label.

If the data obtained do not meet the criteria, the results are considered unreliable and the test should be repeated.

13.4.2. Quantitative Determination

Specialized software for quantitative determination is recommended. Mean OD of calibrators are plotted versus their respective free T_4 concentrations using **4PL or 5PL fit** (Fig. 2 Typical standard curve). Calculate concentration of free T_4 in samples using standard curve.

Any extrapolation of the standard curve to free T_4 concentration above the nominal value of the calibrator 5 (approximately 100 pmol/L) is forbidden.


Fig. 2 Example of typical standard curve
Do not use for evaluation of real assay data!

13.5. Assay scheme

Reagents \ Wells	«Blank»	CAL CONTROL	Samples
CAL CONTROL	–	20 µL	–
Samples	–	–	20 µL
CONJ	–	150 µL	150 µL
Incubation No.1	60 min, +37 °C, 500–800 rpm		
WASH P (diluted)	4 x 300 µL		
SUB	100 µL	100 µL	100 µL
Incubation No.2	15-30 min, +18...+25 °C, in the dark		
	10 min, +37 °C, 500–800 rpm		
STOP	100 µL	100 µL	100 µL
Stirring	1–2 min, +18...+25 °C		
OD measuring	450 nm		
Calculations	Corresponding software		

14. AUTOMATIC TEST (REF 24-03A)

14.1. Material Provided

MP	Microplate: 12 breakable 8-well strips (total 96 wells) coated with anti-T ₄ monoclonal antibodies	1 pcs
CONJ	Conjugate: solution containing T ₄ conjugated with HRP	18 mL, ready to use
0-5 CAL	Free T₄ calibrators: protein-based solutions or lyophilized preparations containing known free T ₄ concentrations- 0; 3; 6; 12; 30; 100 pmol/L. For exact free T ₄ concentrations see vial labels	8 vials, 0.5 mL each: CAL 2, 3 – 2 × 0.5 mL; CAL 0, 1, 4, 5 – 1 × 0.5 mL; ready to use or lyophilized
CONTROL	Free T₄ control: protein-based solution or lyophilized preparation containing known free T ₄ concentration. For exact range of free T ₄ concentration see vial label	2 × 0.5 mL; ready to use or lyophilized
SUB	Substrate (TMB solution): 3,3',5,5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide	14 mL, ready to use
TRIAL 5000X	Trial solution, 5000X concentrated: solution of detergent	It is delivered by separate order
STOP	Stop solution: 1 M HCl solution	19 mL, ready to use
WASH P 20X	Wash solution P, 20X concentrated: surfactant in buffered saline, sufficient for preparation of 1000 mL of solution	50 mL, concentrated

Note: Extra vials of **CAL** 2, 3 and **CONTROL** are provided for recalibration on reference calibration curve. For detailed information see instructions for use for analyzer "Alisei".

14.2. Equipment and Materials Required but not provided

- Automatic analyzer for ELISA kits on microplate (analyzer “Alisei”);
- multipurpose polypropylene tubes 12x75 volume 5.5 mL;
- **TRIAL 5000X** Trial solution, 5000X concentrated*;
- deionized or distilled water;
- latex or plastic gloves.

14.3. Test Procedure

Free T4 kit is designed for 96 tests. This is sufficient for 88 unknowns, 6 calibrators, 1 control and 1 blank (OD of TMB solution) in monoplicates, provided that all the strips are used simultaneously.

14.3.1. Reagent Preparation

Allow all the reagents to reach room temperature, and then thoroughly stir.

Note: Preparation of other reagents see in section 8.

TRIAL Working solution of Trial for washing internal hydraulic lines and needles of automatic analyzer. Prepare required volume of **Trial solution** by dilution of the concentrate 5000-fold with distilled or deionized water before analysis:

2 mL **TRIAL 5000X** + 9998 mL water

Mix thoroughly, avoid foaming.

* Reagent is not included in the kit, it is delivered by separate order.
Astra Biotech

14.3.2. Assay Procedure for Automatic Test

While using for the procedure analyzer “Alisei”, refer to its relative manual.

The analysis on EIA analyzer “Alisei” is completely automatic: pipetting of reagents, washing, incubation, OD measurement, analysis of results. The program of calculation free T4 concentration is hardwired in analyzer memory.

14.4. Data Processing

14.4.1. Data Reliability (OD measured at 450 nm)

See criteria in section 13.4.1.

14.5. Safety precautions

- If kit for automatic analyzer is used in several separate experiments it is necessary to take reagents from automatic analyzer immediately after pipetting them in the wells of all plates because liquid evaporates from vials. Put the reagents in refrigerator;
- Recalibration using calibration curve, obtained with a kit of any other lot, is not permitted.

April, 12, 2018

Astra Biotech GmbH
Rudower Chaussee 29,
12489 Berlin, Germany
Telefon: +49 (0)30 74696509
E-Mail: info@astrabiotech.de
www.astrabiotech.de