# Prolactin kit

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## 2. INTENDED USE

**The Prolactin kit is provided for the quantitative determination of prolactin in human serum.**

This test has 2 complete sets:

**REF 21-04** for manual use,  
**REF 21-04 A** must be used with ELISA automatic instrument "Alisei" manufactured by NEXT Level S.r.l. hereinafter referred to as analyser "Alisei".

Instructions for use:

Section 13 for manual test,  
Section 14 for *analyser* "Alisei".

**Note 1:** *Take into account that calibrators' nominals can be different for manual and automatic test kits.*

**Note 2:** *We guarantee applications of test*

**REF 21-04 A** *only on analyser "Alisei",*  
**REF 21-04** *only for manual use.*

*While using non predefined methods of use, it is under end user responsibility, to make sure that it was appropriately tested for ELISA kits.*


Prolactin is a single-chain protein hormone with a molecular mass of about 23,000 Da. It is synthesized in the frontal lobe of pituitary gland. The determination of prolactin is a valuable tool in diagnostics of testicular and ovarian dysfunctions. Hyperprolactinemia in women may lead to galactorrhea, amenorrhea and other distortions in menstrual cycle. In men it may result in impotence or loss of libido. Measurement of circulating prolactin is used as a primary test for barrenness.

Pathological hyperprolactinemia takes place in the case of hypothyroidism, renal insufficiency and in patients with a pituitary gland tumor - prolactinoma. Physiological increase of prolactin level takes place during the gestation, lactation, in sleep and after physical and emotional stress.

### 3. PRINCIPLE OF TEST

The **Prolactin kit** is a “sandwich” type of solid-phase enzyme immunoassay, based on two monoclonal antibodies that are specific for different epitopes of prolactin molecule. One of these antibodies is conjugated with horseradish peroxidase (HRP); the other is coated onto the inner surface of microwells. Prolactin molecules from the serum sample bind to both immobilized antibody and anti-prolactin-peroxidase conjugate. Then the wells are washed with wash solution to remove any material not bound on the inner surface of the wells. Quantity of the bound conjugate is directly proportional to prolactin level in sample.

During the incubation with TMB solution the colour is developing. The intensity of the colour is directly proportional to prolactin concentration in specimens. Prolactin concentration in the patient sample is read from a standard curve that is processed in each assay.

*Fig. 1 Assay scheme*

#### **4. STORAGE AND STABILITY OF THE KIT**

The expiration date of the kit is printed on the box label; the expiration date for each component is printed on the respective label.

The **Prolactin kit** should be stored at +2...+8 °C upon receipt, preferably in the original kit box, until the expiration date. Storage at +25 °C is allowed but for no more than 5 days.

Shelf life of the kit is 18 months.

After initial opening, the kit is stable for 12 month if stored at +2...+8 °C.

If used in several separate experiments, after initial opening kit contents should be stored as follows but never used longer than the expiration date:

- unused strips: in a firmly closed resealable zipper bag at +2...+8 °C until the expiration date;
- vials with conjugate, calibrators and control (ready-to-use): for 12 month;
- vials with calibrators and control (reconstituted): at +2...+8 °C for no more than 6 months after opening;

- vial with substrate: at +2...+8 °C until the expiration date, protected from light;
- vials with concentrated Trial, concentrated wash solution and stop solution: at +2...+8 °C until the expiration date;
- wash solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days or at +2...+8°C for 4 weeks in a firmly closed bottle;
- Trial solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days, in a firmly closed bottle.

### **Damaged Test Kits**

In case of any severe damage of the test kit or components, Astra Biotech GmbH has to be informed in writing, during one week after receiving the kit. Usage of severely damaged single components for a test run is not recommended.

## **5. SAMPLE COLLECTION AND STORAGE**

Collect blood by venipuncture. After clotting, the serum is separated by centrifugation.

Do not use plasma, haemolyzed (bright red) or lipaemic (milky) serum samples as well as samples containing sodium azide as preservative.

Store serum samples at +2...+8 °C for no more than 5 days. Aliquot and freeze the samples for a longer storage ( $\leq -20$  °C). Avoid repeated freezing.

## 6. EXPECTED VALUES

Serum samples of 120 apparently healthy people were assayed with **Prolactin** kit. The results are shown in the table. These limits should be considered as guidelines only.

Group	n	Prolactin concentration range (mIU/L)
Female	120	67-726
Male		105-540

It is highly recommended for each laboratory to determine its own reference range of prolactin concentrations.

## 7. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

## 8. REAGENT PREPARATION

Allow all the reagents to reach room temperature, and then thoroughly stir.

### **MP** Microplate

Keep the **microplate** at room temperature (+18...+25 °C) for at least 30 minutes before opening the bag. Place required number of strips onto strip holder. Place unused strips into the resealable zipper bag and reseal duly.

### **CAL** **CONTROL** Calibrators and Control

Liquid calibrators and control are ready to use.

Prepare lyophilized calibrators and control as follows. Gently tap on the vial caps to knock off all the dry matter. Open the vials and carefully place the caps upside down on the clean dry

surface. Add 0.5 mL of distilled or deionized water to each vial with lyophilized calibrators and control, close vials with the corresponding caps and leave for 10 min at room temperature without stirring. Then stir gently avoiding foaming, until the dry matter is completely dissolved. Leave for another 10 minutes at room temperature stirring gently periodically. Make sure that no dry matter is left on the caps and walls of the vials.

### **WASH P** Wash solution

Prepare required volume of **wash solution** by dilution of the concentrate 20-fold with distilled or deionized water. For example:

5 mL **WASH P** **20X** + 95 mL water.

Mix thoroughly, avoid foaming.

### **SUB** Substrate

Protect **substrate** from direct light.

## **9. SAMPLE PREPARATION**

Allow samples to reach room temperature. Stir samples gently in order to ensure homogeneity.

## **10. PERFORMANCE CHARACTERISTICS OF THE ASSAY**

### **10.1 Calibration - Traceability**

The **Prolactin kit** was calibrated against the 3<sup>rd</sup> WHO International Standard 84/500.

## 10.2 Specificity

No cross-interaction of both anti-prolactin monoclonal antibodies with LH, FSH, TSH, placental lactogen, and hGH was detected.

## 10.3 Analytical Sensitivity

Analytical sensitivity of **Prolactin kit** i.e. concentration that can be distinguished from zero calibrator, is 50 mIU/L. It is defined as mean OD of 10 replicates of calibrator 0 plus two standard deviations.

## 10.4 Measurement Range

The Prolactin kits (manual and Alisei kit) were validated for measurement of prolactin concentration within the concentration diapason of 50 – 4500 mIU/L.

## 10.5 Hook Effect

The **Prolactin kit** shows no high dose hook effect for prolactin concentrations up to 100,000 mIU/L. High dose hook effect was determined by spiking calibrator 0 matrix with antigen.

## 10.6 Intra- and Inter-Assay Variation

For **Intra-Assay CV** determination, 7 serum samples (manual kit) and 8 serum samples (Alisei kit) were assayed in 9 replicates each. The results are shown below.

## Manual kit

Sample	Mean prolactin concentration, mIU/L	Intra-Assay CV	
		SD	CV, %
1	155	21.7	14.0
2	432	26.8	6.2
3	513	29.2	5.7
4	660	33.7	5.1
5	1234	55.5	4.5
6	1546	94.3	6.1
7	2567	179.7	7.0

## Alisei kit

Sample	Mean prolactin concentration, mIU/L	Intra-Assay CV	
		SD	CV, %
1	125	9.9	7.9
2	362	15.9	4.4
3	498	23.0	4.6
4	772	17.9	2.3
5	1240	76.8	6.2
6	1718	92.9	5.4
7	2004	66.0	3.3
8	2757	112.7	4.1

For **Inter-Assay CV** determination, 8 serum samples were assayed 3 times by different operators with 1-week interval. Each specimen was assayed in 9 replicates. The results are shown below.

## Manual kit

Sample	Mean prolactin concentration, mIU/L			Inter-Assay CV	
	1 assay	2 assay	3 assay	SD	CV, %
1	154	139	144	7.6	5.4
2	218	202	195	11.8	5.9
3	302	278	262	20.1	7.5
4	324	301	281	21.5	7.4
5	348	311	308	22.3	7.2
6	470	432	431	22.2	5.2
7	958	840	851	65.2	7.7
8	2410	2409	2049	208.1	9.3

## Alisei kit

Sample	Mean prolactin concentration, mIU/L			Inter-Assay precision	
	1 assay	2 assay	3 assay	SD	CV, %
1	125	145	140	10.4	7.6
2	362	379	344	17.5	4.8
3	498	439	468	29.5	6.3
4	772	861	769	52.3	6.5
5	1240	1334	1272	47.8	3.7
6	1718	1900	1803	91.1	5.0
7	2004	2175	2184	101.4	4.8
8	2757	2513	2499	145.1	5.6

## 11. LIMITATION OF THE METHOD

Any clinical diagnosis should not be based on the results of in vitro diagnostic method alone. For diagnosis establishment, a physician is supposed to consider all available clinical and laboratory findings.

## 12. SAFETY PRECAUTIONS

- **This kit is for in vitro diagnostic use only.** Operator should follow the manual closely in order to ensure reliable data. The manual is valid for the present kit only, within the listed composition. Any substitution of kit components is not allowed by CE regulations.
- Do not use kit or its components after the expiration date indicated on the label. Take into account stability period for reconstituted reagents.
- Do not mix or use together reagents from different lots of the kit except substrate, stop solution and wash solution.
- Do not use substrate, stop solution and wash solution supplied by other vendors.
- Use only “P”-labeled Astra Biotech wash solution.
- Note that stop solution is 1 N HCl solution. Avoid contacts with skin and mucosa. In case of contact rinse affected area thoroughly with plenty of water and seek medical advice.
- Take into account the following common procedural notes:
  - always pipette reagents into wells immediately after washing procedure;
  - avoid contamination of the solutions;
  - in case of partial use of the kit, dispense only required volume of the reagent into the tray;
  - do not pour unused reagents back into the original vials;
  - avoid exposure to direct sunlight during incubations;
  - always pipette reagents in the same order to minimize reaction time differences between wells; the total dispensing time for the calibrators, control and samples must not exceed 15 min;

- the incubation temperature for all the immunological reactions must be kept at +37 °C;
- do not touch the bottom of the wells.
- calibrators should be measured in each separate assay. It is also recommended to measure each time prolactin concentration in the control.
- remove the adhesive foil carefully to avoid contamination and don't use the adhesive foil repeatedly.
- TMB solution should be colourless. Light colouring of solution is admissible. Avoid direct exposure of substrate to sunlight.
-  Source materials of human origin used for kit components preparation were tested and found negative for HBsAg, anti-HIV and anti-HCV antibodies. However, none of known laboratory test guarantees absence of these viral agents. Therefore, all kit components and patient's samples should be handled as potentially hazardous.
-  After usage strips, calibrators, control, specimens and all consumables which contacted with specimens during handling, storage or assay (tubes, vials, gloves, pipette tips etc.) should be collected separately and sterilized by autoclaving. Instead of autoclaving pipette tips may be sterilized by disinfectant treatment. After sterilization all components and expendable materials may be utilized as non-dangerous garbage. Other components of the kit should be discarded into conventional garbage.
-  During manual washing procedure do not discard the contents of the wells directly to drainage. Use a container with disinfectant solution.

-  As the kit contains potentially hazardous material, the following precautions should be taken:
  - do not smoke, eat or drink while performing the assay;
  - always use protective gloves;
  - never pipette material by mouth;
  - in case of spilling, wipe up the spills promptly and wash affected area thoroughly using decontaminant.
- GLP including all general and individual regulations should be applied for the kit usage.

As the kit contains irritant (**CONJ**, **CAL**, **CONTROL**), the following precautions should be observed:

- P261 - Avoid breathing spray;
- P272 - Contaminated work clothing should not be allowed out of the workplace;
- P280 - Wear protective gloves/protective clothing/eye protection;
- P302+P352 - IF ON SKIN: Wash with plenty of soap and water;
- P333+P313 - If skin irritation or rash occurs: Get medical advice/attention;
- P363 - Wash contaminated clothing before reuse;
- P501 - Dispose of contents/container in accordance with national regulation.

Precautionary statements according to Regulation EC № 1272/2008.

### 13. MANUAL TEST (REF 21-04)

#### 13.1. Material Provided

<b>MP</b>	<b>Microplate:</b> 12 breakable 8-well strips (total 96 wells), coated with anti-prolactin monoclonal antibodies	1 pcs
<b>CONJ</b>	<b>Conjugate:</b> solution contains anti-prolactin monoclonal antibodies conjugated with HRP	14 mL, ready to use
<b>0-5</b> <b>CAL</b>	<b>Prolactin calibrators:</b> protein-based solution or lyophilized preparations containing known prolactin concentrations – 0; 100; 500; 1200; 2500; 4500 mIU/L. (approximate values - do not use for evaluation of real assay data). The concentrations of calibrators may be different for Assay Protocol 1 and 2. For lot specific Prolactin concentrations see values in the Quality Control Sheet.	6 vials, 0.5 mL each, ready to use or lyophilized
<b>CONTROL</b>	<b>Prolactin control:</b> protein-based solution or lyophilized preparation containing known prolactin concentration. The range of concentration may be different for Assay Protocol 1 and 2. For lot specific ranges of Prolactin concentrations see values in the Quality Control Sheet.	0.5 mL, ready to use or lyophilized
<b>SUB</b>	<b>Substrate (TMB solution):</b> 3,3',5,5'-tetramethylbenzidine solution in citrate solution containing hydrogen peroxide	14 mL, ready to use
<b>WASH P</b>	<b>Wash solution P, 20X concentrated:</b> surfactant in buffered	2x14 mL, concentrated

<b>20X</b>	saline, sufficient for preparation of 560 mL of solution	
<b>STOP</b>	<b>Stop solution:</b> 1 N HCl solution	14 mL, ready to use
	<b>Adhesive foil</b>	2x1 foil (optional)

### 13.2. Equipment and Materials Required but not provided

- 1-channel calibrated variable precision pipettes, with disposable tips;
- 8-channel calibrated variable precision pipette, with disposable tips;
- microplate incubator (+37 °C) or microplate incubator-shaker (+37 °C, 500–800 rpm);
- manual or automatic equipment for rinsing wells;
- calibrated microplate reader (450 nm);
- vortex tube mixer;
- deionized or distilled water;
- graduated beaker and cylinder of appropriate volume;
- latex or plastic gloves;
- trays for pipetting reagents with 8-channel pipette;
- disinfectant;
- absorbent material (for manual wash).

### 13.3 Test Procedure

The **Prolactin kit** is designed for 96 tests. This is sufficient for 40 unknowns, 6 calibrators, 1 control and 1 blank (OD of TMB solution) in duplicates, provided that all the strips are used simultaneously.

#### 13.3.1 Assay Procedure

**Note:** *Please check in the Quality Control Sheet supplied with the kit for which Assay Protocol your kit has been validated. If your kit has been validated for Assay Protocol 1, please follow the Assay Protocol 1 described below. If your kit has been additionally validated for Assay Protocol 2, please consider the Supplement 1 to the Instructions for use provided with the kit.*

#### Assay Protocol 1

(see also Assay scheme to Assay Protocol 1, section 13.5)

**Note:** *Consider concentrations of calibrators and control range for Assay Protocol 1 provided in the Quality Control Sheet.*

**All samples should be tested in duplicates.**

- A.** Pipette **100  $\mu$ L conjugate** **CONJ** into each well, **except wells A1-A2 (blank)**;
- B.** Pipette **20  $\mu$ L calibrators** **CAL**, control **CONTROL** and prepared patient's samples in duplicates, into the respective wells.

**Leave wells A1-A2 empty for blank!**

**Note:** *total time of dispensing must not exceed 15 minutes, otherwise the test result may be unreliable, because the time of incubation will substantially vary for different samples.*

**Note:** *If an adhesive foil is used, put it now onto the plate to seal the cavities.*

- C. Incubate for **60 minutes at +37 °C while shaking (500–800 rpm)**.
- D. **Note:** *If an adhesive foil is used, remove it now from the plate.* Wash 5 times, as described in section 13.3.2.
- E. Pipette **100 µL substrate** **SUB** into each well (including blank); incubate strips **at room temperature (+18...+25 C) in the dark for 15-30 minutes**, depending on the colour intensity, or **10 minutes while shaking (500-800 rpm) at +37 °C**.
- F. Pipette **100 µL stop solution** **STOP** into each well (including blank) in the same sequence and at the same speed as used for dispensing TMB substrate. Shake for 1–2 min at room temperature.
- G. Read **OD at 450 nm within 20 min**.

### 13.3.2 Wash Procedure

It is advisable to use an automatic microplate washer set at 5 wash cycles and a volume of 300 µL of wash solution per well per cycle.

If an automatic washer is not available, the wash procedure can be carried out manually as follows:

- remove the contents of the wells into a container with disinfectant;

- dispense 300 µL of wash solution (section 8) into each well, shake the plate carefully for 5–10 sec and remove the contents of the wells; repeat 5 times;
- strike the wells sharply on absorbent material to remove any liquid residue.

### 13.4 Data Processing

If the reader cannot be adjusted to zero using the substrate blank in wells A1-A2, subtract mean OD value of wells A1-A2 from all OD values before further calculations.

*Example:*

OD (Cal 5) measured = 2.28 and OD (blank) = 0.06;

OD (Cal 5) calculated =  $2.28 - 0.06 = \underline{2.22}$

#### 13.4.1 Data Reliability (OD<sub>450 nm</sub>)

The data should meet the following criteria:


- average blank OD (in wells A1-A2)  $\leq 0.100$ ;
- average OD of Cal 5  $\geq 1.5$  (after blank subtraction);
  - control's concentration must fall within the acceptability range that is shown on the Quality Control Sheet.

If the data obtained do not meet the criteria, the results are considered unreliable and the test should be repeated.

### 13.4.2 Quantitative Determination

Specialized software for quantitative determination is recommended. Mean OD values of the calibrators at 450 nm versus their respective prolactin concentrations using **4PL** or **5PL fit** (Fig. 2 Typical standard curve). The arithmetical mean OD of the two «Blank» wells is used for zero setting of device.

Extrapolation of the standard curve to prolactin concentration above the nominal value of calibrator 5 is not permitted.

*Fig. 2 Example of typical standard curve*  
**Do not use for evaluation of real assay data!**

### 13.5 Assay scheme to Assay Protocol 1

Reagents	Wells		
	«Blank»	<b>CAL</b> <b>CONTROL</b>	Samples
<b>CONJ</b>	–	100 µL	100 µL
<b>CAL</b> <b>CONTROL</b>	–	20 µL	–
Samples	–	–	20 µL
Incubation No. 1	60 min, +37 °C, 500–800 rpm		
<b>WASH P</b> (diluted)	5 x 300 µL		
<b>SUB</b>	100 µL	100 µL	100 µL
Incubation No. 2	15–30 min, +18...+25 °C, in the dark		
	10 min, +37 °C, 500–800 rpm		
<b>STOP</b>	100 µL	100 µL	100 µL
Stirring	1–2 min, +18...+25 °C		
OD measuring	450 nm		
Calculations	Corresponding software		

## 14. AUTOMATIC TEST (REF 21-04 A)

### 14.1 Material Provided

<b>MP</b>	<b>Microplate:</b> 12 breakable 8-well strips (total 96 wells) coated with anti-Prolactin monoclonal antibodies	1 pcs
<b>CONJ</b>	<b>Conjugate:</b> solution contains anti-Prolactin monoclonal antibodies conjugated with HRP	14 mL, ready to use
<b>0-5 CAL</b>	<b>Prolactin calibrators:</b> protein-based solution or lyophilized preparations containing known prolactin concentrations – 0; 100; 500; 1200; 2500; 4500 mIU/L. (approximate values - do not use for evaluation of real assay data). For lot specific Prolactin concentrations see values provided in the Quality Control Sheet.	<b>8 vials,</b> 0.5 mL each: <b>CAL</b> 1, 4 – 2 × 0.5 mL; <b>CAL</b> 0, 2, 3, 5 – 1 × 0.5 mL; ready to use or lyophilized
<b>CONTROL</b>	<b>Prolactin control:</b> protein-based solution or lyophilized preparation containing known prolactin concentration. For lot specific ranges of Prolactin concentrations see Quality Control Sheet	2 × 0.5 mL; ready to use or lyophilized
<b>WASH P 20X</b>	<b>Wash solution P, 20X concentrated:</b> surfactant in buffered saline, sufficient for preparation of 1000 mL of solution	50 mL, concentrated
<b>SUB</b>	<b>Substrate (TMB solution):</b> 3,3',5,5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide	14 mL, ready to use
<b>STOP</b>	<b>Stop solution:</b> 1 N HCl solution	19 mL, ready to use
<b>TRIAL 5000X</b>	<b>Trial solution, 5000X concentrated:</b> solution of detergent	It is delivered by separate order.

**Note:** Extra vials of **CAL** 1, 4 and **CONTROL** are provided for recalibration on reference calibration curve. For detailed information see instructions for use for analyser "Alisei".

## 14.2. Equipment and Materials Required but not provided

- Automatic analyser for ELISA kits on microplate (analyser “Alisei”);
- multipurpose polypropylene tubes 12x75 volume 5.5 mL;
- **TRIAL 5000X** Trial solution, 5000X concentrated<sup>2</sup>;
- deionized or distilled water;
- latex or plastic gloves.

## 14.3. Test Procedure

The **Prolactin kit** is designed for 96 tests. This is sufficient for 88 unknowns, 6 calibrators, 1 control and 1 blank (OD of TMB solution) in monoplicates, provided that all the strips are used simultaneously.

### 14.3.1 Reagent Preparation

Allow all the reagents to reach room temperature (+18...+25 °C), and then thoroughly stir.

**Note:** *Preparation of other reagents see in section 8.*

**TRIAL** Working solution of Trial for washing internal hydraulic lines and needles of automatic analyser “Alisei”. Prepare required volume of **Trial solution** by dilution of the concentrate 5000-fold with distilled or deionized water before analysis:

2 mL **TRIAL 5000X** + 9998 mL water

Mix thoroughly, avoid foaming.

---

<sup>2</sup> Reagent is not included in the kit, it is delivered by separate order.

### **14.3.2 Assay Procedure for Automatic Test**

While using for the procedure analyser “Alisei”, refer to its relative manual.

The analysis on analyser “Alisei” is completely automatic: pipetting of reagents, washing, incubation, OD measurement, analysis of results. The program of calculation Prolactin concentration is hardwired in analyser memory.

## **14.4 Data Processing**

### **14.4.1 Data Reliability (OD<sub>450 nm</sub>)**

See criteria in section **13.4.1**.

## **14.5 Safety precautions**

- If kit for automatic analyser is used in several separate experiments it is necessary to take reagents from automatic analyser immediately after pipetting them in the wells of all plates because liquid evaporates from vials. Put the reagents in refrigerator;
- Recalibration using calibration curve, obtained with a kit of any other lot, is not permitted.

August 23, 2019,

Astra Biotech GmbH  
Rudower Chaussee 29,  
12489 Berlin, Germany  
Telefon: +49 (0) 30 74696509  
E-Mail: [info@astrabiotech.de](mailto:info@astrabiotech.de)  
[www.astrabiotech.de](http://www.astrabiotech.de)