# Testosterone kit

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## 2. INTENDED USE

The **Testosterone kit** is provided for the **quantitative** determination of **testosterone** in human serum.

This test has 2 complete sets:

**REF 21-02** for manual use,

**REF 21-02 A** must be used with ELISA automatic instrument "Alisei" manufactured by NEXT Level S.r.l. hereinafter referred to as analyser "Alisei".

Instructions for use:

Section 13 for manual test,

Section 14 for analyser "Alisei".

**Note 1:** *Take into account that calibrators' nominals can be different for manual and automatic test kits.*

**Note 2:** *We guarantee applications of test:*

**REF 21-02 A** only on analyser "Alisei",

**REF 21-02** only for manual use.


*While using non predefined methods of use, it is under end user responsibility, to make sure that it was appropriately tested for ELISA kits.*

Testosterone is a steroid hormone with a molecular mass of 288.4 Da. Testosterone is synthesized mainly in testicles and, in sufficiently less extent, in ovaries and adrenal cortex. The determination of serum testosterone is a valuable tool for investigation of testis function and diagnostics of some adrenal, ovary and testicle tumors as well as female hirsutism.

### 3. PRINCIPLE OF THE TEST

The **Testosterone kit** is a competitive solid phase enzyme immunoassay. During the incubation testosterone of tested samples and horseradish peroxidase (HRP) labeled testosterone bind to the antibodies coated onto the inner surface of the microplate wells until balance between them occurs. Separation of free and bound to antibodies testosterone and conjugate testosterone-peroxidase occurs while extracting the contents of the wells. The amount of bound conjugate is inversely proportional to the quantity of testosterone in the sample (Fig. 1).

During the incubation with TMB solution the colour is developing. The intensity of the colour is inversely proportional to the testosterone concentration in specimens. The testosterone concentration in the patient sample is read from a standard curve that is processed in each assay.

*Fig. 1 Assay scheme*

## 4. STORAGE AND STABILITY OF THE KIT

Expiry date of the kit is printed on the box label; expiry date for each component is printed on the respective label.

The **Testosterone kit** should be stored at +2...+8 °C upon receipt, preferably in the original kit box, until the expiration date. Storage at +25 °C is allowed but for no more than 5 days.

Shelf life of the kit is 18 months.

After initial opening the kit is stable for 12 months if stored at +2...+8 °C.

If used in several separate experiments, after initial opening kit contents should be stored as follows but never used longer than the expiration date:

- unused strips: in a firmly closed resealable zipper bag at +2...+8 °C until expiration date;
- vials with calibrators and control and conjugate (liquid): at +2...+8 °C for 12 months; vials with calibrators and control (reconstituted): at +2...+8 °C for no more than 1 month after opening;
- vial with substrate: at +2...+8 °C until expiration date, protected from light;
- vials with concentrated Trial, concentrated wash solution, and stop solution: at +2...+8 °C until expiry date;
- wash solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days or at +2...+8°C for no more than 4 weeks in a firmly closed bottle;
- Trial solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days in a firmly closed bottle.

## Damaged Test Kits

In case of any severe damage of the test kit or components, Astra Biotech GmbH has to be informed in writing, during one week after receiving the kit. Usage of severely damaged single components for a test run is not recommended.

### 5. SAMPLE COLLECTION AND STORAGE

Collect blood by venipuncture in a tube without anticoagulants. Allow blood to clot. Centrifuge the specimens to separate serum from the blood corpuscles.

Do not use plasma, haemolyzed (bright red) or lipaemic (milky) serum samples as well as samples containing sodium azide as preservative.

Store serum samples at +2...+8 °C for no more than 2 days. Aliquot and freeze the samples for longer storage ( $\leq -20$  °C). Avoid repeated freezing.

### 6. EXPECTED VALUES

Serum samples collected between 9 and 11 a.m. from 70 apparently healthy people, both males and females, between the ages of 21–45, were assayed with **Testosterone kit**. The results are listed below. These limits should be considered as guidelines only.

Group	Number of samples	Mean (nmol/L)	Range (nmol/L)
Female	37	1.8	0.5-4.3
Male	33	19.2	12.1-38.3

It is highly recommended each laboratory to determine its own reference range of testosterone concentrations.

## 7. QUALITY CONTROL

It is recommended to use control samples according to the state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

## 8. REAGENT PREPARATION

Allow all the reagents to reach room temperature (+18...+25 °C), and then thoroughly stir.

### **MP** Microplate

Keep the **microplate** at room temperature for at least 30 minutes before opening the bag. Place required number of strips onto strip holder. Place unused strips into the resealable zipper bag and reseal duly.

### **CAL** **CONTROL** Calibrators and Control

Liquid calibrators and control are ready to use.

Prepare lyophilized calibrators and control as follows. Gently tap on the vial caps to knock off all the dry matter. Open the vials and carefully place the caps upside down on the clean dry surface. Add 0.5 mL of distilled or deionized water to each vial with lyophilized calibrators and control, close vials with the corresponding caps and leave for 10 min at room temperature without stirring. Then stir gently avoiding foaming, until the dry matter is completely dissolved. Leave for another 10 minutes at room temperature stirring gently periodically. Make sure that no dry matter is left on the caps and walls of the vials.

### **WASH P** Wash solution

Prepare required volume of **wash solution** by dilution of the concentrate 20-fold with distilled or deionized water. For example:

5 mL **WASH P 20X** + 95 mL water

Mix thoroughly, avoid foaming.

### **SUB** Substrate

Protect **substrate** from direct light.

## 9. SAMPLE PREPARATION

Allow samples to reach room temperature. Stir samples gently in order to ensure homogeneity.

## 10. PERFORMANCE CHARACTERISTICS OF THE ASSAY

### 10.1 Calibration-Traceability:

The **Testosterone kit** was calibrated against the Working Standard, which had been manufactured by gravimetric method based on weighing purified synthetic testosterone into analyte-free matrix.

### 10.2 Specificity:

Cross-reaction of anti-testosterone monoclonal antibodies with different steroids is shown below.

<b>Steroid</b>	<b>Cross-reaction, %</b>
Testosterone	100.0
5- $\alpha$ -Dihydrotestosterone	9.8
Androstendiol	2.0
Androstenolone	0.2
Androstandione	1.3
Dehydroepiandrosterone	0.2
Cortisol	0.04
Progesterone	0.1
Estrone	0.06
Estradiol	0.9
Estriol	0.02

### 10.3 Analytical Sensitivity:

Analytical sensitivity of **Testosterone kit**, i.e. concentration, that can be distinguished from zero calibrator is 0.2 nmol/L (for manual and Alisei kits). It was defined as mean OD of 10 replicates of calibrator 0 minus 2 SD.

### 10.4 Measurement Range:

The **Testosterone kit** was validated for measurement of testosterone concentration within the concentration diapason of 0.2 - 50 nmol/L (for manual test) and 0.2 - 100 nmol/L (for Alisei kit).

### 10.5 Measurement Units:

In **Testosterone kit** the concentrations of calibrators are specified in nmol/L. To convert into ng/mL, multiply the concentration in nmol/L by 0.288.

### 10.6 Intra- and Inter-Assay Variation (Precision):

For **intra-assay CV** determination, 8 serum samples were run, each in 9 replicates. The results are shown below.

Manual kit

Sample	Mean testosterone concentration, nmol/L	Intra-assay CV	
		SD	CV, %
1	1.14	0.177	15.4
2	1.37	0.209	15.3
3	4.68	0.295	6.3
4	13.0	0.75	5.8
5	16.9	0.42	2.5
6	17.2	0.60	3.5
7	19.6	0.59	3.0
8	26.2	0.74	2.8

## Alisei kit

Sample	Mean testosterone concentration, nmol/L	Intra-assay CV	
		SD	CV, %
1	1.07	0.050	4.7
2	2.83	0.100	3.5
3	5.61	0.183	3.3
4	7.47	0.229	3.1
5	11.34	0.436	3.8
6	15.70	0.557	3.5
7	19.11	0.488	2.6
8	29.56	0.726	2.5

For **inter-assay CV** determination, 10 serum samples were assayed 4 times by different operators with 1-week interval. Each specimen was run in 2 replicates. The results are shown below.

## Manual kit

Sample	Mean testosterone concentration, nmol/L				Inter-assay CV	
	Assay 1	Assay 2	Assay 3	Assay 4	SD	CV, %
1	0.71	0.76	0.79	0.69	0.046	6.2
2	1.34	1.44	0.97	1.07	0.221	18.4
3	1.34	1.33	1.50	1.29	0.093	6.8
4	2.07	1.93	1.93	1.95	0.067	3.4
5	2.33	2.60	2.14	2.25	0.196	8.4
6	12.3	14.2	14.1	12.9	0.92	6.9
7	12.7	14.7	14.1	13.9	0.84	6.0
8	13.8	15.1	17.7	15.6	1.60	10.3
9	14.3	17.0	16.4	15.3	1.20	7.6
10	16.4	18.4	16.6	16.5	0.96	5.7

## Alisei kit

Sample	Mean testosterone concentration, nmol/L				Inter-assay CV	
	Assay 1	Assay 2	Assay 3	Assay 4	SD	CV, %
1	1.12	1.09	1.13	0.99	0.064	5.9
2	2.96	2.54	2.77	2.61	0.187	6.9
3	5.79	5.21	5.68	5.53	0.252	4.5
4	7.85	7.16	7.49	7.01	0.373	5.1
5	11.49	11.06	12.62	11.98	0.670	5.7
6	17.64	16.90	15.31	15.85	1.045	6.4
7	19.00	17.62	20.44	19.87	1.227	6.4
8	23.51	21.19	22.82	23.99	1.223	5.3
9	26.58	27.89	26.54	28.12	0.840	3.1
10	35.06	31.85	33.41	32.79	1.350	4.1

## 11 LIMITATION OF THE METHOD

Any clinical diagnosis should not be based on the results of in vitro diagnostic methods alone. For diagnosis establishment, a physician is supposed to consider all available clinical and laboratory findings.

## 12 SAFETY PRECAUTIONS

- **This kit is for in vitro diagnostic use only.** Operator should follow the manual closely in order to ensure reliable data. The manual is valid for the present kit only, within the listed composition. Any substitution of kit components is not allowed by CE regulations.
- Do not use kit or its components after expiry date indicated on the label. Take into account stability period for reconstituted reagents.

- Do not mix or use together reagents from different lots of the kit except substrate, stop solution and wash solution.
- Do not use substrate, stop solution and wash solution supplied by other vendors.
- Use only “P”-labeled Astra Biotech wash solution.
- Note that stop solution is 1 N HCl solution. Avoid contacts with skin and mucosa. In case of contact rinse affected area thoroughly with plenty of water and seek medical advice.
- Take into account the following common procedural notes:
  - always pipette reagents into wells immediately after washing procedure;
  - avoid contamination of the solutions;
  - in case of partial use of the kit, dispense only required volume of the reagent into the tray;
  - do not pour unused reagents back into the original vials;
  - avoid exposure to direct sunlight during incubations;
  - always pipette reagents in the same order to minimize reaction time differences between wells; the total dispensing time for the calibrators, control and samples must not exceed 15 min;
  - the incubation temperature for all the immunological reactions must be kept at recommended by this instruction levels;
  - do not touch the bottom of the wells;
  - calibrators should be measured in each separate assay. It is also recommended to measure each time testosterone concentration in the control;
  - remove the adhesive foil carefully to avoid contamination and don't use the adhesive foil repeatedly.

- TMB solution should be colourless. Light colouring of solution is admissible. Avoid direct exposure of substrate to sunlight.
-  Source materials of human origin used for kit components preparation were tested and found negative for HBsAg, anti-HIV and anti-HCV antibodies. However, none of known laboratory test guarantees absence of these viral agents. Therefore, all kit components and patient's samples should be handled as potentially hazardous.
-  After usage strips, calibrators, control, specimens and all consumables which contacted with specimens during handling, storage or assay (tubes, vials, gloves, pipette tips etc.) should be collected separately and sterilized by autoclaving. Instead of autoclaving pipette tips may be sterilized by disinfectant treatment. After sterilization all components and expendable materials may be utilized as non-dangerous garbage. Other components of the kit should be discarded into conventional garbage.
-  During manual washing procedure do not discard the contents of the wells directly to drainage. Use a container with disinfectant solution.
-  As the kit contains potentially hazardous material, the following precautions should be taken:
  - do not smoke, eat or drink while performing the assay;
  - always use protective gloves;
  - never pipette material by mouth;

- in case of spilling, wipe up the spills promptly and wash affected area thoroughly using decontaminant.
- GLP including all general and individual regulations should be applied for the kit usage.

As the kit contains irritant (**CONJ**, **CAL**, **CONTROL**), the following precautions should be observed:

- P261 - Avoid breathing spray;
- P272 - Contaminated work clothing should not be allowed out of the workplace;
- P280 - Wear protective gloves/protective clothing/eye protection;
- P302+P352 - IF ON SKIN: Wash with plenty of soap and water;
- P333+P313 - If skin irritation or rash occurs: Get medical advice/attention;
- P363 - Wash contaminated clothing before reuse;
- P501 - Dispose of contents/container in accordance with national regulation.

Precautionary statements according to Regulation EC No 1272/2008.

## 13 MANUAL TEST (REF 21-02)

### 13.1 Material Provided

<b>MP</b>	<b>Microplate:</b> 12 breakable 8-well strips (total 96 wells) coated with anti-testosterone monoclonal antibodies	1 pcs
<b>CONJ</b>	<b>Conjugate:</b> solution containing testosterone conjugated with HRP	18 mL, ready to use
<b>0-5 CAL</b>	<b>Testosterone calibrators:</b> protein-based solutions or lyophilized preparations containing known testosterone concentrations – 0; 0.5; 1.5; 5; 15; 50 nmol/L (approximate values - do not use for evaluation of real assay data). For lot specific testosterone concentrations see values in the Quality Control Sheet.	<b>6 vials,</b> 0.5 mL each; ready to use or lyophilized
<b>CONTROL</b>	<b>Testosterone control:</b> protein-based solution or lyophilized preparation containing known testosterone concentration. For lot specific ranges of testosterone concentration see values in the Quality Control Sheet.	0.5 mL, ready to use or lyophilized
<b>SUB</b>	<b>Substrate (TMB solution):</b> 3,3',5,5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide	14 mL, ready to use
<b>WASH P 20X</b>	<b>Wash solution P, 20X concentrated:</b> surfactant in buffered saline, sufficient for preparation of 560 mL of solution	2x14 mL, concentrated
<b>STOP</b>	<b>Stop solution:</b> 1 N HCl solution	14 mL, ready to use
	<b>Adhesive foil</b>	2x1 foil (optional)

## 13.2 Equipment and Materials required but not provided

- 1-channel calibrated variable precision pipettes, with disposable tips;
- 8-channel calibrated variable precision pipette, with disposable tips;
- microplate incubator (+37 °C) or microplate incubator-shaker (+37 °C, shaking speed 500–800 rpm);
- manual or automatic equipment for rinsing wells;
- calibrated microplate reader (450 nm);
- vortex tube mixer;
- deionized or distilled water;
- graduated beaker and cylinder of appropriate volume;
- latex or plastic gloves;
- trays for pipetting reagents with 8-channel pipette;
- disinfectant;
- absorbent material (for manual wash).

## 13.3 Test Procedure

The **Testosterone kit** is designed for 96 tests. This is sufficient for the **quantitative assay** of 40 unknowns, 6 calibrators, 1 control and 1 blank (OD of TMB solution) in duplicates, provided that all the strips are used simultaneously.

### 13.3.1 Assay Procedure

#### Assay Protocol 1

(see also Assay scheme to Assay Protocol 1, section 13.5)

**All samples should be tested in duplicates.**

- A.** Pipette **50 µL** calibrators **CAL** (0-5), control **CONTROL** and prepared **patient's samples in duplicates** into the respective wells; **Leave wells A1-A2 empty for blank!**
- B.** Pipette **150 µL** conjugate **CONJ** into each well, **except wells A1-A2.**

**Note:** *total time of dispensing must not exceed 15 minutes, otherwise the test result may be unreliable, because the time of incubation will substantially vary for different samples.*

**Note:** *If an adhesive foil is used, put it now onto the plate to seal the cavities.*

- C.** Incubate strips for **90 minutes while shaking (500–800 rpm) at room temperature (+18...+25 °C).**

**Note:** *If an adhesive foil is used, remove it now from the plate.*

- D.** Wash 4 times, as described below.
- E.** Pipette **100 µL** substrate **SUB** into each well (including blank); incubate **at room temperature in the dark for 15-30 minutes** (depending on the colour intensity) **or 10 minutes while shaking (500-800 rpm) at +37 °C.**
- F.** Pipette **100 µL** stop solution **STOP** into each well (including blank) in the same sequence and at the same speed as used for dispensing TMB substrate. **Shake for 1–2 min at room temperature.**
- G.** Read **OD at 450 nm within 20 min.**

## Assay Protocol 2: (without shaking)

(see also Assay scheme to Assay Protocol 2, section 13.6)

**All samples should be tested in duplicates.**

- A. Pipette **50 µL** calibrators **CAL** (0-5), control **CONTROL** and **patient's samples in duplicates** into the respective wells; **Leave wells A1-A2 empty for blank!**
- B. Pipette **150 µL** conjugate **CONJ** into each well, **except wells A1-A2.**

**Note:** *total time of dispensing must not exceed 15 minutes, otherwise the test result may be unreliable, because the time of incubation will substantially vary for different samples.*

**Note:** *If an adhesive foil is used, put it now onto the plate to seal the cavities.*

- C. Incubate strips for **90 minutes at room temperature** (+18...+25 °C) (pre-shake for 1-2 minutes at room temperature).

**Note:** *If an adhesive foil is used, remove it now from the plate.*

- D. Wash 4 times, as described below.
- E. Pipette **100 µL** substrate **SUB** into each well (including blank); incubate **at room temperature in the dark for 15-30 minutes** (depending on the colour intensity).
- F. Pipette **100 µL** stop solution **STOP** into each well (including blank) in the same sequence and at the same speed as used for dispensing TMB substrate. **Shake for 1–2 min at room temperature.**
- G. Read **OD at 450 nm within 20 min.**

### 13.3.2. Wash Procedure

It is advisable to use an automatic microplate washer set at 4 wash cycles and a volume of 300  $\mu\text{L}$  of wash solution per well per cycle.

If an automatic washer is not available, the wash procedure can be carried out manually as follows:

- remove the contents of the wells into a container with disinfectant;
- dispense 300  $\mu\text{L}$  of wash solution (section 8) into each well, shake the plate carefully for 5–10 sec and remove the contents of the wells; repeat 4 times;
- strike the wells sharply on absorbent material to remove any liquid residue.

### 13.4 Data Processing

If the reader cannot be adjusted to zero using the substrate blank in wells A1-A2, subtract mean OD value of wells A1-A2 from all OD values before further calculations.

*Example:*

OD (Cal 0) measured = 2.28 and OD (blank) = 0.06;

OD (Cal 0) calculated =  $2.28 - 0.06 = \underline{2.22}$

#### 13.4.1 Data Reliability (for OD measured at 450 nm)

The data should meet the following criteria:


- average blank OD (in wells A1-A2)  $\leq 0.100$ ;
- average OD of Cal 0  $\geq 1.0$  (after blank subtraction);
- control's concentration must fall within the acceptability range that is shown in Quality Control Sheet.

If the data obtained do not meet the criteria, the results are considered unreliable and the test should be repeated.

### 13.4.2 Quantitative Determination

Specialized software for quantitative determination is recommended. Mean OD of calibrators are plotted versus their respective testosterone concentrations using **4PL** or **5PL** fit (Fig. 2 Typical standard curve). Calculate concentration of testosterone in samples using standard curve.

Any extrapolation of the standard curve to testosterone concentration above the nominal value of the calibrator 5 (approximately 50 nmol/L) is forbidden.

*Fig. 2 Example of typical standard curve*  
**Do not use for evaluation of real assay data!**

### 13.5 Assay scheme to Assay Protocol 1

Reagents	Wells	«Blank»	<b>CAL</b>	Samples
			<b>CONTROL</b>	
<b>CAL</b>		–	50 µL	–
<b>CONTROL</b>				
Samples		–	–	50 µL
<b>CONJ</b>		–	150 µL	150 µL
Incubation No.1	90 min, +18...+25 °C, 500–800 rpm			
<b>WASH P</b> (diluted)	4 x 300 µL			
<b>SUB</b>	100 µL	100 µL	100 µL	100 µL
Incubation No.2	15-30 min, +18...+25 °C, in the dark			
	10 min, +37 °C, 500–800 rpm			
<b>STOP</b>	100 µL	100 µL	100 µL	100 µL
Stirring	1–2 min, +18...+25 °C			
OD measuring	450 nm			
Calculations	Corresponding software			

## 13.6 Assay scheme to Assay Protocol 2 (without shaking)

Wells	«Blank»	<b>CAL</b> <b>CONTROL</b>	Samples
Reagents			
<b>CAL</b> <b>CONTROL</b>	–	50 µL	–
Samples	–	–	50 µL
<b>CONJ</b>	–	150 µL	150 µL
Incubation No.1	90 min, +18...+25 °C (pre-shake for 1-2 minutes at room temperature)		
<b>WASH P</b> (diluted)	4 x 300 µL		
<b>SUB</b>	100 µL	100 µL	100 µL
Incubation No.2	15-30 min, +18...+25 °C, in the dark		
<b>STOP</b>	100 µL	100 µL	100 µL
Stirring	1–2 min, +18...+25 °C		
OD measuring	450 nm		
Calculations	Corresponding software		

## 14. AUTOMATIC TEST (REF 21-02 A)

### 14.1. Material Provided

<b>MP</b>	<b>Microplate:</b> 12 breakable 8-well strips (total 96 wells) coated with anti-testosterone monoclonal antibodies	1 pcs
<b>CONJ</b>	<b>Conjugate:</b> solution containing testosterone conjugated with HRP	18 mL, ready to use
<b>0-5 CAL</b>	<b>Testosterone calibrators:</b> protein-based solutions or lyophilized preparations containing known testosterone concentrations – 0; 1.5; 5; 15; 50, 100 nmol/L (approximate values - do not use for evaluation of real assay data). For lot specific testosterone concentrations see values provided in the Quality Control Sheet	<b>8 vials,</b> 0.5 mL each: <b>CAL</b> 2, 3 – 2 × 0.5 mL; <b>CAL</b> 0, 1, 4, 5 – 1 × 0.5 mL; ready to use or lyophilized
<b>CONTROL</b>	<b>Testosterone control:</b> protein-based solution or lyophilized preparation containing known testosterone concentration. For lot specific ranges of testosterone concentration see Quality Control Sheet	2 × 0.5 mL; ready to use or lyophilized
<b>WASH P 20X</b>	<b>Wash solution P, 20X concentrated:</b> surfactant in buffered saline, sufficient for preparation of 1000 mL of solution	50 mL, concentrated
<b>SUB</b>	<b>Substrate (TMB solution):</b> 3,3',5,5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide	14 mL, ready to use
<b>STOP</b>	<b>Stop solution:</b> 1 N HCl solution	19 mL, ready to use
<b>TRIAL 5000X</b>	<b>Trial solution, 5000X concentrated:</b> solution of detergent	It is delivered by separate order

**Note:** *Extra vials of **CAL** 2, 3 and **CONTROL** are provided for recalibration on reference calibration curve. For detailed information see instructions for use for analyser “Alisei”.*

## 14.2. Equipment and Materials Required but not provided

- Automatic analyser for ELISA kits on microplate (analyser “Alisei”);
- multipurpose polypropylene tubes 12x75 volume 5.5 mL;
- **TRIAL** **5000X** Trial solution, 5000X concentrated\*;
- deionized or distilled water;
- latex or plastic gloves.

## 14.3. Test Procedure

The **Testosterone kit** is designed for 96 tests. This is sufficient for 88 unknowns, 6 calibrators, 1 control and 1 blank (OD of TMB solution) in monoplicates, provided that all the strips are used simultaneously.

### 14.3.1. Reagent Preparation

Allow all the reagents to reach room temperature, and then thoroughly stir.

**Note:** *Preparation of other reagents see in section 8.*

**TRIAL** Working solution of Trial for washing internal hydraulic lines and needles of automatic analyser “Alisei”.

Prepare required volume of **Trial solution** by dilution of the concentrate 5000-fold with distilled or deionized water before analysis:

2 mL **TRIAL** **5000X** + 9998 mL water

Mix thoroughly, avoid foaming.

---

\* Reagent is not included in the kit, it is delivered by separate order.

### **14.3.2. Assay Procedure for Automatic Test**

While using for the procedure analyser “Alisei”, refer to its relative manual.

The analysis on analyser “Alisei” is completely automatic: pipetting of reagents, washing, incubation, OD measurement, analysis of results. The program of calculation testosterone concentration is hardwired in analyser memory.

### **14.4. Data Processing**

#### **14.4.1. Data Reliability (OD 450 nm)**

See criteria in section 13.4.1.

### **14.5. Safety precautions**

- If kit for analyser “Alisei” is used in several separate experiments it is necessary to take reagents from analyser “Alisei” immediately after pipetting them in the wells of all plates because liquid evaporates from vials. Put the reagents in refrigerator;
- Recalibration using calibration curve, obtained with a kit of any other lot, is not permitted.

August, 22, 2019

Astra Biotech GmbH  
Rudower Chaussee 29,  
12489 Berlin, Germany  
Telefon: +49 (0)30 74696509  
E-Mail: [info@astrabiotech.de](mailto:info@astrabiotech.de)  
[www.astrabiotech.de](http://www.astrabiotech.de)