# TG kit

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## 2. INTENDED USE

The **TG** kit is provided for the quantitative determination of thyroglobulin (TG) in human serum.

This test has 2 complete sets:

**REF 24-08** for manual use,

**REF 24-08 A** must be used with ELISA automatic instrument "Alisei" manufactured by NEXT Level S.r.l. hereinafter referred to as analyser "Alisei".

Instructions for use:

Section 13 for manual test,

Section 14 for *analyser*.

**Note 1:** *Take into account that calibrators' nominals can be different for manual and automatic test kits.*

**Note 2:** *We guarantee applications of test:*

**REF 24-08 A** only on analyser "Alisei",

**REF 24-08** only for manual use.

*While using non predefined methods of use, it is under end user responsibility, to make sure that it was appropriately tested for ELISA kits.*

Thyroglobulin (TG) is the main protein produced by thyroid gland. It is the precursor and the deposit form of thyroid hormones thyroxin ( $T_4$ ) and triiodothyronine ( $T_3$ ). TG is an iodinated glycoprotein with a molecular mass of 660 KDa that consists of two subunits.

TG concentration in serum depends on mass of thyroid gland, its functional status and degree of inflammation.

In healthy people TG concentration varies from 2 to 50 ng/mL and may increase about twofold during the pregnancy.

After surgical operation or treatment with radioactive iodine TG concentration remains increased for several weeks. In the population suffering from environmental iodine deficiency mean TG concentration in blood is increased. Since 1994 TG has been regarded by WHO as one of the indicators of endemic goiter.

Serum TG measurement is performed mostly for monitoring of thyroid carcinoma patients. Dynamics of TG concentration is used as a marker for evaluation of treatment efficiency, revealing of possible recurrence and metastases in patients after thyrectomia. This parameter is not used for primary diagnostics of thyroid cancer.

TG concentration in serum may be useful for exact diagnostics in the case of congenital hypothyrosis. Increased TG concentration permits to differentiate subacute thyroiditis from pharmacological thyrotoxicosis (when TG level remains stable). At patients with Grave's disease increase in TG concentration forewarns of recurrence after canceling suppressive therapy.


TG measurement can be distorted by presence of anti-TG autoantibodies in the serum sample. In this case the results may be false-negative. To obtain the reliable results it is recommended always to combine TG and anti-TG autoantibodies determination.

Serial TG determinations for patient's monitoring should be always done by the same method. The physician should be always informed not only of the results of TG determination but also of the type of assay used. When the type of diagnostic kit is changed, the laboratory must perform the appropriate clinical evaluation and compare the range of the results obtained with two types of kits.

### 3. PRINCIPLE OF TEST

The **TG kit** is a “sandwich” type of solid-phase enzyme immunoassay, based on two monoclonal antibodies that are specific for different epitopes of TG molecule. One of these antibodies is conjugated with horseradish peroxidase (HRP); the other is coated onto the inner surface of microwells. TG molecules from the serum sample bind to both immobilized antibody and anti-TG-peroxidase conjugate (Fig. 1). Then the wells are washed with wash solution to remove any material not bound to the inner surface of the wells. Quantity of the bound conjugate is directly proportional to the TG level in sample.

During incubation with substrate the colour is developing. The intensity of the colour is directly proportional to the concentration of TG in specimens or calibrators. TG concentration in the patient sample is read from a standard curve that is processed in each assay.

*Fig. 1 Assay scheme*

## 4. STORAGE AND STABILITY OF THE KIT

The expiration date of the kit is printed on the box label; the expiration date for each component is printed on the respective label.

The **TG kit** should be stored at +2...+8 °C upon receipt, preferably in the original kit box, until the expiration date. Storage at +25 °C is allowed but for no more than 5 days.

Shelf life of the kit is 18 months. After initial opening the kit is stable for 12 months if stored at +2...+8 °C.

If used in several separate experiments, after initial opening kit contents should be stored as follows but never used longer than the expiration date:

- unused strips: in a firmly closed resealable zipper bag, at +2...+8 °C until the expiration date;
- vials with conjugate, calibrators and control (ready-to-use): at +2...+8 °C for no more than 12 months after opening; vials with calibrators and control (reconstituted): at +2...+8 °C for no more than 1 month after opening;
- vial with substrate: at +2...+8 °C until the expiration date, protected from light;
- vials with concentrated Trial, concentrated wash solution, sample diluent, stop solution : at +2...+8 °C until the expiration date;
- wash solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days or at +2...+8°C for no more than 4 weeks in a firmly closed bottle;
- Trial solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days in a firmly closed bottle.

## **Damaged Test Kits**

In case of any severe damage of the test kit or components, Astra Biotech GmbH has to be informed in writing, during one week after receiving the kit. Usage of severely damaged single components for a test run is not recommended.

## **5. SAMPLE COLLECTION AND STORAGE**

Collect blood by venipuncture. After clotting, the serum is separated by centrifugation.

Do not use plasma, haemolyzed (bright red) or lipaemic (milky) serum samples as well as samples containing sodium azide as preservative.

Store serum samples at +2...+8 °C for no more than 2 days. Aliquot and freeze the samples for longer storage ( $\leq -20$  °C). Avoid repeated freezing.

## **6. EXPECTED VALUES**

TG concentrations range determined with **TG kit** in serum samples collected between 9 and 11 a.m. from 223 apparently healthy people (both males and females) at the age of 21–45 did not exceed 55 ng/mL in 99.1% of cases. These limits should be considered as guidelines only.

It is highly recommended for each laboratory to determine its own reference range of TG concentrations.

## **7. QUALITY CONTROL**

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

## 8. REAGENT PREPARATION

Allow all the reagents to reach room temperature (18...25 °C), and then thoroughly stir.

### **MP** Microplate

Keep the **microplate** at room temperature for at least 30 minutes before opening the bag. Place required number of strips into strip holder. Place unused strips onto the resealable zipper bag and reseal duly.

### **CAL** **CONTROL** Calibrators and Control

Liquid calibrators and control are ready to use.

Prepare lyophilized calibrators and control as follows. Gently tap on the vial caps to knock off all the dry matter. Open the vials and carefully place the caps upside down on the clean dry surface. Add 0.5 mL of distilled or deionized water to each vial with lyophilized calibrators and control, close vials with the corresponding caps and leave for 10 min at room temperature (+18...+25 °C) without stirring. Then stir gently avoiding foaming, until the dry matter is completely dissolved. Leave for another 10 minutes at room temperature stirring gently periodically. Make sure that no dry matter is left on the caps and walls of the vials.

**WASH P** Prepare required volume of **wash solution** by dilution of the concentrate 20-fold with distilled or deionized water. For example:

5 mL **WASH P** **20X** + 95 mL water

Mix thoroughly, avoid foaming.

**SUB** Protect **substrate** from direct light.

## 9. SAMPLE PREPARATION

Allow samples to reach room temperature (+18...+25 °C). Stir samples gently in order to ensure homogeneity.

If expected TG concentration in the samples is higher than in **calibrator 5**, the samples should be diluted 50-fold with sample diluent **DIL** in concordance with instructions for use for analyser “Alisei” or manually test before analysis.

The example of manual sample dilution as follows:  
490 µL sample diluent **DIL** + 10 µL serum sample  
vortex or mix thoroughly

## 10. PERFORMANCE CHARACTERISTICS OF THE ASSAY

### 10.1. Calibration – Traceability:

The **TG kit** was calibrated against the International TG Reference Material CRM 457 (*Community Bureau of references, BCR, Brussels, Belgium*).

### 10.2. Specificity:

For the monoclonal antibodies used in the assay, no detectable cross-reactivities were found with the following substances: AFP, FSH, TBG, T<sub>2</sub>, T<sub>3</sub>, T<sub>4</sub>.

Nevertheless, the falsification of any TG determination by specific (e.g. anti-TG antibodies) or non-specific serum factors can not in principle be excluded.

### 10.3. Analytical Sensitivity:

Analytical sensitivity of manual and Alisei **TG assay**, or the lowest detectable concentration that can be distinguished from zero calibrator, is 1.0 ng/mL. It is defined as mean OD of 10 replicates of calibrator 0 plus two standard deviations.

#### 10.4. Measurement Range:

The **TG kits (manual and Alisei kit)** were validated for measurement of TG concentration within the concentration diapason (without dilution) of 1.0 - 300 ng/mL.

#### 10.5 Hook Effect:

For **TG kit** high dose hook effect was not detected for concentrations up to 30 000 ng/mL. High dose hook effect was determined by spiking calibrator 0 with antigen.

#### 10.6 Intra- and Inter-Assay Variation (Precision):

For **inter-assay CV** determination 8 serum samples were assayed in 9 replicates each. The results are shown below.

##### Manual kit

Sample	Mean TG concentration, ng/mL	Intra-assay CV	
		SD	CV, %
1	1.6	0.12	7.5
2	5.0	0.12	2.4
3	10.0	0.26	2.6
4	20.4	0.64	3.1
5	34.6	0.76	2.2
6	37.1	1.47	4.0
7	48.3	1.30	2.7
8	228	6.6	2.9

##### Alisei kit

Sample	Mean TG concentration, ng/mL	Intra-assay CV	
		SD	CV, %
1	2.0	0.11	5.5
2	7.5	0.38	5.0
3	12.8	0.56	4.3
4	21.6	1.21	5.6

5	35.1	1.90	5.4
6	85.9	2.52	2.9
7	143.9	2.93	2.0
8	268.2	5.87	2.2

For **inter-assay CV** determination 8 samples were assayed 3 times by different operators with 1-week interval. Each specimen was assayed in 9 replicates. The results are shown below.

#### Manual kit

Sample	Mean TG concentration, ng/mL			Inter-assay precision	
	1 assay	2 assay	3 assay	SD	CV, %
1	1.11	1.51	1.36	0.200	14.1
2	2.02	2.45	2.06	0.240	10.5
3	6.46	5.96	6.92	0.480	7.5
4	13.0	12.3	12.8	0.36	2.9
5	38.8	38.9	41.6	1.59	3.9
6	47.8	46	50.2	2.11	4.4
7	111	128	140	14,6	10.9
8	286	275	273	7.0	2.6

#### Alisei kit

Sample	Mean TG concentration, ng/mL			Inter-assay precision	
	1 assay	2 assay	3 assay	SD	CV, %
1	2.0	2.3	2.1	0.15	7.2
2	7.5	7.0	7.9	0.45	6.0
3	12.8	14.3	15.0	1.12	8.0
4	21.6	19.4	19.0	1.40	7.0
5	35.1	39.6	37.8	2.26	6.0
6	85.9	79.1	89.5	5.28	6.2
7	143.9	159.7	155.5	8.18	5.3
8	268.2	257.9	280.6	11.37	4.2

## 11. LIMITATION OF THE METHOD

Any clinical diagnosis should not be based on the results of in vitro diagnostic methods alone. To state a diagnosis, the physician is supposed to consider all the available clinical and laboratory findings.

TG measurement may be affected by presence of anti-TG autoantibodies in the serum sample. In this case the results may be false-negative. To obtain the reliable results it is recommended always to combine TG and anti-TG autoantibodies determination.

Serial TG determinations for patient monitoring should be always done by the same method. The physician should be always informed not only of the results of TG determination but also of the type of assay used. When the type of diagnostic kit is changed, the laboratory must perform the appropriate clinical evaluation and compare the range of the results obtained with two types of kits.

## 12. SAFETY PRECAUTIONS

- **This kit is for in vitro diagnostic use only.** Operator should follow the manual closely in order to ensure reliable data. The manual is valid for the present kit only, within the listed composition. Any substitution of kit components is not allowed by CE regulations.
- Do not use kit or its components after the expiration date indicated on the label. Take into account stability period for reconstituted reagents.
- Do not mix or use together reagents from different lots of the kit except substrate, stop solution and wash solution.
- Do not use substrate, stop solution and wash solution supplied by other vendors.
- Use only “P”-labeled Astra Biotech wash solution.
- Note that stop solution is 1 N HCl solution. Avoid contacts with skin and mucosa. In case of contact rinse affected area thoroughly with plenty of water and seek medical advice.
- Take into account the following common procedural notes:
  - always pipette reagents into wells immediately after washing procedure;
  - avoid contamination of the solutions;
  - in case of partial use of the kit, dispense only required volume of the reagent into the tray;
  - do not pour unused reagents back into the original vials;
  - avoid exposure to direct sunlight during incubations;
  - always pipette reagents in the same order to minimize reaction time differences between wells; the total dispensing time for the calibrators, control and samples must not exceed 15 min;

- the incubation temperature for all the immunological reactions must be kept at +37 °C;
- do not touch the bottom of the wells;
- calibrators should be measured in each separate assay. It is also recommended to measure each time TG concentration in the control.
- remove the adhesive foil carefully to avoid contamination and don't use the adhesive foil repeatedly.
- TMB solution should be colourless. Light colouring of solution is admissible. Avoid direct exposure of substrate to sunlight.
-  Source materials of human origin used for kit components preparation were tested and found negative for HBsAg, anti-HIV and anti-HCV antibodies. However, none of known laboratory test guarantees absence of these viral agents. Therefore, all kit components and patient's samples should be handled as potentially hazardous.
-  After usage strips, calibrators, control, specimens and all consumables which contacted with specimens during handling, storage or assay (tubes, vials, gloves, pipette tips etc.) should be collected separately and sterilized by autoclaving. Instead of autoclaving pipette tips may be sterilized by disinfectant treatment. After sterilization all components and expendable materials may be utilized as non-dangerous garbage. Other components of the kit should be discarded into conventional garbage.
-  During manual washing procedure do not discard the contents of the wells directly to drainage. Use a container with disinfectant solution.

-  As the kit contains potentially hazardous material, the following precautions should be taken:
  - do not smoke, eat or drink while performing the assay;
  - always use protective gloves;
  - never pipette material by mouth;
  - in case of spilling, wipe up the spills promptly and wash affected area thoroughly using decontaminant.
- GLP including all general and individual regulations should be applied for the kit usage.

- As the kit contains irritant (**CONJ**, **DIL**, **CAL**, **CONTROL**), the following precautions should be observed:
- P261 - Avoid breathing spray;
  - P272 - Contaminated work clothing should not be allowed out of the workplace;
  - P280 - Wear protective gloves/protective clothing/eye protection;
  - P302+P352 - IF ON SKIN: Wash with plenty of soap and water;
  - P333+P313 - If skin irritation or rash occurs: Get medical advice/attention;
  - P363 - Wash contaminated clothing before reuse;
  - P501 - Dispose of contents/container in accordance with national regulation.

Precautionary statements according to Regulation EC № 1272/2008.

## 13. MANUAL TEST (REF 24-08)

### 13.1. Material Provided

<p><b>MP</b></p>	<p><b>Microplate:</b> 12 breakable 8-well strips (total 96 wells) coated with anti-TG monoclonal antibodies</p>	<p>1 pcs</p>
<p><b>CONJ</b></p>	<p><b>Conjugate:</b> solution contains anti-TG monoclonal antibodies conjugated with HRP</p>	<p>14 mL, ready to use</p>
<p><b>0-5 CAL</b></p>	<p><b>TG calibrators:</b> protein-based solution or lyophilized preparations containing known TG concentrations - 0; 3; 10; 30; 100; 300 ng/mL (approximate values - do not use for evaluation of real assay data). The concentrations of calibrators may be different for Assay Protocol 1 and 2. For lot specific TG concentrations see values for respective Assay Protocol provided in the Quality Control Sheet.</p>	<p>6 vials, 0.5 mL each; ready to use or lyophilized</p>
<p><b>CONTROL</b></p>	<p><b>TG control:</b> protein-based solution or lyophilized preparation containing known TG concentration The range of TG concentration may be different for Assay Protocol 1 and 2. For lot specific ranges of TG concentration see values for respective Assay Protocol provided in the Quality Control Sheet.</p>	<p>0.5 mL, ready to use or lyophilized</p>
<p><b>WASH P 20X</b></p>	<p><b>Wash solution P, 20X concentrated:</b> surfactant in buffered saline, sufficient for preparation of 560 mL of solution</p>	<p>2x14 mL, concentrated</p>
<p><b>SUB</b></p>	<p><b>Substrate (TMB solution):</b> 3,3',5,5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide</p>	<p>14 mL, ready to use</p>

<b>STOP</b>	<b>Stop solution:</b> 1 N HCl solution	14 mL, ready to use
<b>DIL</b>	<b>Sample diluent:</b> protein-based solution	3.0 mL, ready to use
	<b>Adhesive foil</b>	2x1 foil (optional)

### 13.2. Equipment and Materials Required but not provided

- 1-channel calibrated variable precision pipettes, with disposable tips;
- 8-channel calibrated variable precision pipette, with disposable tips;
- microplate incubator (+37 °C) or microplate incubator/shaker (+37 °C, shaking speed 500–800 rpm);
- manual or automatic equipment for rinsing wells;
- calibrated microplate reader (450 nm);
- vortex tube mixer;
- deionized or distilled water;
- graduated beaker and cylinder of appropriate volume;
- latex or plastic gloves;
- trays for pipetting reagents with 8-channel pipette;
- disinfectant;
- absorbent material (for manual wash).

### 13.3. Test Procedure

The **TG kit** is designed for 96 tests. This is sufficient for 40 unknowns, 6 calibrators, 1 control and 1 blank (OD of TMB solution) in duplicates, provided that all the strips are used simultaneously.

### 13.3.1. Assay Procedure

**Note:** *Please check in the Quality Control Sheet supplied with the kit for which Assay Protocol your kit has been validated. If your kit has been validated for Assay Protocol 1, please follow the Assay Protocol 1 described below. If your kit has been additionally validated for Assay Protocol 2, please consider the Supplement 1 to the Instructions for use provided with the kit.*

#### Assay Protocol 1

(see also Assay scheme to Assay Protocol 1, section 13.5)

**Note:** *Consider concentrations of calibrators and control range for Assay Protocol 1 provided in the Quality Control Sheet.*

*All samples should be tested in duplicates.*

- A.** Pipette **100 µL conjugate** **CONJ** into each well, **Leave wells A1-A2 empty for blank!**
- B.** Pipette **50 µL of calibrators** **CAL**, **control** **CONTROL** and prepared **patient's samples** in duplicates; **except wells A1-A2.**

**Note:** *total time of dispensing must not exceed 15 minutes, otherwise the test result may be unreliable, because the time of incubation will substantially vary for different samples.*

**Note:** *If an adhesive foil is used, put it now onto the plate to seal the cavities.*

- C.** Incubate for **60 minutes at +37 °C while shaking (500–800 rpm)**.

**Note:** *If an adhesive foil is used, remove it now from the plate.*

- D.** Wash 5 times, as described below.

- E.** Pipette **100 µL substrate** **SUB** into each well (including blank); incubate strips **at room temperature (+18...+25 °C) in the dark for 15-30 minutes**, depending on the colour intensity, or **10 minutes while shaking (500-800 rpm) at +37 °C**.

- F.** Pipette **100 µL stop solution** **STOP** to all the wells in the same sequence and at the same speed as used for dispensing TMB substrate. **Shake for 1–2 min at room temperature**.

- G.** Read OD at **450 nm within 20 min**.

### **13.3.2. Wash Procedure**

It is advisable to use an automatic microplate washer set at 5 wash cycles and a volume of 300 µL of wash solution per well per cycle. If an automatic washer is not available, the wash procedure can be carried out manually:

- remove the contents of the wells into a container with disinfectant;
- dispense 300 µL of wash solution (section 8) into each well, shake the plate carefully for 5–10 sec and remove the contents of the wells; repeat 5 times;

- strike the wells sharply on absorbent material to remove any liquid residue.

### 13.4. Data Processing

If the reader cannot be adjusted to zero using the substrate blank in wells A1-A2, subtract mean OD value of wells A1-A2 from all OD values before further calculations.

*Example:*

OD (Cal 5) measured = 2.28 and OD (blank) = 0.06;

OD (Cal 5) calculated =  $2.28 - 0.06 = \underline{2.22}$

#### 13.4.1. Data Reliability (for OD measured at 450 nm)

The data should meet the following criteria:


- average blank OD (in wells A1-A2)  $\leq 0.100$ ;
- average OD of Cal 5  $\geq 1.500$  (after blank subtraction);
- control's concentration must fall within the acceptability range that are shown on the Quality Control Sheet.

If the data obtained do not meet the criteria, the results are considered unreliable and the test should be repeated.

#### 13.4.2. Quantitative Determination

Specialized software for quantitative determination is recommended. Mean OD values of the calibrators at 450 nm versus their respective TG concentrations using **4PL** or **5PL fit** (Fig. 2 Typical standard curve). Calculate concentration of TG in samples using standard curve.

Any extrapolation of the standard curve to TG concentration above the nominal value of calibrator 5 is forbidden. In this case the sample should be diluted 50-fold with sample diluent and re-tested. Multiply the measured concentration of pre-diluted samples by dilution factor (50-fold).

*Fig. 2 Example of typical standard curve*  
**Do not use for evaluation of real assay data!**

### 13.5. Assay scheme to Assay Protocol 1

Reagents	Wells		
	«Blank»	<b>CAL</b> <b>CONTROL</b>	Samples
<b>CONJ</b>	–	100 µL	100 µL
<b>CAL</b> <b>CONTROL</b>	–	50 µL	–
Samples	–	–	50 µL
Incubation №1	60 min, +37 °C, 500–800 rpm		
<b>WASH P</b> (diluted)	5 x 300 µL		
<b>SUB</b>	100 µL	100 µL	100 µL
Incubation №2	15–30 min, +18...+25 °C, in the dark		
	10 min, +37 °C, 500–800 rpm		
<b>STOP</b>	100 µL	100 µL	100 µL
Stirring	1–2 min, +18...+25 °C		
OD measuring	450 nm		
Calculations	Corresponding software		

## 14. AUTOMATIC TEST (REF 24-08 A)

### 14.1. Material Provided

<b>MP</b>	<b>Microplate:</b> 12 breakable 8-well strips (total 96 wells) coated with anti-TG monoclonal antibodies	1 pcs
<b>CONJ</b>	<b>Conjugate:</b> solution contains anti-TG monoclonal antibodies conjugated with HRP	14 mL, ready to use
<b>0-5 CAL</b>	<b>TG calibrators:</b> protein-based solution or lyophilized preparations containing known TG concentrations - 0; 3; 10; 30; 100; 300 ng/mL (approximate values - do not use for evaluation of real assay data). For lot specific TG concentrations see values provided in the Quality Control Sheet	8 vials, 0.5 mL each: <b>CAL</b> 1, 4 – 2 × 0.5 mL; <b>CAL</b> 0, 2, 3, 5 – 1 × 0.5 mL; ready to use or lyophilized
<b>CONTROL</b>	<b>TG control:</b> protein-based solution or lyophilized preparation containing known TG concentration. For lot specific ranges of TG concentration see Quality Control Sheet	2 × 0.5 mL, ready to use or lyophilized
<b>WASH P 20X</b>	<b>Wash solution P, 20X concentrated:</b> surfactant in buffered saline, sufficient for preparation of 1000 mL of solution	50mL, concentrated
<b>DIL</b>	<b>Sample diluent:</b> protein-based solution	3.0 mL, ready to use
<b>SUB</b>	<b>Substrate (TMB solution):</b> 3,3',5,5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide	14 mL, ready to use
<b>STOP</b>	<b>Stop solution:</b> 1 N HCl solution	19 mL, ready to use
<b>TRIAL 5000X</b>	<b>Trial solution, 5000X concentrated:</b> solution of detergent	It is delivered by separate order

**Note:** Extra vials of **CAL** 1, 4 and **CONTROL** are provided for recalibration on reference calibration curve. For detailed information see instructions for use for analyser “Alisei”.

## 14.2. Equipment and Materials Required but not provided

- Automatic analyser for ELISA kits on microplate (analyser “Alisei”);
- multipurpose polypropylene tubes 12x75 volume 5.5 mL;
- **TRIAL 5000X** Trial solution, 5000X concentrated\*;
- deionized or distilled water;
- latex or plastic gloves.

## 14.3. Test Procedure

The **TG kit** is designed for 96 tests. This is sufficient for 88 unknowns, 6 calibrators, 1 control and 1 blank (OD of TMB solution) in monoplicates, provided that all the strips are used simultaneously.

### 14.3.1. Reagent Preparation

Allow all the reagents to reach room temperature (+18...+25 °C), and then thoroughly stir.

**Note:** Preparation of other reagents see in section 8.

**TRIAL** Working solution of Trial for washing internal hydraulic lines and needles of automatic analyser “Alisei”.

Prepare required volume of **Trial solution** by dilution of the concentrate 5000-fold with distilled or deionized water before analysis:

---

\* Reagent is not included in the kit, it is delivered by separate order.

2 mL **TRIAL** **5000X** + 9998 mL water

Mix thoroughly, avoid foaming.

### **14.3.2. Assay Procedure for Automatic Test**

While using for the procedure analyser “Alisei”, refer to its relative manual.

The analysis on analyser “Alisei” is completely automatic: pipetting of reagents, washing, incubation, OD measurement, analysis of results. The program of calculation TG concentration is hardwired in analyser memory.

## **14.4. Data Processing**

### **14.4.1. Data Reliability (OD measured at 450 nm)**

See criteria in section **13.4.1.**

## **14.5. Safety precautions**

- If kit for analyser “Alisei” is used in several separate experiments it is necessary to take reagents from analyser “Alisei” immediately after pipetting them in the wells of all plates because liquid evaporates from vials. Put the reagents in refrigerator;  
Recalibration using calibration curve, obtained with a kit of any other lot, is not permitted.

August 26, 2019

Astra Biotech GmbH  
Rudower Chaussee 29  
12489 Berlin, Germany  
Telefon: +49 (0)30 74696509

E-Mail: [info@astrabiotech.de](mailto:info@astrabiotech.de)  
[www.astrabiotech.de](http://www.astrabiotech.de)