Treponema-Total Antibodies kit

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2. INTENDED USE

Treponema-Total Antibodies kit is provided for the **qualitative** determination of IgG, IgA and IgM (total Ab) to *Treponema pallidum* (*T. pallidum*) in human serum, plasma and liquor.

Treponema-Total Antibodies kit has 4 complete sets for qualitative determination both manually and with automatic analyser «Alisei»:

REF 60-01-96	96 tests
REF 60-01-192	192 tests
REF 60-01-480	480 tests
REF 60-01-960	960 tests

Instructions for use are described in:

Paragraph 13	for manual test,
Paragraph 14	for analyser.

Syphilis is a chronic infection caused by *Treponema pallidum*. *T. pallidum* is an anaerobic gram-negative obligate bacteria of the Spirochaeta family.

It usually transmitted from human to human during sexual acts via the mucosa. Although the disease may be transmitted by transfusion of infected blood. Intrauterine infection also occurs.

The infection progresses through several distinct stages – primary, secondary, early latent, late latent, and tertiary.

The serological diagnosis of syphilis is performed demonstrating the presence of significant levels of specific anti-treponemal antibodies detectable within approximately 2-4 weeks post infection.

Note 1: Take into account that *Ccrit* can be different for manual and automatic test kits.

Note 2: In case of automated application we guarantee applications of test only on analyser “Alisei”. While using non predefined methods of use, it is under end user responsibility.

3. PRINCIPLE OF THE TEST

The **Treponema-Total Antibodies** kit is a one-step solid phase enzyme immunoassay. Anti-*T. pallidum* total Ab from a sample or a control bind to the recombinant antigens TpN15, TpN17, TpN47 coated onto the inner surface of the microplate wells. The same recombinant antigens TpN15, TpN17, TpN47 conjugated with horseradish peroxidase (HRP) bind to the “antigen-antibody” complex (Fig. 1). During the incubation with TMB solution the colour is developing. The intensity of the colour correlates to the concentration of anti- *T. pallidum* total antibodies in specimens or control. After OD measurement, diagnoses of the specimens are determined on the basis of the Cut-off value.


Fig. 1. Assay scheme

4. STORAGE AND STABILITY OF THE KIT

The expiration date of the kit is printed on the box label; expiration date for each component is printed on the respective label.

Treponema-Total Antibodies kit should be stored at +2...+8 °C, preferably in the original kit box, until the expiration date. Storage at +25 °C is allowed but for no more than 7 days.

Shelf life of the kit is 12 months.

After initial opening the kit is stable for 2 months if stored at +2...+8 °C.

If used in several separate experiments, kit contents should be stored as follows:

MP unused strips: in a firmly closed ziplock pouch at +2...+8 °C for no more than 2 months after opening;

NEG, **POS**, **CONJ**, **SUB** opened vials with controls, conjugate and substrate: at +2...+8 °C for no more than 2 months after opening;

TRIAL **5000X**, **WASH H** **20X**, **STOP** opened vials with concentrated Trial, concentrated wash solution and stop solution: at +2...+8 °C until the expiration date;

TRIAL Trial solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days, in a firmly closed bottle.

WASH H Wash solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days or at +2...+8 °C for no more 4 weeks in a firmly closed bottle.

Damaged Test Kits

In a case of any severe damage of the test kit or components Astra Biotech GmbH has to be informed in writing, during one week after receiving the kit. Usage of severely damaged single components for a test run is not recommended.

5. SAMPLE COLLECTION AND STORAGE

Collect blood by venipuncture in a tube without anticoagulants (for serum) or evacuated tube containing EDTA, citrate or heparin (for plasma). Allow blood to clot for serum samples. Centrifuge the specimens to separate serum or plasma from the cells.

Do not use haemolyzed (bright red) or lipaemic (milky) blood samples as well as samples containing sodium azide as preservative.

Store serum and plasma samples at +2...+8 °C for no more than 5 days, and aliquot and freeze samples for longer storage (-20 °C and lower). Avoid repeated freezing and thawing.

Collect liquor in a tube without dust, EDTA or fluoride.

Store liquor at +2...+8 °C for no more than 7 days, do not freeze liquor samples.

6. QUALITY CONTROL

It is recommended to use control samples according to the state and federal regulations. The use of control samples is advised to assure the day-to-day validity of results.

7. REAGENT PREPARATION

Allow all the reagents to reach room temperature, and then thoroughly stir.

MP Keep **microplate** at room temperature (+18...+25 °C) for at least 30 minutes before opening the bag. Place required number of strips onto strip holder. Place unused strips into the resealable zipper bag and reseal duly.

NEG POS Controls are ready to use.

CONJ Conjugate is ready to use.

WASH H Prepare required volume of **wash solution** by dilution of the concentrate 20-fold with distilled or deionized water. For example:

5 mL of **WASH H 20X** + 95 mL of water.

Mix thoroughly, avoid foaming. Keep ready-to-use wash solution firmly closed.

TRIAL Working solution of Trial for washing internal hydraulic lines and needles of automatic analyser "Alisei". Prepare required volume of **Trial solution** by dilution of the concentrate 5000-fold with distilled or deionized water before analysis:

2 mL of **TRIAL 5000X** + 9998 mL of water.

Mix thoroughly, avoid foaming.

SUB

Substrate is ready to use. Protect **substrate** from direct light.

STOP

Stop solution is ready to use.

8. SAMPLE PREPARATION

Allow samples to reach room temperature. Stir samples gently in order to ensure homogeneity

9. PERFORMANCE CHARACTERISTICS OF THE ASSAY

9.1. Diagnostic Specificity

450 samples classified negative for anti- *T. pallidum* total Ab by reference EIA kit were tested using the **Treponema-Total Antibodies** kit. Diagnostic specificity was found to be 99,1%.

9.2. Diagnostic Sensitivity

47 samples positive for anti- *T. pallidum* total Ab from commercial panels were tested using the **Treponema-Total Antibodies** kit. Diagnostic sensitivity was found to be 100%.

206 samples classified positive for anti- *T. pallidum* total Ab by reference EIA kit were tested using the **Treponema-Total Antibodies** kit. Diagnostic sensitivity was found to be 98,5%.

9.3. Intra- and inter-assay variation (Precision)

For intra-assay CV determination, 6 serum samples were run using one lot of the **Treponema-Total Antibodies** kit accordingly NCCLS EP05-A2 «Evaluation of Precision Performance of Quantitative Measurement Methods». The results are shown below.

Sample	Ratio	Intra-assay CV	
		SD	CV, %
1	1.9	0.068	3.6
2	2.5	0.090	3.6

3	3.8	0.122	3.2
4	4.9	0.181	3.7
5	6.1	0.244	4.0
6	13.7	0.342	2.5

For inter-assay CV determination, 6 serum samples were run using 3 different lots of the **Treponema-Total Antibodies** kit accordingly NCCLS EP05-A2 «Evaluation of Precision Performance of Quantitative Measurement Methods». The results are shown below.

Sample	Ratio			Inter-assay CV	
	lot 1	lot 2	lot 3	SD	CV, %
1	1.7	1.9	1.8	0.171	9.5
2	2.6	2.4	2.2	0.252	10.5
3	3.6	3.9	3.3	0.331	9.2
4	4.6	4.7	4.8	0.446	9.5
5	5.7	5.9	6.1	0.496	8.4
6	12.8	13.0	12.6	1.152	9.0

9.4. High dose hook effect

For **Treponema-Total Antibodies** kit no **high dose hook effect** was detected.

10. LIMITATION OF THE METHOD

Any clinical diagnosis should not be based on the results of in vitro diagnostic methods alone. For diagnosis establishment, a physician is supposed to consider all available clinical and laboratory findings.

11. DATA RELIABILITY

For OD measured at 450 nm

The data should meet the following criteria:

- blank OD (in well A1) ≤ 0.1 ;
- OD negative control **NEG** ≤ 0.2 (after blank subtraction);
- OD positive control **POS** ≥ 0.5 (after blank subtraction).

If the data obtained do not meet the criteria, the results are considered unreliable and the test should be repeated.

12. SAFETY PRECAUTIONS

- **This kit is for in vitro diagnostic use only.** Operator should follow the manual closely in order to ensure reliable data. The manual is valid for the present kit only, within the listed composition. Any substitution of kit components is not allowed by CE regulations.
- Do not use the kit or its components after the expiration date indicated on the label. Take into account stability period for reconstituted reagents.
- Do not mix or use together reagents from different lots of the kit except substrate, stop solution and wash solution.
- Do not use substrate, stop solution and wash solution supplied by other vendors.
- Use only “H”-labeled wash solution for REF 200-50, 200-51, 200-52 and 200-53.
- Note that stop solution is 1N HCl solution. Avoid contacts with skin and mucosa. In case of contact rinse affected area thoroughly with plenty of water and seek medical advice.

As the kit contains irritant (**CONJ**, **NEG**, **POS**) the following precautions should be observed:

- P261 - Avoid breathing spray;
- P272 - Contaminated work clothing should not be allowed out of the workplace;
- P280 - Wear protective gloves/protective clothing/eye protection;
- P302+P352 - IF ON SKIN: Wash with plenty of soap and water;
- P333+P313 - If skin irritation or rash occurs: Get medical advice/attention;
- P363 - Wash contaminated clothing before reuse;
- P501 - Dispose of contents/container in accordance with national regulation.

Precautionary statements according to Regulation EC № 1272/2008.

- Take into account the following common procedural notes:
 - always pipette reagents into wells immediately after washing procedure;
 - avoid contamination of the solutions;
 - in case of partial use of the kit, dispense only required volume of the reagent into the tray; do not pour unused reagents back into the original vials;
 - avoid exposure to direct sunlight during incubations;
 - always pipette reagents in the same order to minimize reaction time differences between wells; the total dispensing time for the calibrators and samples must not exceed 15 min;
 - the incubation temperature for all the immunological reactions must be kept at 37 °C;

- do not touch the bottom of the wells;
- controls should be measured in each separate assay.
- TMB solution should be colourless. Light colouring of solution is admissible. Avoid direct exposure of substrate to sunlight.
-  Source materials of animal origin used for kit components preparation were tested and found negative for BSE; Vesicular stomatitis virus and Bluetongue virus. However, none of known laboratory test guarantees absence of these viral agents. Therefore, all kit components should be handled as potentially hazardous.
-  Source materials of human origin used for kit components preparation were tested and found negative for HBsAg, anti-HIV and anti-HCV antibodies. However, none of known laboratory test guarantees absence of these viral agents. Therefore, all kit components and patient's samples should be handled as potentially hazardous.
-  After usage strips, controls, conjugate, specimens and all consumables which contacted with specimens during handling, storage or assay (tubes, vials, gloves, pipette tips etc.) should be collected separately and sterilized by autoclaving. Instead of autoclaving pipette tips may be sterilized by disinfectant treatment. After sterilization all components and expendable materials may be utilized as non-dangerous garbage. Other components of the kit should be discarded into conventional garbage.
-  During manual washing procedure do not discard the contents of the wells directly to drainage. Use a container with disinfectant solution.

-  As the kit contains potentially hazardous material, the following precautions should be taken:
 - do not smoke, eat or drink while performing the assay;
 - always use protective gloves;
 - never pipette material by mouth;
 - in case of spilling, wipe up the spills promptly and wash affected area thoroughly using decontaminant.
- GLP including all general and individual regulations should be applied for the kit usage.

13. MANUAL TEST

13.1 Materials Provided

see paragraph 14.1

13.2. Equipment and materials required but not provided

- 1-channel calibrated variable precision pipettes, with disposable tips;
- 8-channel calibrated variable precision pipette, with disposable tips;
- microplate incubator/shaker (+37 °C, 400–800 rpm);
- manual or automatic equipment for rinsing wells;
- microplate calibrated reader (450 nm);
- vortex tube mixer;
- deionized or distilled water;
- graduated beaker and cylinder of appropriate volume;
- latex or plastic gloves;
- trays for pipetting reagents using 8-channel pipette;
- disinfectant;
- absorbent material (for manual wash).

13.3. Test procedure

The **Treponema-Total Antibodies kit** is supplied in four bundling:

- bundling for **96 determinations** is sufficient for the qualitative assay of 45 unknowns in duplicates or 90 unknowns in monoplicates*

- bundling for **192 determinations** is sufficient for the qualitative assay of 90 unknowns in duplicates or 180 unknowns in monoplicates*
- bundling for **480 determinations** is sufficient for the qualitative assay of 225 unknowns in duplicates or 450 unknowns in monoplicates*
- bundling for **960 determinations** is sufficient for the qualitative assay of 450 unknowns in duplicates or 900 unknowns in monoplicates*

*provided that all the strips are used simultaneously.

Note 1: *Pipetting in duplicates is recommended by manufacturer.*

Note 2: *Take into account it is necessary to determine OD of negative control , positive control and TMB on each plate.*

13.3.1. Assay Procedure

(See assay scheme, page 43)

- A.** Pipette **90 µL** conjugate **CONJ**; **except well A1**(blank).
- B.** Pipette **10 µL** negative control **NEG** in duplicate, positive control **POS** in triplicate and **patient's samples in duplicates or monoplicates**. Mix thoroughly the components in wells by triple pipetting; except well A1.
- C.** Incubate for **30 minutes at +37°C while shaking (400–800 rpm)** or for **30 minutes at +37°C without shaking**.
- D.** Wash 8 times, as described below.
- E.** Pipette **100 µL substrate SUB** into each well (including blank); incubate for **10 minutes while shaking (400–800 rpm) at +37 °C** or for **15 minutes at +18...+25°C without shaking**.
- F.** Pipette **100 µL stop solution STOP** into each well (including blank) in the same sequence and at the same speed as used for dispensing TMB substrate. **Shake for 1–2 min at room temperature**.
- G.** Read the absorbance at **450 nm** within 20 min.

13.3.2. Wash Procedure

It is advisable to use an automatic microplate washer set at 8 wash cycles with 30 seconds delay between cycles and a volume of 300 µL of wash solution per well per cycle. Check, that the wells are completely filled in, and make sure that the aspiration is carried out fully.

If an automatic washer is not available, the wash procedure can be carried out manually as follows:

- remove the contents of the wells into container with disinfectant;
- dispense 300 μL of wash solution, prepared according to section 7, into each well, wait for 30 sec and remove the contents of the wells; repeat 8 times;
- strike the wells sharply on absorbent material to remove any liquid residue.

13.4. Data processing

If the reader cannot be adjusted to zero using the substrate blank in well A1, subtract OD value of well A1 (blank) from all OD values before further calculations.

Example:

OD (sample 1) measured = 2.28 and OD (blank) = 0.06;

OD (sample 1) calculated = $2.28 - 0.06 = 2.22$.

Calculate the mean OD of the negative control **NEG**.

Only **NEG** values with $\text{OD} < 0.2$ should be used for the calculation of mean OD **NEG**.

Calculate Cut-off using the formula 1:

Cut-off = mean OD **NEG** + C_{crit}

C_{crit} – is a coefficient defined by the manufacturer. C_{crit} for each lot of the kit is specified in the QUALITY CONTROL SHEET (QCS) inserted into the kit.

Calculate Ratio using the formula 2:

$$\text{Ratio} = \frac{\text{Sample OD}}{\text{Cut-off}}$$

13.5. Interpretation of Ratio

Ratio	Result
< 0.9	negative
≥ 0.9 to ≤ 1.1	borderline
>1.1	positive

In cases of borderline test results, an additional patient sample should be taken 7-14 days later and retested in parallel with the first patient sample.

The results of both samples allow proper evaluation of titer changes.

13.6. Assay scheme

Reagents	Wells	«Blank»	NEG	POS	Samples
	CONJ	–	–	90 µL	90 µL
NEG POS	–	–	10 µL triple pipetting	–	–
Samples	–	–	–	–	10 µL triple pipetting
Incubation No.1	30 min, +37 °C, 400–800 rpm or 30 min, +37 °C, without shaking				
WASH H (diluted)	8 x 300 µL				
SUB	100 µL	100 µL	100 µL	100 µL	100 µL
Incubation No.2	10 min, +37 °C, 400–800 rpm or 15 min, +18...+25 °C, without shaking				
STOP	100 µL	100 µL	100 µL	100 µL	100 µL
Stirring	1–2 min, +18...+25 °C				
OD measuring	450 nm				
Calculations	Corresponding software				

14. AUTOMATIC TEST

14.1. Material Provided

		REF 60-01-96	REF 60-01-192	REF 60-01-480	REF 60-01-960
MP	Microplate coated with recombinant antigens TpN15, TpN17, TpN47	12 breakable 8-well strips (total 96 wells)			
		1 pcs	2 pcs	5 pcs	10 pcs
NEG	Negative control protein-based solution	0,8 mL, ready to use	0,8 mL, ready to use	2 x 0,8 mL, ready to use	4 x 0,8 mL, ready to use
POS	Positive control protein-based solution containing anti-T.pallidum IgG and IgM	0,8 mL, ready to use	0,8 mL, ready to use	2 x 0,8 mL, ready to use	4 x 0,8 mL, ready to use
CONJ	Conjugate solution containing recombinant antigens TpN15, TpN17, TpN47 conjugated with HRP	16 mL, ready to use	2 x 16 mL, ready to use	4 x 16 mL, ready to use	8 x 16 mL, ready to use
WASH H 20X	Wash solution H, 20X concentrated: surfactant in buffered saline	50 mL, sufficient for preparation of 1000 mL of solution	2 x 50 mL, sufficient for preparation of 2 x 1000 mL of solution	3 x 100 mL, sufficient for preparation of 3 x 2000 mL of solution	5 x 100mL, sufficient for preparation of 5 x 2000 mL of solution
SUB	Substrate (TMB solution): 3,3',5,5'-	14 mL, ready to	2 x 14 mL, ready to	100 mL, ready to	2 x 100 mL, ready

	tetramethyl benzidine solution in citrate buffer containing hydrogen peroxide	use	use	use	to use
STOP	Stop solution: 1 M HCl solution	50 mL, ready to use	2 x 50 mL, ready to use	3 x 100 mL, ready to use	5 x 100 mL, ready to use
TRIAL 5000X	Trial solution, 5000X concentrated: solution of detergent	Delivered by separate order			

14.2. Equipment and Materials Required but not provided

- Automatic analyser for ELISA kits on microplate (analyser “Alisei”);
- multipurpose polypropylene tubes 12x75 volume 5.5 mL;
- **TRIAL 5000X** Trial solution, 5000X concentrated;
- deionized or distilled water;
- latex or plastic gloves.

14.3. Test procedure

Treponema-Total Antibodies kit is supplied in four bundling:

- bundling for **96 determinations** is sufficient for the qualitative assay of 90 unknowns in monoplicates*
- bundling for **192 determinations** is sufficient for the qualitative assay of 180 unknowns in monoplicates*

- bundling for **480 determinations** is sufficient for the qualitative assay of 450 unknowns in monoplicates*
- bundling for **960 determinations** is sufficient for the qualitative assay of 900 unknowns in monoplicates*

* provided that all the strips are used simultaneously.

14.3.1. Assay Procedure for Automatic Test

While using for the procedure analyser “Alisei”, refer to its relative manual.

The analysis on EIA analyser “Alisei” is completely automatic: pipetting of reagents, washing, incubation, OD measurement, analysis of results. The program of anti-*T. pallidum* total Ab determination is hardwired in analyser memory.

14.4. Data Reliability

See criteria in paragraph 11.

14.5. Safety precautions

See paragraph 12.

If the kit is used in several separate experiments it is necessary to take reagents from analyser “Alisei” immediately after pipetting them in the wells of all plates because liquid evaporates from vials. Put the reagents into refrigerator.

Mai, 29, 2019

Astra Biotech GmbH
Rudower Chaussee 29,
12489 Berlin, Germany
Telefon: +49 (0)30 74696509
E-Mail: info@astrabiotech.de
www.astrabiotech.de