TSH 3rd generation kit

3.Generation ELISA Testkit zur quantitativen Bestimmung des Thyreoidea-stimulierenden Hormons in humanem Blutserum

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2. INTENDED USE

The **TSH (third generation)** kit is provided for the quantitative determination of low concentrations of thyroid stimulating hormone (TSH) in human serum.

TSH is a glycoprotein with a molecular mass of about 30 000 Da, that consists of two subunits – alpha and beta. TSH is secreted by frontal lobe of pituitary gland and stimulates thyroxin and triiodothyronine synthesis in thyroid gland. TSH assay is important in evaluation for thyroid status.

Determination of TSH level by third generation kit is used for the differential diagnosis of euthyroidism and hyperthyroidism.

3. PRINCIPLE OF THE TEST

The **TSH (third generation) kit** is a “sandwich” type of solid-phase enzyme immunoassay, based on two monoclonal antibodies that are specific for different epitopes of TSH β -subunit. One of these antibodies is conjugated with horseradish peroxidase (HRP); the other is coated onto the inner surface of microwells. TSH molecules from the serum sample bind to both immobilized antibody and anti-TSH-peroxidase conjugate (Fig. 1). Then the wells are washed with wash solution to remove any material not bound to the inner surface of the wells. Quantity of the bound conjugate is directly proportional to TSH level in the sample.

During the incubation with TMB solution the colour is developing. The intensity of the colour is directly proportional to TSH concentration in specimens. TSH concentration in the patient sample is read from a standard curve that is processed in each assay.


Fig. 1 Assay scheme

4. MATERIALS PROVIDED

MP	Microplate: 12 breakable 8-well strips (total 96 wells) coated with anti-TSH monoclonal antibodies	1 pcs
CONJ	Conjugate: solution contains anti-TSH monoclonal antibodies conjugated with HRP	14 mL, ready to use
0-5 CAL	TSH calibrators: protein-based solution or lyophilized preparations containing known TSH concentrations - 0; 0.02; 0.075; 0.25; 0.75; 2.5 μ U/mL. For exact TSH concentrations see vial labels	6 vials, 0.5 mL each; ready to use or lyophilized
CONTROL	TSH control: protein-based solution or lyophilized preparation containing known TSH concentration. For exact range of TSH concentration see vial label	0.5 mL, ready to use or lyophilized
SUB	Substrate (TMB solution): 3,3',5,5'-tetramethylbenzidine solution in citrate buffer containing hydrogen peroxide	14 mL, ready to use

WASH P 20X	Wash solution P, 20X concentrated: surfactant in buffered saline, sufficient for preparation of 280 mL of solution	14 mL, concentrated
STOP	Stop solution: 1 M HCl solution	14 mL, ready to use

5. EQUIPMENT AND MATERIALS REQUIRED BUT NOT PROVIDED

- 1-channel calibrated variable precision pipettes, with disposable tips;
- 8-channel calibrated variable precision pipette, with disposable tips;
- microplate incubator/shaker (+37 °C, shaking speed 400–800 rpm);
- manual or automatic equipment for rinsing wells;
- calibrated microplate reader (450 nm);
- vortex tube mixer;
- deionized or distilled water;
- graduated beaker and cylinder of appropriate volume;
- latex or plastic gloves;
- trays for pipetting reagents with 8-channel pipette;
- disinfectant;
- absorbent material (for manual wash).

6. STORAGE AND STABILITY OF THE KIT

The expiration date of the kit is printed on the box label; the expiration date for each component is printed on the respective label.

The **TSH (third generation) kit** should be stored at +2...+8 °C upon receipt, preferably in the original kit box, until the expiration date. Storage at +25 °C is allowed but for no more than 5 days.

Shelf life of the kit is 12 months.

After initial opening, the kit is stable until the expiration date if stored at +2...+8 °C.

If used in several separate experiments, after initial opening kit contents should be stored as follows:

- unused strips: in a firmly closed resealable zipper bag at +2...+8 °C until the expiration date;
- vials with conjugate, calibrators and control (ready-to-use): at +2...+8 °C until the expiration date; vials with calibrators and control (reconstituted): at +2...+8 °C for no more than 1 month after opening;
- vial with substrate: at +2...+8 °C until the expiration date, protected from light;
- vials with concentrated wash solution and stop solution: at +2...+8 °C until the expiration date;
- wash solution prepared for use: at room temperature (+18...+25 °C) for no more than 5 days or at +2...+8 °C for maximal 4 weeks in a firmly closed bottle.

Damaged Test Kits

In case of any severe damage of the test kit or components, Astra Biotech GmbH has to be informed in writing, during one week after receiving the kit. Usage of severely damaged single components for a test run is not recommended.

7. SAMPLE COLLECTION AND STORAGE

Collect blood by venipuncture. After clotting, the serum is separated by centrifugation.

Do not use plasma, haemolyzed (bright red) or lipaemic (milky) serum samples as well as samples containing sodium azide as preservative.

Store serum samples at +2...+8 °C for no more than 2 days. Aliquot and freeze the samples for a longer storage (≤ -20 °C). Avoid repeated freezing.

8. TEST PROCEDURE

The **TSH (third generation) kit** is designed for 96 tests. This is sufficient for 40 unknowns, 6 calibrators, 1 control and 1 blank (OD of TMB solution) in duplicates, provided that all the strips are used simultaneously.

8.1. Reagent Preparation

Allow all the reagents to reach room temperature, and then thoroughly stir.

MP Keep the **microplate** at room temperature (+18...+25 °C) for at least 30 minutes before opening the bag. Place required number of strips into strip holder. Place unused strips onto the resealable zipper bag and reseal duly.

CAL CONTROL Calibrators and Control

Liquid calibrators and control are ready to use.

Prepare lyophilized calibrators and control as follows. Gently tap on the vial caps to knock off all the dry matter. Open the vials and carefully place the caps upside down on the clean dry surface. Add 0.5 mL of distilled or deionized water to each vial with lyophilized calibrators and control, close vials with the corresponding caps and leave for 10 min at room temperature without stirring. Then stir gently avoiding foaming, until the dry matter is completely dissolved. Leave for another 10 minutes at room temperature stirring gently periodically. Make sure that no dry matter is left on the caps and walls of the vials.

WASH P Prepare required volume of **wash solution** by dilution of the concentrate 20-fold with distilled or deionized water. For example:

5 mL **WASH P 20X** + 95 mL water

Mix thoroughly, avoid foaming.

SUB Protect **substrate** from direct light.

8.2. Sample Preparation

Allow samples to reach room temperature. Stir samples gently in order to ensure homogeneity.

8.3. Assay Procedure

(See assay scheme section 15)

- A.** Pipette **50 μ L** of **calibrators** **CAL** (0-5), **control** **CONTROL** and **patient's samples in duplicates** into the respective wells; **Leave wells A1-A2 empty for blank!**
- B.** Pipette **100 μ L** of **conjugate** **CONJ** into each well, **except wells A1-A2.**
- C.** Incubate strips for **60 minutes while shaking (500–800 rpm) at +37 °C.**
- D.** Wash 5 times, as described below.
- E.** Pipette **100 μ L substrate** **SUB** into each well (including blank); incubate **at room temperature in the dark for 15-30 minutes** (depending on the colour intensity) or **10 minutes while shaking (500-800 rpm) at +37 °C.**
- F.** Pipette **100 μ L stop solution** **STOP** into each well (including blank) in the same sequence and at the same speed as used for dispensing TMB substrate. Shake for 1–2 min at room temperature.
- G.** Read OD at **450 nm within 20 min.**

Wash Procedure

It is advisable to use an automatic microplate washer set at 5 wash cycles and a volume of 300 μ L of wash solution per well per cycle.

If an automatic washer is not available, the wash procedure can be carried out manually as follows:

- remove the contents of the wells into a container with disinfectant;

- dispense 300 µL of wash solution (section 8.1) into each well, shake the plate carefully for 5–10 sec and remove the contents of the wells; repeat 5 times;
- strike the wells sharply on absorbent material to remove any liquid residue.

9. DATA PROCESSING

If the reader cannot be adjusted to zero using the substrate blank in wells A1-A2, subtract mean OD value of wells A1-A2 from all OD values before further calculations.

Example:

OD (Cal 5) measured = 2.28 and OD (blank) = 0.06;

OD (Cal 5) calculated = $2.28 - 0.06 = \underline{2.22}$

9.1. Data Reliability (for OD measured at 450 nm)

The data should meet the following criteria:

- average blank OD (in wells A1-A2) ≤ 0.100 ;
- average OD of Cal 5 ≥ 1.5 (after blank subtraction);
- control's concentration must fall within the acceptability range that is shown on the vial label.

If the data obtained do not meet the criteria, the results are considered unreliable and the test should be repeated.

9.2. Quantitative Determination

Specialized software for quantitative determination is recommended. Mean OD values of the calibrators at 450 nm versus their respective TSH concentrations using **4PL** or **5PL fit** (Fig. 2 Typical standard curve).

Calculate concentration of TSH in samples using standard curve.

Any extrapolation of the standard curve to TSH concentration above the nominal value of the calibrator 5 (approximately 2.5 $\mu\text{IU/mL}$) is forbidden.


Fig. 2 Example of typical standard curve
Do not use for evaluation of real assay data!

10. EXPECTED VALUES

Level of TSH in serum samples collected from 140 apparently healthy people at the age 21–45 was 0.23–3.4 $\mu\text{IU/mL}$.

These limits should be considered as guidelines only.

It is highly recommended for each laboratory to determine its own reference range of TSH concentrations.

11. QUALITY CONTROL

It is recommended to use control samples according to the state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

12. PERFORMANCE CHARACTERISTICS OF THE ASSAY

12.1. Calibration – Traceability:

The **TSH (third generation)** kit is calibrated against the WHO 2nd International Reference Preparation 80/558.

12.2 Specificity:

According to Product specification provided by the supplier, cross-reactivity of both monoclonal antibodies, used for kit manufacture, with intact hCG; β -hCG; LH and FSH does not exceed 1.0%.

12.3 Analytical Sensitivity:

Analytical sensitivity of **TSH (third generation) assay**, or the lowest detectable concentration that can be distinguished from zero calibrator, is 0.01 $\mu\text{IU/mL}$. It is defined as mean OD of 10 replicates of calibrator 0 plus two standard deviations.

12.4. Measurement Range:

The **TSH (third generation) kit** was validated for measurement of TSH concentration within the concentration diapason (without dilution) of 0.01-2.5 $\mu\text{IU/mL}$.

12.5. Intra- and Inter-Assay Variation:

For **intra-assay CV** determination 8 serum samples were assayed in 9 replicates each. The results are shown below.

Sample	Mean TSH concentration, $\mu\text{IU/mL}$	Intra-assay CV	
		SD	CV, %
1	0.02	0.004	20.9
2	0.13	0.020	15.0
3	0.2	0.022	11.2
4	0.5	0.041	8.3
5	1.2	0.067	5.6
6	1.8	0.058	3.2
7	2.2	0.023	1.0
8	2.3	0.019	0.8

For **inter-assay CV** determination 8 serum samples were assayed 3 times by different operators with 1-week interval.

Each specimen was assayed in 9 replicates. The results are shown below.

Sample	Mean TSH concentration, $\mu\text{IU/mL}$			Inter-assay CV	
	1 assay	2 assay	3 assay	SD	CV, %
1	0.21	0.20	0.25	0.026	12.0
2	0.42	0.38	0.42	0.023	5.7
3	1.04	1.12	1.08	0.040	3.7
4	1.32	1.26	1.30	0.031	2.4
5	1.24	1.28	1.30	0.031	2.4
6	1.68	1.65	1.58	0.051	3.1
7	2.05	2.00	2.14	0.071	3.4
8	2.10	2.12	2.16	0.031	1.4

13. LIMITATION OF THE METHOD

Any clinical diagnosis should not be based on the results of in vitro diagnostic methods alone. To state a diagnosis, the physician is supposed to consider all the available clinical and laboratory findings.

This **TSH (third generation) kit** is specially intended for differential diagnostics of euthyrosis and hyperthyrosis.

If TSH concentration in the sample is higher than in calibrator 5, the sample should be assayed with **TSH kit** (REF 24-01) from Astra Biotech with detectable concentration range of 0.25-15 $\mu\text{IU/mL}$.

14. SAFETY PRECAUTIONS

- **This kit is for in vitro diagnostic use only.** Operator should follow the manual closely in order to ensure reliable data. The manual is valid for the present kit only, within the listed composition. Any substitution of kit components is not allowed by CE regulations.
- Do not use kit or its components after the expiration date indicated on the label. Take into account stability period for reconstituted reagents.
- Do not mix or use together reagents from different lots of the kit except substrate, stop solution and wash solution.
- Do not use substrate, stop solution and wash solution supplied by other vendors.
- Use only “P”-labeled Astra Biotech wash solution.
- Note that stop solution is 1 M HCl solution. Avoid contacts with skin and mucosa. In case of contact rinse affected area thoroughly with plenty of water and seek medical advice.
- Take into account the following common procedural notes:
 - always pipette reagents into wells immediately after washing procedure;
 - avoid contamination of the solutions;
 - in case of partial use of the kit, dispense only required volume of the reagent into the tray;
 - do not pour unused reagents back into the original vials;
 - avoid exposure to direct sunlight during incubations;
 - always pipette reagents in the same order to minimize reaction time differences between wells; the total dispensing time for the calibrators, control and samples must not exceed 15 min;

- the incubation temperature for all the immunological reactions must be kept at +37 °C;
- do not touch the bottom of the wells;
- calibrators should be measured in each separate assay. It is also recommended to measure each time TSH concentration in the control.
- TMB solution should be colourless. Light colouring of solution is admissible. Avoid direct exposure of substrate to sunlight.
-  Source materials of human origin used for kit components preparation were tested and found negative for HBsAg, anti-HIV and anti-HCV antibodies. However, none of known laboratory test guarantees absence of these viral agents. Therefore, all kit components and patient's samples should be handled as potentially hazardous.
-  After usage strips, calibrators, control, specimens and all consumables which contacted with specimens during handling, storage or assay (tubes, vials, gloves, pipette tips etc.) should be collected separately and sterilized by autoclaving. Instead of autoclaving pipette tips may be sterilized by disinfectant treatment. After sterilization all components and expendable materials may be utilized as non-dangerous garbage. Other components of the kit should be discarded into conventional garbage.
-  During manual washing procedure do not discard the contents of the wells directly to drainage. Use a container with disinfectant solution.

-  As the kit contains potentially hazardous material, the following precautions should be taken:
 - do not smoke, eat or drink while performing the assay;
 - always use protective gloves;
 - never pipette material by mouth;
 - in case of spilling, wipe up the spills promptly and wash affected area thoroughly using decontaminant.
- GLP including all general and individual regulations should be applied for the kit usage.

 As the kit contains irritant (**CONJ**, **CAL**, **CONTROL**), the following precautions should be observed:

- P261 - Avoid breathing spray;
 - P272 - Contaminated work clothing should not be allowed out of the workplace;
 - P280 - Wear protective gloves/protective clothing/eye protection;
 - P302+P352 - IF ON SKIN: Wash with plenty of soap and water;
 - P333+P313 - If skin irritation or rash occurs: Get medical advice/attention;
 - P363 - Wash contaminated clothing before reuse;
 - P501 - Dispose of contents/container in accordance with national regulation.
- Precautionary statements according to Regulation EC № 1272/2008.

15. ASSAY SCHEME

Wells	«Blank»	CAL CONTROL	Samples
Reagents			
CAL CONTROL	–	50 µL	–
Samples	–	–	50 µL
CONJ	–	100 µL	100 µL
Incubation №1	60 min, +37 °C, 500–800 rpm		
WASH P (diluted)	5 x 300 µL		
SUB	100 µL	100 µL	100 µL
Incubation №2	15–30 min, +18...+25 °C, in the dark		
	10 min, +37 °C, 500–800 rpm		
STOP	100 µL	100 µL	100 µL
Stirring	1–2 min, +18...+25 °C		
OD measuring	450 nm		
Calculations	Corresponding software		

January, 26, 2016

Astra Biotech GmbH
Rudower Chaussee 29
12489 Berlin, Germany
Telefon: +49 (0)30 74696509
E-Mail: info@astrabiotech.de
www.astrabiotech.de