Multi DNA Extraction kit

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

Getrocknete Blutropfen (DBS)

2. INTENDED USE

Multi DNA Extraction kit is designed for fast and effective DNA extraction from blood serum/plasma, whole blood, dried blood spots (DBS), saliva, urine, urethral swabs, endocervical and vaginal smears for further nucleic acid amplification testing including Real-Time PCR. The Kit is designed for 96 extractions.

3. PRINCIPLE OF THE TEST

Multi DNA Extraction kit is based on guanidine thiocyanate lysis followed by alcohol precipitation.

4. MATERIALS PROVIDED

LB	Lysis buffer	1 bottle - 50 mL
PS	DNA-precipitation solution	2 bottles - 2x 35 mL
WASH	Washing solution	2 bottles - 2x 50 mL
EB	Elution buffer	1 bottle – 10 mL
NC	Negative control	1 microcentrifuge tube - 1.5 mL

Note: *If extracted DNA is intended for the subsequent analysis with Astra Biotech Infectious PCR kits, internal control **IC** contained in these PCR kits are needed for the extraction procedure with the **Multi DNA Extraction kit**.*

5. EQUIPMENT AND MATERIALS REQUIRED BUT NOT SUPPLIED

- *Internal Control **IC** from the respective Infectious PCR diagnostic kit (in case of subsequent analyses with Astra Biotech Infectious PCR kits);*
- Biosafety cabinet;
- Heating block or thermo-shaker (temperature up to +75 °C) suitable for 1.5 mL tubes;
- Vortex mixer;
- Microcentrifuge (max. speed 13,000 rpm) suitable for 1.5 mL tubes;
- Medical aspirator;
- Non-filter pipette tips (up to 200 µL) for use with medical aspirator;
- 2 separate sets of pipettes for DNA-free and DNA-containing components (20-200 µL; 100-1,000 µL);
- Filter tips (100; 200 and 1,000 µL);
- Disposable safe-seal 1.5 mL microcentrifuge tubes;
- DNA-zone-only working area with separate protective laboratory clothing;
- Waste bin with disinfectant;
- 0.9 % sodium chloride solution and transport medium;
- **DBSb** Dried blood spot pretreatment buffer (REF 80-03) and Proteinase K (REF 80-04) for isolation of genomic DNA from dried blood spots.

6. STORAGE CONDITIONS AND STABILITY OF THE KIT

The expiry date of the kit is stated on the box label, expiry date for each component is indicated on the respective label. **Multi DNA Extraction kit** can be stored in the original kit box at +2...+8 °C during the entire shelf life, or at +25 °C for no more than 7 days.

DBSb is stable for 2 years if stored at +2...8 C.

Store Proteinase K, lyo. for 3 years at -18...-22°C.

Prot.K (ready to use) can be stored for one year from the date of preparation at -18...-22°C.

7. SAMPLE COLLECTION AND STORAGE

Type	Conditions	Stability time
Dried blood spots (DBS)	+18...+25 °C	3 years
Whole blood	+2...+8 °C	2 days (quantification assays)
		1 month (genomic DNA extraction only)
Blood serum/ plasma	+2...+8 °C	5 days
	-18...-22 °C	1 year
	-68...-72 °C	Longtime stability

Swabs, smears and urine samples	+18...+25 ° C	2 days
	+2...+8 °C	4 days (2 weeks for urine)
	-18...-22 °C	Longtime stability
Saliva	+18...+25 ° C	6 hours
	+2...+8 °C	1 week
	-18...-22 °C	1 month
	-68...-72 °C	Longtime stability

8. SAMPLE PREPARATION

- **Urine samples**

Shake the urine specimen cup. Transfer 1 mL of urine into new 1.5 mL microcentrifuge tube. Centrifuge at 13,000 rpm for 5 min. Remove the supernatant using medical aspirator. Add 0.9 % NaCl solution to obtain overall volume of 200 µL. Vortex thoroughly until pellet dissolves.

- **Saliva-Mouthwash samples**

Rinse mouth thoroughly with 10 mL of 0.9 % sodium chloride solution during 10-15 s. Collect fluid in sterile 15 mL microcentrifuge tube. Centrifuge at 3,000 rpm for 3 min. Remove approximately 9 mL of the supernatant using medical aspirator. Residual volume should be approximately 0.5-1 mL

(pellet and liquid phase). Add 200 μ L of transport medium. Vortex thoroughly until pellet dissolves.

- **Swabs and scrapes samples**

Collect the swabs with a suitable tool. After specimen collection place swab into a capped 1.5 mL microcentrifuge tube with transport medium and rotate it for 15 s. Avoid splashing the solution. Squeeze out residual liquid from the swab and discard the swab.

- **Blood**

- **Whole blood**

Collect blood into the anticoagulant-treated plastic tube. EDTA and trisodium citrate could be used as the anticoagulant.

Note: *Utilizing heparin as the anticoagulant is not allowed.*

- **Blood plasma**

To obtain blood plasma centrifuge whole blood with anticoagulant at 3,000 rpm for 20 min at room temperature (+18...+25 °C). Transfer the supernatant into 1.5 mL microcentrifuge tube.

- **Blood serum**

Collect whole blood into the anticoagulant-free plastic tube. Allow the whole blood to clot at room temperature for 30 min. Centrifuge at 3,000 rpm for 10 min. Transfer the supernatant into 1.5 mL microcentrifuge tube.

- **Dried blood spots (DBS)**

Apply a few drops of blood, drawn by lancet onto specially manufactured absorbent filter paper. Allow the blood to

thoroughly saturate the paper and dry it on air for several hours. Touching of blood collection card with bare hands is prohibited for everyone except examinee.

Store blood collection cards in zip lock plastic bag with 3-4 desiccant bags and humidity indicator card. Cut three 3 mm (1/8 inch) diameter discs from a dried blood spot with a single-hole paper puncher and place them into new 1.5 mL microcentrifuge tube. Then mark the tubes.

Note: *It is strongly recommended to cut discs immediately before DNA extraction procedure.*

9. ASSAY PROCEDURE FOR DNA EXTRACTION FROM BLOOD SERUM/PLASMA, SALIVA, URINE, SWABS, SMEARS AND WHOLE BLOOD

Note: *In case of no subsequent analysis with Astra Biotech Infectious PCR kits is followed; skip steps D, E and G.*

- A.** Preheat a heating block to +75 °C
- B.** Check the Lysis buffer **LB** for precipitation. If **LB** forms a precipitate, it should be heated until the precipitate dissolves. Shake the bottle briefly and allow cooling to room temperature before use.
- C.** Prepare and label required quantity of 1.5 mL microcentrifuge tubes equal to the number of samples.

Note: *In case of subsequent analyses with Astra Biotech Infectious PCR kits: take one 1.5 mL microcentrifuge tube for negative extraction control and label as “NEC”. **NC** must*

be extracted simultaneously with clinical samples in “NEC” tube.

- D.** Vortex thoroughly Internal control **IC** tube then centrifuge briefly for 2-3 s at 1,500-2,400 rpm.
- E.** Add 20 μL **IC** into each prepared microcentrifuge tube (including “NEC”).
- F.** Add 500 μL **LB** into each prepared microcentrifuge tube.
- G.** Add **NC** to the microcentrifuge tube marked as “NEC”:

Material	NC
Blood serum/ plasma, saliva, urine, swabs, smears	100 μL
Whole blood	25 μL

- H.** Add respective clinical sample to microcentrifuge tubes:

Material	Sample
Blood serum/plasma, saliva, urine, swabs, smears	100 μL
Whole blood	25 μL

- I.** Vortex for 3-5 s.
- J.** Incubate at +75 °C for 15 min. Vortex tube every 3 min for 1-2 s. If thermo-shaker is used, set up 1,300 rpm.
- K.** Set up heating block/ thermo-shaker at +60 °C.
- L.** Cool the tubes at room temperature for 2 min.
- M.** Briefly centrifuge the 1.5 mL microcentrifuge tubes to remove drops from the inside of the lid.
- N.** Add DNA-precipitation solution **PS**. Vortex for 5-10 s.

Material	PS
Blood serum/plasma, saliva, urine, swabs, smears	600 μ L
Whole blood	550 μ L

- O.** Centrifuge at 13,000 rpm for 5 min. Remove the supernatant with medical aspirator. In case of precipitate is barely discernible reduce supernatant volume to the level of 3 mm above tube bottom.
- P.** Add 1 mL **WASH** to the pellet. Vortex tube for 5-10 s.
- Q.** Centrifuge at 13,000 rpm for 3 min. Carefully remove all supernatant with medical aspirator.
- R.** Incubate at room temperature for 3-5 min with tube caps opened.

Note: *Do not allow the pellet to overdry!*

- S.** Add 60 μ L Elution buffer **EB**. Vortex thoroughly.

Note: *Volume of elution could be increased up to 100 μ L if needed.*

- T.** Incubate at +60 °C for 10 min. Vortex tube every 2 min for 5-7 s. If thermo-shaker is used set up 1,300 rpm.
- U.** Centrifuge at 13,000 rpm for 1 min.
- V.** The eluate contains pure genomic DNA. Store DNA at -18...- 22 °C continuously, or at +2...+8 °C for up to 1 week.

Note: *Vortex tube thoroughly and centrifuge it for 1 min at 13,000 rpm before repeated use.*

10. ASSAY PROCEDURE FOR DNA EXTRACTION FROM DRIED BLOOD SPOTS

- A.** Check the Lysis buffer **LB** for precipitation. If **LB** forms a precipitate, it should be heated until the precipitate dissolves. Shake the bottle briefly and allow cooling to room temperature before use.
- B.** Prepare **Prot.K**: add 1 ml of Dilution buffer to the tube with Proteinase K, lyo (10 mg). Mix **Prot.K** and briefly centrifuge to remove drops from the inside of the lid.
- C.** Preheat a heating block/ thermo-shaker to +56°C
- D.** Take the tubes with three 3 mm (1/8 inch) diameter disks from a dried blood spots (see preparation paragraph 8) and add 100 µL **DBSb**.
- E.** Add 10 µL **Prot.K**, mix by vortexing, and incubate at 56°C for 60 min, vortex tubes every 10 min for 1-2 s. Briefly centrifuge to remove drops from inside the lid.
- F.** Set up heating block/ thermo-shaker at +75 °C.
- G.** Add 500 µL **LB**.
- H.** Vortex thoroughly for 3-5 s.
- I.** Incubate at +75 °C for 15 min. Vortex tube every 3 min for 1-2 s.
- J.** Set up heating block/ thermo-shaker at +60 °C.
- K.** Centrifuge at 13,000 rpm for 3 min, then transfer supernatant into new 1.5 mL microcentrifuge tube.
- L.** Cool the tubes at room temperature for 2 min.
- M.** Briefly centrifuge the 1.5 mL microcentrifuge tubes to remove drops from the inside of the lid.
- N.** Add 600 µL DNA-precipitation solution **PS**. Vortex for 5-10s.
- O.** Centrifuge at 13,000 rpm for 5 min. Remove the supernatant with medical aspirator. In case of precipitate is barely

discernible reduce supernatant volume to the level of 3 mm above tube bottom.

P. Add 1 mL **WASH** to the pellet. Vortex tube for 5-10 s.

Q. Centrifuge at 13,000 rpm for 3 min. Carefully remove all supernatant with medical aspirator.

R. Incubate at room temperature for 3-5 min with tube caps opened.

Note: *Do not allow the pellet to overdry!*

S. Add 60 μ L Elution buffer **EB**. Vortex thoroughly.

Note: Volume of elution could be increased up to 100 μ L if needed.

T. Incubate at +60°C for 10 min. Vortex tube every 2 min for 5-7 s. If thermo-shaker is used set up 1,300 rpm.

U. Centrifuge at 13,000 rpm for 1 min.

V. The eluate contains pure genomic DNA. Store DNA at -18...-22 °C continuously, or at +2...8 °C for up to 1 week.

Note: *Vortex tube thoroughly and centrifuge it for 1 min at 13,000 rpm before repeated use.*

11. TROUBLESHOOTING

1. Low yield of DNA

- Volume of added specimen is too large. Do not exceed volumes designated in this manual.
- Specimen may be old or degraded. Store specimens appropriately or use fresh specimens.
- Specimen inefficiently disrupted. Warm Lysis buffer until the precipitate dissolves completely.

- Pellet is dissolved incompletely. Do not overdry pellets.
2. Downstream applications are inhibited
- Residual Washing buffer in purified DNA: Carefully remove all Washing buffer as much as possible then dry pellets as consisted with procedure.
3. DNA is degraded
- Sample was undergone repeated frozen-thaw cycles. Avoid repeated frozen-thaw cycles, use fresh preparations immediately.
 - Inappropriate storage conditions. Store samples and components as consisted with protocols.

12. PERFORMANCE CHARACTERISTICS OF THE ASSAY

Expected yields of genomic DNA will vary depending on the amount and type of starting material used. The efficiency of extraction is not less than 80 %.

13. QUALITY CONTROL

The quality of obtained DNA with the kit can be assessed:

- by analysis of internal control in case of subsequent analyses with Astra Biotech Infectious PCR kits;
- by DNA gel electrophoresis on a 1 % agarose.

14. SAFETY PRECAUTIONS

- The operator should thoroughly follow the manual to obtain the reliable data.

- Lysis buffer **LB** is harmful and irritant, the following precautions should be observed:
 - P261 - Avoid breathing spray;
 - P272 - Contaminated work clothing should not be allowed out of the workplace;
 - P280 - Wear protective gloves/protective clothing/eye protection;
 - P363 - Wash contaminated clothing before reuse;
 - P501 - Dispose of contents/container in accordance with national regulation;
 - P302+P352 - IF ON SKIN: Wash with plenty of soap and water;
 - P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

- DNA-precipitation solution **PS** and Washing solution **WASH** are highly flammable and can cause serious eye damage: ,
 - P210 - Keep away from heat/ sparks/ open flames/ hot surfaces.-No smoking;
 - P233 - Keep container tightly closed;
 - P240 - Ground/ bond container and receiving equipment.
 - P241 - Use explosion-proof electrical/ ventilating/light equipment;
 - P242 - Use only non-sparking tools;
 - P243 - Take precautionary measures against static discharge;

- P280 - Wear protective gloves/protective clothing/eye protection/face protection;
- P310 - Immediately call a POISON CENTER or doctor/physician;
- P305+351+338- IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do – continue rinsing.

- Elution buffer **EB** and Negative control **NC** contain 0.05 % of sodium azide which is very toxic:
 - P260 - Do not breathe dust/fume/gas/mist/vapours/spray;
 - P262 - Do not get in eyes, on skin, or on clothing;
 - P264 – Wash hands thoroughly after handling;
 - P270 - Do not eat, drink or smoke when using this product;
 - P271 - Use only outdoors or in a well-ventilated area;
 - P273 - Avoid release to the environment;
 - P280 - Wear protective gloves/protective clothing/eye protection/face protection;
 - P284 - Wear respiratory protection;
 - P310 - Immediately call a POISON CENTER or doctor/physician;
 - P330 - Rinse mouth;
 - P361 - Remove/Take off immediately all contaminated clothing;
 - P363 - Wash contaminated clothing before reuse;
 - P391 - Collect spillage;
 - P405 - Store locked up;
 - P501 - Dispose of contents/container in accordance with national regulation;

- P301 +P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician;
- P302 +P350 - IF ON SKIN: Gently wash with plenty of soap and water;
- P304 +P340 - IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing;
- P403 +P233 - Store in a well-ventilated place. Keep container tightly closed.

If assay procedure for DNA extraction from dried blood spots is used, the following precautions should be observed:

- P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P342 + P311 If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician.


Precautionary statements according to Regulation EC № 1272/2008.

-  When using Multi DNA Extraction kit for DNA extraction of clinical materials, **biosafety level 2 practices and procedures** have to be followed. Human blood, blood products, body fluids and tissues should be treated as potential infectious.
- To avoid the contamination:
 - use two separate pipette sets for DNA-free and DNA-containing components,
 - use filter tips,


- do not use the same tip for two different components, neither DNA-free nor DNA-containing.
- Do not use the kit after its expiration date.
- Do not pool reagents from different lots or from different bottles of the same lot. Immediately after use, close all bottles in order to avoid leakage. After opening, store all bottles and vials in an upright position.
- In case of using several Astra Biotech Infectious Real-Time PCR kits be sure of using **IC** from latest lot.

15. ASSAY SCHEME

Swabs, smears, whole blood, blood serum/ plasma, urine and saliva samples

Dried blood spots

June, 30, 2015

Astra Biotech GmbH
Rudower Chaussee 29
12489 Berlin, Germany
Telefon: +49 (0)30 74696509
E-Mail: info@astrabiotech.de
www.astrabiotech.de