# ***Neisseria gonorrhoeae* PCR kit**

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## 2. INTENDED USE

The ***Neisseria gonorrhoeae* PCR kit** is provided for fast and qualitative detection of *Neisseria gonorrhoeae* DNA in samples extracted from urethral swab specimens, endocervical and vaginal smears, urine samples, prostatic fluid specimens and other human body fluids and tissue samples by means of Real-Time PCR for diagnosis of *Neisseria gonorrhoeae* presence.

## 3. PRINCIPLE OF TEST

The detection of the ***Neisseria gonorrhoeae* DNA** includes three stages: DNA extraction, PCR amplification and data analysis.

At first stage DNA extraction from specimen with addition of Internal Control **IC** is performed. **IC** is provided with ***Neisseria gonorrhoeae* PCR kit** and is added to each sample controlling both DNA extraction and PCR amplification stages.

The second stage is the Real-Time PCR to detect the DNA of *Neisseria gonorrhoeae*. During the Real-Time PCR the fragments of *Neisseria gonorrhoeae* genomic DNA and **IC** DNA are both amplified and detected simultaneously. During PCR increases quantity of amplified fragments of *Neisseria gonorrhoeae* DNA and **IC** DNA which hybridize to fluorescent reporter probes.

The exonuclease activity of the Taq DNA-polymerase **Taq** degrades oligonucleotide probes and separates the fluorophore from the quencher, thus cancelling quenching effect and allowing emission of fluorescence detected and measured by Real-Time PCR cyclers.

Therefore, resulting quantification cycle C<sub>q</sub> (threshold cycle (C<sub>t</sub>) or crossing point (C<sub>p</sub>)) is used to determine the presence of target DNA in the sample.

The “hot start” Taq DNA-polymerase **Taq** improves sensitivity and specificity of the test.

## 4. MATERIALS PROVIDED

**4.1 Form S** (liquid reaction mix aliquoted into 0.2 mL low profile 8-tube strips and separate DNA polymerase)

Label	Component	Volume, $\mu\text{L}$	Quantity
<b>Ng-Mix-S</b>	PCR-Mix <i>N. gonorrhoeae</i> , Form S	10	14x 8 ×0.2 mL tubes
<b>Taq</b>	Taq DNA polymerase	1120	1 tube
<b>PC</b>	Positive Control	250	1 tube
<b>IC</b>	Internal Control	1120	2 tubes
<b>NC</b>	Negative Control	1500	1 tube

The format is compatible with the following amplifiers: CFX96/CFX96 Touch, SaCycler, Light Cycler 96 and other Real-Time PCR cyclers compatible with low profile PCR tubes.

#### 4.2 Form T (liquid reaction mix aliquoted into 0.2 mL regular profile reaction tubes and separate DNA polymerase)

Label	Component	Volume, $\mu\text{L}$	Quantity
<b>Ng-Mix-T</b>	PCR-Mix <i>N. gonorrhoeae</i> , Form T	10	112× 0.2 mL tubes
<b>Taq</b>	Taq DNA - polymerase	1120	1 tube
<b>PC</b>	Positive Control	250	1 tube
<b>IC</b>	Internal Control	1120	2 tubes
<b>NC</b>	Negative Control	1500	1 tube

The format is compatible with the following amplifiers: Rotor-Gene Q/ 3000/ 6000 and other Real-Time PCR cyclers compatible with regular profile PCR tubes.

### 5. MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

- Biosafety cabinet;
- Real-Time PCR cyclers;
- Vortex mixer;
- Microcentrifuge (suitable for 0.2 mL tubes, max. 2400 rpm);
- 2 separate pipette sets for DNA-free and DNA-containing components;
- Filter tips (10, 100; 200 and 1000  $\mu\text{L}$ );
- Disposable safe-seal 1.5 mL microcentrifuge tubes;
- Separate working areas for DNA extraction, PCR preparation and PCR amplification each with set of pipettes, instruments and laboratory clothing.

## 6. STORAGE CONDITIONS AND STABILITY OF THE KIT

The kit is designed to perform 112 reactions, including controls.

***Neisseria gonorrhoeae* PCR kit** should be stored at -18...-30°C, preferably in the original kit box until the expiration date. The expiry date of the kit is stated on the box label, expiry date for each component is indicated on the respective label.

Storage at +25 °C is allowed for no more than 7 days (any kit format).

## 7. SAMPLES

Total DNA can be extracted from urethral swab specimens, endocervical and vaginal smears, urine samples, prostatic fluid specimens and other human body fluids and tissue samples. All methods of extraction can be used.

## 8. ASSAY PROCEDURE

Positive and Negative Controls must be included in each run.

The complete procedure consists of three stages:

1. DNA Extraction
2. Real-Time PCR
3. Data analysis and interpretation

## 1. DNA extraction

**NOTE:** Application of Astra Biotech DNA extraction kits for specimen preparation is strongly recommended.

- 1.1 Take **IC** to thaw at room temperature and place the rest components back in freezer. Gently vortex **IC** and centrifuge for 2-3 s.
- 1.2 Prepare and mark required quantity of 1.5 mL microcentrifuge tubes equal to the number of samples + 1 tube for negative extraction control and label as “NEC”.

**NOTE:** Each extraction procedure should be performed using “NEC”.

- 1.3 If Astra Biotech DNA extraction kit is used, ignore steps 1.4-1.6 and follow the Astra Biotech DNA extraction kit manual; otherwise continue with step 1.4.
- 1.4 Place **20 µL IC** into each 1.5 mL microcentrifuge tube.
- 1.5 Add **100 µL NC** into the microcentrifuge tube marked “NEC”.

**NOTE:** The Elution buffer contained in DNA extraction kit could be used instead of **NC**.

- 1.6 Follow the manual of respective DNA extraction kit.

**NOTE:** Elution volume should not exceed 200 µL.

## 2. Real-Time PCR

### 2.1 Kit components preparation:

- Take the kit out from the refrigerator, thaw all liquid components.
- Gently vortex all liquid components and centrifuge briefly.

### 2.2 Microcentrifuge tubes preparation:

- Prepare and mark required quantity of 0.2 mL microcentrifuge tubes/ wells **Ng-Mix-S/T** equal to the number of samples + 1 tube/ well for “PC” + 1 tube/ well for “NEC”.

**NOTE:** *Different positions for marking the tubes are necessary because of differences in the fluorescence detection: When using CFX96/ 96 Touch, LC 96 or SaCycler Real-Time PCR cyclers mark tubes outside the centre of tube cap or on the tube side. When using Rotor-Gene Q/ 3000/ 6000 instruments mark the tubes on the caps.*

- Place **10 µl Taq** into each tube avoiding tip contact with paraffin.
- Add **10 µl DNA** (extracted from specimens and “NEC”) avoiding tip contact with paraffin.
- Add **10 µL positive control PC** in the PCR tube marked “PC”. Use separate filter tip for each sample.
- Close the tubes.
- Vortex, then centrifuge at 1500-2400 rpm for 2-3 s.

### 2.3 Real-Time PCR Amplification:

- Place tubes/ plate into the Real-Time PCR cycler.
- Fill in the location of the samples. Set up fluorescent detection from channels ROX/ Orange and FAM/ Green.
- Run the PCR protocol.

Stage	Temperature, °C	Data collection	Time	Number of cycles
Hold	94	-	3 min	1
Cycling	94	-	10 s	5
	60		20 s	
Cycling	94	-	10 s	45
	60	ROX/ Orange FAM/ Green	20 s	

## 3. Data Analysis

### 3.1 General conditions:

ROX dye is used to detect the *Neisseria gonorrhoeae* DNA and FAM is used to detect **IC** DNA.

Analysis of the results can be performed by means of a Real-Time PCR thermal cycler's software in accordance to its instruction. For the validated cyclers detailed instructions concerning protocol set up and results analysis are given in « Astra Biotech: Real-Time PCR kits Guidelines».

Analysis of the results can be performed only by using following parameters: Cq **Ng**, Cq **PC** and Cq **IC**, specified in the Quality control Sheet enclosed in every kit.

Cq **Ng** critical value of quantification cycle for *Neisseria gonorrhoeae* in the ROX channel, crucial for detection positive results;

Cq **PC** critical value of quantification cycle for *Positive Control* (PC) in the ROX channel is used to determine whether kit components work properly (general control of the reagents effectiveness), crucial for detection negative results

Cq **IC** critical value of quantification cycle for *Internal Control* (IC) in the FAM channel is used to control all stages of the analysis (overall effectiveness) in each tube, including inhibition and DNA loss.

### **The results of single run are valid if:**

Cq PC  $\leq$  Cq **PC** in the ROX channel.

Nevertheless, if Cq value for PC is determined as  $>$  Cq **PC**, in that case all the specimens defined as positive in this run should be considered as valid, whereas all the specimens defined as negative should be considered as invalid.

Cq NEC  $>$  Cq **Ng** in the ROX channel and

Cq NEC  $\leq$  Cq **IC** in the FAM channel,

but there are some exclusions, for more information see 3.2.

### 3.2 Data interpretation

#### **Specimen is considered as positive on DNA of *Neisseria g.*:**

$Cq \leq Cq_{Ng}$  in the ROX channel,

even if  $Cq > Cq_{IC}$  in the FAM channel or is not determined.

#### **Specimen is considered as negative on DNA *Neisseria g.*:**

- $Cq > Cq_{Ng}$  in the ROX channel or is not determined
- $Cq \leq Cq_{IC}$  in the FAM channel
- $Cq_{PC} \leq Cq_{PC}$  in the ROX channel

#### **Specimen is considered as invalid:**

$Cq > Cq_{Ng}$  in the ROX channel and at least one of this conditions are satisfied:

- $Cq > Cq_{IC}$  in the FAM channel or is not determined.
- $Cq_{PC} > Cq_{PC}$  in the ROX channel or is not determined.

#### **In case of contamination:**

$Cq_{NEC} \leq Cq_{Ng}$  in the ROX channel

Measures must be taken to eradicate the source of contamination. In case of contamination, all positive specimens satisfying condition  $Cq \geq (Cq_{NEC} - 5)$  in the ROX channel are considered as invalid, the other results are considered as valid. For all specimens with invalid results It is recommended to perform the analysis from the DNA extraction stage.

Data interpretation, summary table:

Cq ROX Detection of <i>N.g.</i> <b>DNA for sample</b>	Cq FAM Detection of IC DNA <b>for sample</b>	Cq <b>ROX</b> Detection of <i>N.g.</i> DNA <b>for PC</b>	Cq <b>ROX</b> Detection of <i>N.g.</i> DNA <b>for NEC</b>	Result
$Cq \leq Cq_{Ng}$	any	any	$Cq > Cq_{Ng}$ or N/A	Positive
$Cq > Cq_{Ng}$ or N/A	$Cq \leq Cq_{IC}$	$Cq \leq Cq_{PC}$	any	Negative
$Cq > Cq_{Ng}$ or N/A	$Cq > Cq_{IC}$ or N/A	any	any	Invalid
$Cq > Cq_{Ng}$ or N/A	any	$Cq > Cq_{PC}$ or N/A	any	
$Cq \leq Cq_{Ng}$ and $Cq > Cq_{NEC} - 5$	any	any	$Cq \leq Cq_{Ng}$	Positive
$Cq \leq Cq_{Ng}$ and $Cq \leq C_{NEC} - 5$	any	any	$Cq \leq Cq_{Ng}$	

$Cq_{Ng}$  = *critical value of quantification cycle for Neisseria gonorrhoeae*

$Cq_{IC}$  = *critical value of quantification cycle for Internal Control (IC)*

$Cq_{PC}$  = *critical value of quantification cycle for Positive Control (PC)*

NEC = *negative extraction control*

$Cq_{NEC}$  = *quantification cycle for NEC in the ROX channel is determined in case of contamination*

N/A = *Cq value is not determined*

## 9. PERFORMANCE CHARACTERISTICS OF THE KIT

### Analytical sensitivity:

The detection limit of the kit evaluated via probit analysis is 180 copies of *Neisseria gonorrhoeae* genomic DNA per 1 mL (3 copies of genomic DNA per reaction). This limit is defined for confidence level of 95%. It is assumed that 100 µl of specimen was sampled for DNA extraction and elution volume was 60 µl.

### Analytical specificity:

A false-positive reaction is absent with standard human DNA samples and standard DNA samples of *Ureaplasma urealyticum*, *Ureaplasma parvum*, *Mycoplasma hominis*, *Mycoplasma genitalium*, *Gardnerella vaginalis*, *Lactobacillus spp*, *Escherichia coli*, *Staphylococcus aureus*, *Streptococcus pyogenes*, *Streptococcus agalactiae*, *Trichomonas vaginalis*, *Chlamydia trachomatis*, *Chlamydophila pneumoniae*, *Chlamydophila psittaci*, *Neisseria flava*, *Neisseria subflava*, *Neisseria mucosa*, *Neisseria sicca*, *Treponema pallidum*, *Toxoplasma gondii*, *Candida albicans*, *HSV type 1 and 2*, *CMV* and *HPV*.

## 10. LIMITATIONS OF THE METHOD

Any clinical diagnosis should not be based on the results of *in vitro* diagnostic methods alone. For diagnosis establishment a physician is supposed to consider all the available clinical and laboratory findings.

## 11. SAFETY PRECAUTIONS

- **This kit is for in vitro diagnostic use only.** The operator should thoroughly follow the manual to obtain reliable data. This instruction manual is valid only for the present kit with the listed contents. Any exchange of the kit components is not allowed by CE regulations.
- Test should only be performed by skilled personnel considering GLP (Good Laboratory Practice) guidelines.
- Don't use the kit after its expiration date.
- Do not pool reagents from different lots or from different tubes of the same lot (the only exception for **IC** is described below). Immediately after use, close all tubes in order to avoid leakage.
- In case of simultaneous using of two or more Astra Biotech Infectious Real-Time PCR kits for detection of different infectious agents, **IC** with the latest lot number can be used for all of these kits.
- We recommend setting up blank protocols of PCR before the analysis.
- PCR technology is extremely sensitive. The amplification of a single DNA molecule generates millions of identical copies. Therefore, set up three separate working areas for a) sample preparation b) PCR reagent preparation and c) PCR amplification. For each working area a different set of pipettes and protective clothing should be reserved.
- Use sterile filter tips for pipetting and use special PCR pipettes for aerosol free pipetting.
- Do not use the same tip for two different components, neither DNA-free nor DNA containing.

- Do not eat, drink or smoke in the laboratory work area.
- Be careful while handling Positive Control, avoid its splashing.
- Kits contaminated with *Neisseria gonorrhoeae* DNA/ **PC** should be discarded in a biological waste container.
- Routinely decontaminate your pipettes and the laboratory benches.
- After usage, all the reagents and test components should be discarded as consisted with local statements.

September 24, 2018

Astra Biotech GmbH  
Rudower Chaussee 29  
12489 Berlin, Germany  
Telefon: +49 (0)30 74 696 509  
E-Mail: [info@astrabiotech.de](mailto:info@astrabiotech.de)  
[www.astrabiotech.de](http://www.astrabiotech.de)