Rapid DNA Extraction kit

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2. INTENDED USE

Rapid DNA Extraction kit is designed to provide fast and effective DNA extraction from saliva, urine, spinal fluid, prostatic fluid, swabs and scrapes from genitourinary system, oropharynx, conjunctiva, erosive and ulcerative lesions for further biomolecular applications including Real Time PCR.

Kit is designed for 100 extractions and for *in vitro* diagnostic use only.

3. PRINCIPLE OF THE TEST

Rapid DNA Extraction kit is based on rapid thermal lysis. Clinical specimen is placed into Transport Medium **TM-E** wherein it could be transported and stored. Sample aliquot in **TM-E** is transferred into Lysis Buffer **LB**. Afterwards DNA is released due to destruction of cell membranes and various biopolymer complexes during heat treatment. Subsequent centrifugation causes precipitation of insoluble components. Hence, the supernatant contains ready-to-PCR DNA.

4. MATERIALS PROVIDED

LB	Lysis Buffer	100 tubes x 300 µL
TM-E	Transport Medium	100 tubes x 500 µL
NC	Negative Control	1 tube - 1.5 mL

Note: *If extracted DNA is intended for the subsequent analysis with Astra Biotech Infectious PCR kits, Internal controls **IC** contained in these PCR kits are needed for the extraction procedure with the Rapid DNA Extraction kit.*

5. MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

- Biosafety cabinet;
- Internal Control **IC** from the respective PCR diagnostic kit (*optional*);
- Heating block or thermo-shaker (temperature up to +100 °C) suitable for 1.5 mL tubes;
- Vortex mixer (maximum speed 2,400 rpm);
- Microcentrifuge (maximum speed 13,000 rpm) suitable for 1.5 mL tubes;
- Medical aspirator (for sample preparation only);
- Non-filter pipette tips (up to 200 µL) for use with medical aspirator;
- 2 separate pipette sets for DNA-free and DNA-containing components;
- Filter tips (200 and 1,000 µL);
- 1.5 mL microcentrifuge tubes;
- 15 mL conical centrifuge tubes;
- DNA-zone-only working area with separate laboratory coat, hat, shoes and disposable powder-free gloves;
- Waste bin with disinfectant for used tips.
- 0.9 % sodium chloride solution (for mouthwash sample preparation only).

6. STORAGE CONDITIONS AND STABILITY OF THE KIT

The expiry date of the kit is stated on the box label, expiry date for each component is indicated on the respective label. **Rapid DNA Extraction kit** can be stored in the original kit box at +2...+8 °C during the entire shelf life or at +25 °C for no more than 7 days.

7. SAMPLE COLLECTION AND STORAGE

Specimen	Temperature range	Stability
Swabs and scrapes from mucous membranes, erosive and ulcerative lesions	+18...+25 °C	48 hours
	+2...+8 °C	4 days
	-18...-22 °C	Longtime
Urine Spinal fluid	+18...+25 °C	48 hours
	+2...+8 °C	2 weeks
	-18...-22 °C	Longtime
Saliva	+18...+25 °C	6 hours
	+2...+8 °C	1 week
	-18...-22 °C	1 year
	-68...-72 °C	Longtime
Prostatic fluid	+18...+25 °C	6 hours
	+2...+8 °C	1 day
	-18...-22 °C	1 week
	-68...-72 °C	Longtime

8. SAMPLE PREPARATION

Clinical specimen collection should be performed as consisted with local or state department of health regulatory guidelines.

*NOTE: If there is no special instruction, clinical specimens should be collected into **TM-E** tubes provided with kit.*

Clinical specimens transportation conditions:

Temperature range	Stability
+18...+25 °C	48 hours
+2...+8 °C	14 days
-20 °C and below	longtime

If extracted DNA is intended for the subsequent analysis with Astra Biotech Infectious PCR kits, only one freeze-thaw cycle is allowed.

- **Swabs and scrapes from mucous membranes, erosive and ulcerative lesions**

Clinical specimen should be placed into **TM-E** tube. Use disposable medical tools only.

- **Urine**

Shake the urine specimen cup. Transfer 1 mL of urine into new 1.5 mL microcentrifuge tube. Centrifuge at 13,000 rpm for 5 min. Remove the supernatant using medical aspirator. Add 500 µL of **TM-E** solution to pellet. Vortex thoroughly.

- **Spinal fluid**

Collect at least 1 mL of spinal fluid into new 1.5 mL microcentrifuge tube using disposable needle. There is no special preparation needed.

In case of some infectious diseases such as WNV, TBEV, leptospirosis and Lyme disease, concentrate virus-containing cells by centrifuging at 10,000 – 11,000 rpm for 5 min. Reduce the supernatant volume to 100 µL approximately using medical aspirator. Vortex residual supernatant and pellet thoroughly.

- **Saliva**

Mouthwash sampling. Rinse mouth thoroughly with 8-10 mL of 0.9 % sodium chloride solution in 10-15 s. Collect fluid in sterile 15 mL conical centrifuge tube. Centrifuge at 3,000 rpm for 3 min. Reduce the supernatant volume to 0.5-1 mL approximately using medical aspirator. Add 200 µL of Transport Medium **TM-E**. Vortex thoroughly until pellet dissolves.

- **Buccal swabs sampling.**

Clinical specimen should be placed into **TM-E** tube. Use disposable medical tools only.

- **Prostatic fluid**

Collect specimen into new 1.5 mL microcentrifuge tube. Add 100 µL of **TM-E**. Vortex thoroughly.

9. ASSAY PROCEDURE

Preparation Instructions

- A. Preheat a heating block to +95 °C
- B. Take **Lysis Buffer** **LB** and **Negative Control** **NC** out from the refrigerator. Incubate them at room temperature for 15 min.
- C. Prepare and mark an appropriate quantity of **Lysis Buffer** **LB** 1.5 mL microcentrifuge tubes equal to the number of samples.

NOTE: *In case of subsequent analyses with Astra Biotech Real Time PCR diagnostic kits take one 1.5 mL microcentrifuge tube for negative extraction control and label as “NEC”. **NC** must be extracted simultaneously with clinical samples.*

Extraction procedure

NOTE: *If no Real Time PCR-analysis using Astra Biotech Real Time PCR diagnostic kit is followed, ignore steps A-C.*

- A. Briefly centrifuge the **IC**, **NC** and **LB** tubes to remove drops from the inside of the lid.
- B. Add 20 µL of **Internal control** **IC** into each prepared microcentrifuge tube.
- C. Add 100 µL of **Negative Control** **NC** to the microcentrifuge tube marked as «NEC».
- D. Add 100 µL of prepared samples (paragraph **Fehler! Verweisquelle konnte nicht gefunden werden.**) to the **LB** microcentrifuge tubes.
- E. Close the tubes tightly. Vortex tubes for 3-5 s.

NOTE: To avoid opening during the heat treatment check on tubes are closed tight!

- F. Incubate at +95 °C for 5 min.
- G. Cool the tubes at room temperature for 2 min.
- H. Centrifuge at 13,000 rpm for 1 min.
- I. The supernatant contains ready to molecular biology applications (including Real Time PCR) total DNA. It can be stored at -18...-22 °C for one year or at 2-8 for one week. If repeated use is needed vortex tubes for 3-5 s and centrifuge them at 13,000 rpm for 1 min.

10. TROUBLESHOOTING

Problem	Cause and solution
PCR application is inhibited	The amount of loaded specimen is too high. Do not exceed volumes designated in this manual.
DNA is degraded	<ol style="list-style-type: none"> 1. The sample undergone repeated freeze-thaw cycles. Avoid repeated freeze-thaw cycles, use fresh preparations immediately. 2. Inappropriate storage conditions. Store samples and components as consistent with instructions.

11. PERFORMANCE CHARACTERISTICS OF THE ASSAY

Expected yields of genomic DNA will vary depending on the amount and type of starting material used. The efficiency of extraction is not less than 80 %.

12. QUALITY CONTROL

The quality of DNA obtained with the kit can be assessed PCR using an internal control.

13. SAFETY PRECAUTIONS

- The operator should thoroughly follow the manual to obtain the reliable data.

- Components **LB**, **TM-E**, **NC** of **Rapid-DNA Extraction kit** contain 0.05% of sodium azide which is very toxic:

- P260 - Do not breathe dust/ fume/ gas/ mist/ vapours/ spray;
- P262 - Do not get in eyes, on skin, or on clothing;
- P264 – Wash hands thoroughly after handling;
- P270 - Do not eat, drink or smoke when using this product;
- P271 - Use only outdoors or in a well-ventilated area;
- P280 - Wear protective gloves/ protective clothing/ eye protection/ face protection;
- P284 - Wear respiratory protection;
- P310 - Immediately call a POISON CENTER or doctor/physician;

- P361 - Remove/Take off immediately all contaminated clothing;
- P363 - Wash contaminated clothing before reuse;
- P301 +P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician;
- P302 +P350 - IF ON SKIN: Gently wash with plenty of soap and water;
- P304 +P340 - IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing;
- P405 - Store locked up;
- P403 +P233 - Store in a well-ventilated place. Keep container tightly closed;
- P 273 - Avoid release to the environment;
- P 391 - Collect spillage;
- P501 - Dispose of contents/container in accordance with national regulation.

Precautionary statements according to Regulation EC No 1272/2008.

-  When using Rapid DNA Extraction Kit for DNA extraction of clinical materials, **biosafety level 2 practices and procedures** have to be followed. Human blood, blood products, body fluids and tissues should be treated as potentially infectious.
- To avoid the contamination:
 - use two separate pipette sets for DNA-free and DNA-containing components,
 - use filter tips,
 - do not use the same tip for two different components, neither DNA-free nor DNA- containing.

- Do not use the kit after its expiration date.
- In case of using several Astra Biotech Real Time PCR diagnostic kits be sure of using **IC** from latest lot.
- In case of using third-party vendor's diagnostical PCR kits be sure of using internal control compatible with actual kit.

March, 19, 2019

Astra Biotech GmbH
Rudower Chaussee 29
12489 Berlin, Germany
Telefon: +49 (0)30 74 696 509
E-Mail: info@astrabiotech.de
www.astrabiotech.de