***Trichomonas vaginalis* PCR kit**

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2. INTENDED USE

The ***Trichomonas vaginalis* PCR kit** is provided for fast and qualitative detection of *Trichomonas vaginalis* (TV) DNA in urethral swab specimens, endocervical and vaginal smears, urine samples, prostatic fluid specimens and other human body fluids and tissue samples by means of Real-Time PCR for diagnosis of trichomoniasis.

3. PRINCIPLE OF TEST

The performance of the ***Trichomonas vaginalis* PCR kit** includes two stages: DNA extraction and PCR amplification.

At first DNA extraction from specimen with addition of Internal Control **IC** is performed. **IC** is provided with ***Trichomonas vaginalis* PCR kit** and is added to each reaction controlling both DNA extraction and PCR amplification stages.

The second stage is the Real-Time PCR to detect the DNA of *Trichomonas vaginalis*. During the Real-Time PCR the fragments of *Trichomonas vaginalis* genomic DNA and **IC** DNA are both amplified and detected simultaneously. The amplification causes the increase of hybridization of fluorescent reporter probes with *Trichomonas vaginalis* DNA and **IC** DNA amplicons.

The exonuclease activity of the Taq DNA-polymerase degrades oligonucleotide probes and separates the fluorophore from the quencher, thus cancelling quenching effect and allowing emission of fluorescence detected and measured by Real-Time PCR machine. Therefore, resulting quantification cycle, C_q, (threshold cycle (C_t) or crossing point (C_p)) is used to determine the presence of target DNA in the sample.

The “hot start” Taq DNA-polymerase **Taq** improves sensitivity and specificity of the test.

4. MATERIALS PROVIDED

4.1 Format S (aliquoted into 0.2 mL 8-tube strips liquid reaction mix and separate DNA polymerase)

Label	Component	Volume, μL	Quantity
S PCR-TV	PCR-Mix S-TV	10	14x8x0.2 mL tubes
Taq	Taq DNA polymerase	1120	1 tube
PC	Positive Control	250	1 tube
IC	Internal Control	1120	2 tubes
NC	Negative Control	1500	1 tube

The format is compatible with the following amplifiers: CFX96/CFX90 Touch and SaCycle, Light Cycler 96 and other Real-Time PCR cyclers that are compatible with low profile tubes.

4.2 Format T (aliquoted into 0.2 mL reaction tubes liquid reaction mix and separate DNA polymerase)

Label	Component	Volume, μL	Quantity
T PCR-TV	PCR-Mix T-TV	10	112x0.2 mL tubes
Taq	Taq DNA polymerase	1120	1 tube
PC	Positive Control	250	1 tube
IC	Internal Control	1120	2 tubes
NC	Negative Control	1500	1 tube

The format is compatible with the following amplifiers: Rotor-Gene 3000/6000 and Rotor-Gene Q and other Real-Time PCR cyclers that are compatible with regular profile tubes.

5. MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

- Biosafety cabinet;
- Real-Time PCR cycler;
- Vortex mixer;
- Microcentrifuge (suitable for 0.2 mL tubes, max. 2400 rpm);
- 2 separate pipette sets for DNA-free and DNA-containing components;
- Filter tips (10, 100; 200 and 1000 μ L);
- Sterile optically clear 0.2 mL microcentrifuge tubes/ strip tubes/ PCR plates suitable for Real-Time PCR system used;
- Disposable safe-seal 1.5 mL microcentrifuge tubes
- Separate working areas for DNA extraction, PCR preparation and PCR amplification each with set of pipettes, instruments and laboratory clothing.

6. STORAGE CONDITIONS AND STABILITY OF THE KIT

The kit is designed to perform 112 reactions, including controls.

***Trichomonas vaginalis* PCR kit** should be stored at +2...+8°C, preferably in the original kit box until the expiration date.

The expiry date of the kit is stated on the box label, expiry date for each component is indicated on the respective label.

Storage at +25 °C is allowed for no more than 7 days (any kit format).

7. SAMPLES

Total DNA can be extracted from urethral swab specimens, endocervical and vaginal smears, urine samples, prostatic fluid specimens and other human body fluids and tissue samples. All methods of extraction can be used.

8. ASSAY PROCEDURE

Negative Controls must be included in each run.

The complete procedure consists of three stages:

1. DNA Extraction
2. Real-Time PCR
3. Data analysis and interpretation

1. DNA extraction

NOTE: Application of Astra Biotech DNA extraction kits for specimen preparation is strongly recommended.

1.1 Take **IC** and place the rest components back in refrigerator. Gently vortex **IC** and centrifuge for 2-3 s.

1.2 Prepare and mark required quantity of 1.5 mL microcentrifuge tubes equal to the number of samples + 1 tube for negative extraction control and label as “NEC”.

NOTE: Each extraction procedure should be performed using “NEC”.

1.3 If Astra Biotech DNA extraction kit is used, ignore steps 1.4 –1.6 and follow the Astra Biotech DNA extraction kit manual; otherwise continue from step 1.4.

1.4 Place **20 µL IC** into each 1.5 mL microcentrifuge tube.

1.5 Add **100 µL NC** into the microcentrifuge tube marked “**NEC**”.

NOTE: The Elution buffer contained in DNA extraction kit could be used instead of **NC**.

1.6 Follow the manual of respective DNA extraction kit.

NOTE: Elution volume should not exceed 200 µL.

2. Real-Time PCR

2.1 Kit components preparation:

- Take the kit out from the refrigerator. Gently vortex all liquid components and centrifuge briefly.

2.2 Microcentrifuge tubes preparation:

- Prepare and mark required quantity of 0.2 mL microcentrifuge tubes/ wells equal to the number of samples + 1 tube/ well for “PC” + 1 tube/ well for “NEC”.

NOTE: Different positions for marking the tubes are necessary because of differences in the fluorescence detection: When using CFX96, CFX96Touch or SaCycler Real-Time PCR cyclers mark microcentrifuge tubes outside the centre of tube cap or on the tube side. When using Rotor-Gene 3000/6000 or Rotor-Gene Q instruments mark the tubes on the caps.

- Place **10 µL Taq** into each reaction tube avoiding tip contact with paraffin.

- Add **10 µL DNA** (extracted from specimens and “NEC”) avoiding tip contact with paraffin. Use separate filter tip for each sample.
- Add **10 µL positive control** **PC** in the **PCR tube marked “PC”**
- Close the tubes.
- Centrifuge at **1500-2400 rpm for 2-3 s.**

2.3 Real-Time PCR Amplification:

- Place tubes/ plate into the Real-Time PCR cycler.
- Complete the plate information. Set up fluorescent detection from channels ROX/ Orange and FAM/ Green.
- Run PCR protocol for all thermal cyclers:

Stage	Temperature, °C	Data collection	Time	Number of cycles
Hold	94	-	3 min	1
Cycling	94	-	10 s	5
	60		20 s	
Cycling	94	-	10 s	45
	60	ROX/ Orange FAM/ Green	20 s	

3. Data Analysis

3.1 General conditions:

Quantification cycles (= Cq) are determined for both FAM and ROX dyes in all samples. ROX dye is used to detect the *Trichomonas vaginalis* DNA and FAM is used to detect **IC** DNA.

The results of single run are valid if Cq_{ROX} "PC" is within the range determined in the quality control sheet to the product.

A contamination is present if Cq_{ROX} "NEC" < Cq TV (see quality control sheet). In that case all the specimens tested as positive in this run should be considered invalid. Process all these specimens again from the DNA extraction stage.

3.2 Data interpretation:

Cq _{ROX/ Orange} Detection of TV DNA	Cq _{FAM/Green} Detection of IC DNA	Interpretation
> Cq TV or N/A	≤ Cq IC	Negative
≤ Cq TV	ANY	Positive
> Cq TV or N/A	> Cq IC or N/A	Invalid, run the original specimen again or process new collected one

Abbreviations:

Cq TV= critical value of quantification cycle for Trichomonas vaginalis

*Cq **IC** = critical value of quantification cycle for Internal Control*

N/A = Cq value cannot be determined

9. PERFORMANCE CHARACTERISTICS OF THE KIT

Analytical sensitivity:

The detection limit of the kit evaluated via probit analysis is 480 copies of *Trichomonas vaginalis* genomic DNA per 1 mL. This limit is defined for confidence level of 99%.

Analytical specificity:

A false-positive reaction is absent with standard human DNA samples and standard DNA samples of *Trichomonas tenax*, *Trichomonas hominis*, *Trichomonas foetus*, *Neisseria gonorrhoeae*, *Toxoplasma gondii*, *Chlamydia trachomatis*, *Mycoplasma* sp., *Ureaplasma* sp., *Candida albicans*, *Enterococcus faecalis*, *Gardnerella vaginalis*, *Staphylococcus aureus*, *Lactobacillus acidophilus*.

10. LIMITATIONS OF THE METHOD

Any clinical diagnosis should not be based on the results of *in vitro* diagnostic methods alone. For diagnosis establishment a physician is supposed to consider all the available clinical and laboratory findings.

11. SAFETY PRECAUTIONS

- **This kit is for in vitro diagnostic use only.** The operator should thoroughly follow the manual to obtain reliable data. This instruction manual is valid only for the present kit with the listed contents. Any exchange of the kit components is not allowed by CE regulations.

- Test should only be performed by skilled personnel considering GLP (Good Laboratory Practice) guidelines.
- Don't use the kit after its expiration date.
- Do not pool reagents from different lots or from different bottles of the same lot (the only exception for **IC** is described below). Immediately after use, close all bottles in order to avoid leakage. After opening, store all bottles and vials in an upright position.
- In case of simultaneous using of two or more Astra Biotech Infectious Real-Time PCR kits for detection of different infectious agents, **IC** with the latest lot number can be used for all of these kits.
- We recommend setting up blank protocols of PCR before the analysis.
- PCR technology is extremely sensitive. The amplification of a single DNA molecule generates millions of identical copies. Therefore, set up three separate working areas for a) sample preparation b) PCR reagent preparation and c) PCR amplification. For each working area a different set of pipettes and protective clothing should be reserved.
- Use sterile filter tips for pipetting and use special PCR pipettes for aerosol free pipetting.
- Do not use the same tip for two different components, neither DNA-free nor DNA containing.
- Do not eat, drink or smoke in the laboratory work area.

- All components of ***Trichomonas vaginalis* PCR kit** contain 0.05 % of sodium azide, which is very toxic and dangerous to environment:
 - P260- Do not breathe dust/fume/gas/mist/vapours/spray;
 - P262- Do not get in eyes, on skin, or on clothing;
 - P264- Wash hands thoroughly after handling;
 - P270- Do not eat, drink or smoke when using this product;
 - P271- Use only outdoors or in a well-ventilated area;
 - P273- Avoid release to the environment;
 - P280- Wear protective gloves/protective clothing/eye protection/ face protection;
 - P284- Wear respiratory protection;
 - P310- Immediately call a POISON CENTER or doctor/physician;
 - P320- Specific treatment is urgent (see instruction for use);
 - P321- Specific treatment (see instruction for use);
 - P322- Specific measures (see instruction for use);
 - P330- Rinse mouth;
 - P361- Remove/Take off immediately all contaminated clothing;
 - P363- Wash contaminated clothing before reuse;
 - P 391- Collect spillage;
 - P405- Store locked up;
 - P501- Dispose of contents/container in accordance with national regulation;
 - P301 +P310- IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician;

- P302 +P350- IF ON SKIN: Gently wash with plenty of soap and water;
 - P304 +P340- IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing;
 - P403 +P233- Store in a well-ventilated place. Keep container tightly closed.
-
- Be careful while handling Positive Control, avoid its splashing.
 - Kits contaminated with *Trichomonas vaginalis* DNA/ **PC** should be discarded in a biological waste container.
 - Routinely decontaminate your pipettes and the laboratory benches.
 -  During DNA extraction of clinical materials, **biosafety level 2 practices and procedures** have to be followed. Human blood, blood products, body fluids and tissues should be treated as potential infectious.
 - After usage, all the reagents and test components should be discarded as consisted with local statements.

June 22, 2017

Astra Biotech GmbH
Rudower Chaussee 29
12489 Berlin, Germany
Telefon: +49 (0)30 74 696 509
E-Mail: info@astrabiotech.de
www.astrabiotech.de