

PRODUCT DATA SHEET

Description:	Gastrointestinal Tumor Antigen CA19-9 (Cancer antigen 19-9), Human
REF	11-03
Source	SW 1116 colon adenocarcinoma supernatants
Format	Purified, liquid
Purification	Gel-filtration chromatography after ammonium sulfate precipitation
Concentration	Lot specific (U/mL)
Buffer	0.15 M phosphate buffered saline, pH 7.4
Preservative	0.01 % NaN ₃
Handling	Avoid freeze/thaw cycles
Long term storage	-18...-22°C
Expiration date	Specified on label