

**PRODUCT DATA SHEET**

<b>Description</b>	<b>Monoclonal antibody to hCG</b>	
<b>REF</b>	<b>10-06</b>	
Specificity	Reacts with human Chorionic Gonadotropin	
Clone	11C9, IgG1	
Host Animal	Mouse	
Immunogen	Human hCG	
Cross-reactivity %	TG	<0.1%
	TSH	< 0.5%
	LH does not exceed	<1.3%
Source	Ascitic fluid	
Format	Purified, Liquid	
Purification	Chromatography on Protein G Sepharose, purity is tested by electrophoresis	
Concentration	lot-specific	
Buffer	PBS, pH 7.4	
Preservative	0.02% Sodium azide	
Long term storage	-18...-22°C	
Handling	Avoid freeze/thaw cycles	
Application	ELISA, use as detection antibody in ELISA together with Mab 4B4	